

PRESTO® W80

Heating a 5 liters reactor from -60 °C to +20 °C

Objective

This case study tests the heating power of PRESTO® W80 with a 5 l glass reactor. The PRESTO® W80 is connected to the reactor via two 1 m metal tubings. The PRESTO® W80 is programmed to heat up from -60 °C to +20 °C.



Room temperature +20 °C Humidity 45%

Voltage 230 V / 50 Hz



Test Conditions

JULABO unit PRESTO® W80
Cooling power +20 °C 1.2 kW

0 °C 1.2 kW -20 °C 1.1 kW

Heating capacity

Band limit

Flow pressure

Bath fluid

1.8 kW

without

0.4 bar

Thermal

Bath fluid Thermal HL 80
Reactor 5 I glass reactor (Rettberg)

filled with E I Thermal III. Of

filled with 5 I Thermal HL 80

Jacket volume 2.5 l

Control External (ICC)








Test Results

The PRESTO® W80 heating process from -60 °C to +20°C in 55 min without overshoot.



Setpoint
Temperature in reactor's interior
Temperature in reactor's jacket



