

PRESTO[®] A40 Heating a 6 liters reactor from -20 °C to +20 °C

Objective

0

This case study tests the heating power of PRESTO[®] A40 with a 6 liters glass reactor. The PRESTO[®] A40 is connected to the reactor via two 2 m metal tubings. The PRESTO[®] A40 is programmed to heat up from -20 °C to +20 °C.



Room temperature	+20 °C
Humidity	45%
Voltage	230 V / 50 Hz


Test Conditions

JULABO unit Cooling power

Heating capacity Band limit Flow pressure Bath fluid Reactor

Jacket volume Control +20 °C 1.2 kW 0 °C 0.9 kW -20 °C 0.6 kW 2.7 kW without 0.5 bar Thermal HL60 6 liters glass reactor (QVF) filled with 5 I Thermal HL60 4.5 I External (ICC)

PRESTO® A40






Test Results

0

The PRESTO[®] A40 heating process from -20 °C to +20 °C in 30 min without overshoot.



Setpoint
Temperature in reactor's interior
Temperature in reactor's jacket

Тір

KEY

Use our tube adapters and your tubing will no longer kink.



Tip You can also use the robust Pt100 with PTFE coating.

