

PRESTO® A40

Heating a 6 liters reactor from +20 °C to +100 °C

Objective

This case study tests the heating power of PRESTO® A40 with a 6 liters glass reactor. The PRESTO® A40 is connected to the reactor via two 2 m metal tubings. The PRESTO® A40 is programmed to heat up from +20 °C to +100 °C.

Environment

Room temperature +20 °C Humidity 45%

Voltage 230 V / 50 Hz



Test Conditions

JULABO unit PRESTO® A40 Cooling power +20 °C 1.2 kW 0 °C 0.9 kW

-20 °C 0.6 kW


Heating capacity 2.7 kW
Band limit without
Flow pressure 0.5 bar
Bath fluid Thermal HL60

Reactor 6 liters glass reactor (QVF)

filled with 5 I Thermal HL60

Jacket volume 4.5 l

Control External (ICC)








Test Results

The PRESTO® A40 heating process from +20 °C to +100°C in 1 h without overshoot.



Setpoint

Temperature in reactor's interior

Temperature in reactor's jacket

Tip
Protect your reactor.
The function "band limit" (see above) permits setting the max. temperature difference between jacket and internal vessel.

Profile of reactor

I have been profile of reactor

