

PRESTO® A40

Cooling a 6 liters reactor from +100 °C to +20 °C

Objective

This case study tests the cooling power of PRESTO® A40 with a 6 liters glass reactor. The PRESTO® A40 is connected to the reactor via two 2 m metal tubings. The PRESTO® A40 is programmed to cool down from +100 °C to +20 °C.

Environment

Room temperature +20 °C Humidity 45 %

Voltage 230 V / 50 Hz



Test Conditions

JULABO unit PRESTO® A40 Cooling power +20 °C 1.2 kW

0 °C 0.9 kW -20 °C 0.6 kW

Heating capacity 2.7 kW
Band limit without
Flow pressure 0.5 bar
Bath fluid Thermal HL60

Reactor 6 l glass reactor (QVF)

filled with 5 I Thermal HL60

Jacket volume 4.5 l

Control External (ICC)








Test Results

The PRESTO® A40 cooling process from +100 °C to +20 °C in 1 h without overshoot.





Take advantage of our wide range of accessories. The M+R adapter enables you to display and record an additional temperature.



Tip

You can also use the robust Pt100 with PTFE coating.

