

PRESTO® A40

Heating a 6 liters reactor from +20 °C to +150 °C

Objective

This case study tests the heating power of PRESTO® A40 with a 6 liters glass reactor. The PRESTO® A40 is connected to the reactor via two 2 m metal tubings. The PRESTO® A40 is programmed to heat up from +20 °C to +150 °C.



Room temperature +20 °C Humidity 45%

Voltage 230 V / 50 Hz



Test Conditions

JULABO unit PRESTO® A40 Cooling power +20 °C 1.2 kW 0 °C 0.9 kW

-20 °C 0.6 kW

Heating capacity

Band limit

Flow pressure

Bath fluid

2.7 kW

without

0.5 bar


Thermal HL60

Reactor 6 liters glass reactor (QVF)

filled with 5 I Thermal HL60

Jacket volume 4.5 l

Control External (ICC)








Test Results

The PRESTO® A40 heating process from +20 °C to +150°C in 1 h 20 min without overshoot.



Setpoint
Temperature in reactor's interior
Temperature in reactor's jacket

