

## PRESTO® A40

# Cooling a 6 liters reactor from +150 °C to +20 °C

### **Objective**

This case study tests the cooling power of PRESTO® A40 with a 6 liters glass reactor. The PRESTO® A40 is connected to the reactor via two 2 m metal tubings. The PRESTO® A40 is programmed to cool down from +150 °C to +20 °C.

#### **Environment**

Room temperature +20 °C Humidity 45 %

Voltage 230 V / 50 Hz



#### **Test Conditions**

JULABO unit PRESTO® A40 Cooling power +20 °C 1.2 kW

0 °C 0.9 kW -20 °C 0.6 kW

Heating capacity 2.7 kW
Band limit without
Flow pressure 0.5 bar
Bath fluid Thermal HL60

Reactor 6 l glass reactor (QVF)

reactor 0 i glass reactor (QVI)

filled with 5 I Thermal HL60

Jacket volume 4.5 l

Control External (ICC)








#### **Test Results**

The PRESTO® A40 cooling process from +150 °C to +20 °C in 1 h 15 min without overshoot.





nal temperature.

