

# **PRESTO W50 Cooling a 50 liters reactor** from +150 °C to +20 °C

### **Objective**

0

This case study tests the cooling power of PRESTO W50 with a 50 liters glass reactor. The PRESTO W50 is connected to the reactor via two 2 m metal tubings. The PRESTO W50 is programmed to cool down from +150 °C to +20 °C.

### Environment

+20 °C Room temperature Humidity Voltage

45 % 400 V / 50 Hz

#### **Test Conditions**

JULABO unit Cooling power	PRESTO W50 +20 °C 7.5 kW 0 °C 65 kW	
	-20 °C 3.0 kW	
Heating capacity	6 kW	
Band limit	without	
Flow pressure	0.5 bar	
Bath fluid	Thermal HL60	
Reactor	50 l glass reactor (QVF)	
	filled with 35 I Thermal HL60	
Jacket volume	26.5 l	
Control	External (ICC)	

### **Control Parameters**

Хр	0.2 K	
Tn	695 s	
Tv	85 s	
Хри	15 K	
How is a		
ſ,¥¥-		
~~~	20.1	
	The states	
	P AN	#






## www.julabo.com



#### **Test Results**


0

The PRESTO W50 cooling process from +150 °C to +20 °C in 1 h 10 min without overshoot.



#### Tip

Protect your reactor. The function "band limit" (see above) permits setting the max. temperature difference between jacket and internal vessel.



# Тір

Use our tube adapters and your tubing will no longer kink.

