

PRESTO A85t

Cool-down a 50 liters reactor from +20 °C to lowest possible temperature

Objective

This case study tests the lowest possible temperature of the PRESTO A85t with a 50 liters glass reactor. The PRESTO A85t is connected to the reactor via 2 m metal tubings. The PRESTO A85t cools down from +20 °C to the lowest possible temperature.

Environment

Room temperature +20 °C Humidity 45 % Voltage 400 V / 50 Hz

Test Conditions

JULABO unit PRESTO A85t
Cooling power +20 °C 2.5 kW

0 °C 2.4 kW -20 °C 2.4 kW

Heating capacity 15 kW
Band limit without
Flow pressure 0.5 bar
Bath fluid Thermal HL80

Reactor 50 l glass reactor (Büchiglas)

filled with 50 l Thermal HL80

Jacket volume 11 l


Control External (ICC)

Control Parameters

Xp 0.2 K Tn 720 s Tv 100 s Xpu 24 K









Test Results

The PRESTO $^{\text{TM}}$ A85t cooled the reactor from +20 °C down to the lowest possible temperature in 6 h 30 min. Within these test conditions the lowest possible temperature is -73 °C.

