

PRESTO® A30

Cool-down a 6 liters reactor from +20 °C to lowest possible temperature

Objective

This case study tests the lowest possible temperature of the PRESTO® A30 with a 6 liters glass reactor. The PRESTO® A30 is connected to the reactor via two 2 m metal tubings. The PRESTO® A30 cools down from +20 °C to the lowest possible temperature.

Environment

Room temperature +20 °C Humidity 45 %

Voltage 230 V / 50 Hz



Test Conditions

JULABO unit PRESTO® A30 Cooling power +20 °C 0.5 kW

0 °C 0.4 kW -20 °C 0.2 kW

Heating capacity 2.7 kW
Band limit without
Flow pressure 0.5 bar
Bath fluid Thermal HL60

Reactor 6 l glass reactor (QVF)

filled with 5 l Thermal HL60

Jacket volume 4.5 l

Control External (ICC)








Test Results

The PRESTO® A30 cooled the reactor from +20 °C down to the lowest possible temperature in 4 h. Within these test conditions the lowest possible temperature is -12 °C.



Setpoint
Temperature in reactor's interior
Temperature in reactor's jacket

Tip Take advantage of our wide range of accessories. The M+R adapter enables you to display and record an additional temperature.



