

PRESTO[®] A30 Cooling a 6 liters reactor from +50 °C to +20 °C

Objective

0

This case study tests the cooling power of PRESTO® A30 with a 6 liters glass reactor. The PRESTO® A30 is connected to the reactor via two 2 m metal tubings. The PRESTO® A30 is programmed to cool down from +50 °C to +20 °C.

Environment

Room temperature	+20 °C
Humidity	45 %
Voltage	230 V / 50 Hz

Test Conditions

JULABO unit Cooling power

Heating capacity Band limit Flow pressure Bath fluid Reactor

Jacket volume Control +20 °C 0.5 kW 0 °C 0.4 kW -20 °C 0.2 kW 2.7 kW without 0.5 bar Thermal HL60 6 I glass reactor (QVF) filled with 5 I Thermal HL60 4.5 I External (ICC)

PRESTO® A30








Test Results

0


The PRESTO® A30 cooling process from +50 °C to +20 °C in 50 min without overshoot.



Tip

an and

Use the free of charge EasyTEMP software to control the units with the PC and to show the temperature curves graphically.



Tip

Take advantage of our wide range of accessories. The M+R adapter enables you to display and record an additional temperature.

