

PRESTO® A30

Cooling a 6 liters reactor from +100 °C to +20 °C

Objective

This case study tests the cooling power of PRESTO® A30 with a 6 liters glass reactor. The PRESTO® A30 is connected to the reactor via two 2 m metal tubings. The PRESTO® A30 is programmed to cool down from +100 °C to +20 °C.

Environment

Room temperature +20 °C Humidity 45 %

Voltage 230 V / 50 Hz



Test Conditions

JULABO unit PRESTO® A30
Cooling power +20 °C 0.5 kW

0 °C 0.4 kW -20 °C 0.2 kW

Heating capacity 2.7 kW
Band limit without
Flow pressure 0.5 bar
Bath fluid Thermal HL60

Reactor 6 l glass reactor (QVF)

filled with 5 I Thermal HL60

Jacket volume 4.5 l

Control External (ICC)







Test Results

The PRESTO® A30 cooling process from +100 °C to +20 °C in 1 h 20 min without overshoot.





