

PRESTO® A30

Cooling a 6 liters reactor from +150 °C to +20 °C

Objective

This case study tests the cooling power of PRESTO® A30 with a 6 liters glass reactor. The PRESTO® A30 is connected to the reactor via two 2 m metal tubings. The PRESTO® A30 is programmed to cool down from +150 °C to +20 °C.

Environment

Room temperature +20 °C Humidity 45 %

Voltage 230 V / 50 Hz



Test Conditions

JULABO unit PRESTO® A30
Cooling power +20 °C 0.5 kW

0 °C 0.4 kW -20 °C 0.2 kW

Heating capacity 2.7 kW
Band limit without
Flow pressure 0.5 bar
Bath fluid Thermal HL60

Reactor 6 l glass reactor (QVF)

filled with 5 I Thermal HL60

Jacket volume 4.5 l

Control External (ICC)








Test Results

The PRESTO® A30 cooling process from +150 °C to +20 °C in 1h 35 min without overshoot.



Tip

Take advantage of our wide range of accessories. The M+R adapter enables you to display and record an additional temperature.



Tip
Protect your reactor.
The function "band limit" (see above) permits setting the max. temperature difference between jacket and internal vessel.

Profile of reactor

I a setting the max. temperature difference between jacket and internal vessel.