

PRESTO® A40

Cooling a 6 liters reactor from +20 °C to 0 °C

Objective

This case study tests the cooling power of PRESTO® A40 with a 6 liters glass reactor. The PRESTO® A40 is connected to the reactor via two 2 m metal tubings. The PRESTO® A40 is programmed to cool down from +20 °C to 0 °C.

Environment

Room temperature +20 °C Humidity 45 %

Voltage 230 V / 50 Hz



Test Conditions

JULABO unit PRESTO® A40 Cooling power +20 °C 1.2 kW

0 °C 0.9 kW -20 °C 0.6 kW


Heating capacity 2.7 kW
Band limit without
Flow pressure 0.5 bar
Bath fluid Thermal HL60

Reactor 6 l glass reactor (QVF)

filled with 5 I Thermal HL60

Jacket volume 4.5 l

Control External (ICC)








Test Results

The PRESTO® A40 cooling process from +20 °C to 0 °C in 38 min without overshoot.



Tip

Use the free of charge *EasyTEMP* software to control the units with the PC and to show the temperature curves graphically.



Tip

Take advantage of our wide range of accessories. The M+R adapter enables you to display and record an additional temperature.

