

PRESTO[®] A40 Heating a 6 liters reactor from 0 °C to +20 °C

Objective

0

This case study tests the heating power of PRESTO[®] A40 with a 6 liters glass reactor. The PRESTO[®] A40 is connected to the reactor via two 2 m metal tubings. The PRESTO[®] A40 is programmed to heat up from 0 °C to +20 °C.



Room temperature	+20 °C
Humidity	45%
Voltage	230 V / 50 Hz

Test Conditions

JULABO unit Cooling power

Heating capacity Band limit Flow pressure Bath fluid Reactor

Jacket volume Control +20 °C 1.2 kW 0 °C 0.9 kW -20 °C 0.6 kW 2.7 kW without 0.5 bar Thermal HL60 6 liters glass reactor (QVF) filled with 5 I Thermal HL60 4.5 I External (ICC)

PRESTO® A40







Test Results

0

The PRESTO® A40 heating process from 0 °C to +20°C in 18 min without overshoot.



Setpoint
Temperature in reactor's interior
Temperature in reactor's jacket

Tip

KUN

Protect your reactor. The function "band limit" (see above) permits setting the max. temperature difference between jacket and internal vessel.



Tip

Take advantage of our wide range of accessories. The M+R adapter enables you to display and record an additional temperature.

