

## PRESTO<sup>®</sup> A40 Heating a 6 liters reactor from -10 °C to +20 °C

## Objective

0

This case study tests the heating power of PRESTO<sup>®</sup> A40 with a 6 liters glass reactor. The PRESTO<sup>®</sup> A40 is connected to the reactor via two 2 m metal tubings. The PRESTO<sup>®</sup> A40 is programmed to heat up from -10 °C to +20 °C.



Room temperature	+20 °C
Humidity	45%
Voltage	230 V / 50 Hz


## **Test Conditions**

JULABO unit Cooling power

Heating capacity Band limit Flow pressure Bath fluid Reactor

Jacket volume Control +20 °C 1.2 kW 0 °C 0.9 kW -20 °C 0.6 kW 2.7 kW without 0.5 bar Thermal HL60 6 liters glass reactor (QVF) filled with 5 I Thermal HL60 4.5 I External (ICC)

PRESTO® A40








## **Test Results**

0

The PRESTO® A40 heating process from -10 °C to +20°C in 24 min without overshoot.



Setpoint
Temperature in reactor's interior
Temperature in reactor's jacket

Tip

ったい

You can also use the robust Pt100 with PTFE coating.



**Tip** Use our tube adapters and your tubing will no longer kink.

