# Case Study JULABO F33-MA

Cooling Rate Determination of a JULABO F33-MA Refrigerated Circulator

Case Study
12-2009
Page 1 from 2

# **Objective**

Performance test of a JULABO F33-MA refrigerated circulator. In addition a horizontal uniformity measurement was conducted in the bath.

## **Test Details**

A cooling process from +50 °C to -10 °C was performed to determine the cooling performance.

To conduct the uniformity measurement two identical temperature sensors were placed at the same depth in the center of the bath in order to record the temperature difference between the two sensors. In addition a uniformity comparison was conducted between one sensor in the middle and one in each corner of the bath.

The measurement results were recorded on a notebook using the EasyTemp software.

# **Test Conditions**

Equipment F33-MA
Mains Voltage 230 Volt / 50 Hz
Ambient temperature 20-22 °C (room temperature)

Bath fluid Silicon-Oil "Thermal HY"
Specific Settings Temperature control via
JULABO EasyTemp Software



### **RESULT**

The test result shows that the cooling process was completed within 58 minutes.


The uniformity in comparison to the sensor in the middle of the bath is +0.09K (maximum).

# **Test Results**

TEST		TEMPERATURE RANGE	TIME
Test 1	Cool-down	+50 °C10 °C	58 Min.

### Test 1

Cool-down from +50°C to -10°C in 58 minutes.







# Case Study JULABO F33-MA




**Test 2**Horizontal uniformity measurement.
Measurement setup and calibration of the two sensors.



The measurements were taken in all 4 corners of the bath. The deviations of the temperatures in the corners were stated in relation to the sensor in the middle of the bath. Both sensos had the same immersion depth during themeasurement.

Sensor position	Deviation	
Back left	-0,02k	
Front left	+0,08k	
Back right	+0,08k	
Back front	+0,09k	



### JULABO GmbH

Eisenbahnstraße 45 • 77960 Seelbach / Germany Tel. +49 (0) 7823 51-0 • Fax +49 (0) 7823 24 91 info@julabo.de • www.julabo.de



