Juliubo Case Study

JULABO FP50-HL

Cooling a 10 liters reactor from +80 °C to +20 °C



Objective

This case study tests the cooling power of **JULABO FP50-HL** with a 10 liters glass reactor. The FP50-HL is connected to the reactor via two 2 m metal tubings. The FP50-HL is programmed to cool down from +80 °C to +20°C.

Test Conditions

JULABO unit JULABO FP50-HL Cooling power $+20~^{\circ}\text{C}$ 0.9 kW $0~^{\circ}\text{C}$ 0.8 kW

-20 °C 0.5 kW

Heating capacity 2 kW
Band limit without
Flow pressure 0,4 bar

Bath fluid JULABO Thermal H10

Reactor 10 liters glass reactor (Normag)

filled with 10 liter JULABO Thermal H10

Jacket volume 5

Control External (ICC)

Environment

Room temperature $20 \, ^{\circ}\text{C}$ Humidity $45 \, \%$

Voltage 230 V /50 Hz



Test Results

See chart on back page: The FP50-HL cooling process from +80 °C to +20°C in 95 min without overshoot.

Tip


You can also use the robust Pt100 with PTFE coating.

More tips on back page >>




JULABO GmbH Eisenbahnstraße 45 77960 Seelbach / Germany Tel. +49 (0) 7823 51-0





SetpointTemperature in reactor's interiorTemperature in reactor's jacket

TipUse the free of charge *EasyTEMP* software to control the units with the PC and to show the temperature curves graphically.



JULABO GmbH Eisenbahnstraße 45 77960 Seelbach / Germany Tel. +49 (0) 7823 51-0

