

Julabo Case Study

JULABO FP50-HL

Cooling a 10 liters reactor
from +20°C to -20 °C



Objective


This case study tests the cooling power of **JULABO FP50-HL** with a 10 liters glass reactor. The FP50-HL is connected to the reactor via two 2 m metal tubings. The FP50-HL is programmed to cool down from +20 °C to -20 °C.

Test Conditions

JULABO unit	JULABO FP50-HL
Cooling power	+20 °C 0.9 kW
	0 °C 0.8 kW
	-20 °C 0.5 kW
Heating capacity	2 kW
Band limit	without
Flow pressure	0,4 bar
Bath fluid	JULABO Thermal H10
Reactor	10 liters glass reactor (Normag) filled with 10 liter JULABO Thermal H10
Jacket volume	5 l
Control	External (ICC)

Environment

Room temperature	20 °C
Humidity	45 %
Voltage	230V / 50 Hz



Test Results

See chart on back page: The FP50-HL cooling process from +20°C to -20 °C in 2 h 30 min without overshoot.


Tip

You can also use the robust Pt100 with PTFE coating.

More tips on back page >>




JULABO GmbH
Eisenbahnstraße 45
77960 Seelbach / Germany
Tel. +49 (0) 7823 51-0



Tip

Use the free of charge *EasyTEMP* software to control the units with the PC and to show the temperature curves graphically.



JULABO GmbH
Eisenbahnstraße 45
77960 Seelbach / Germany
Tel. +49 (0) 7823 51-0