Julium Case Study

JULABO FP50-HL

Heating a 10 liters reactor from -20 °C to +20 °C



Objective

This case study tests the heating power of **FP50-HL** with a 10 liters glass reactor. The FP50-HL is connected to the reactor via two 2 m metal tubings. The FP50-HL is programmed to heat up from -20 $^{\circ}$ C to +20 $^{\circ}$ C.

Test Conditions

JULABO unit JULABO FP50-HL Cooling power +20 °C 0.9 kW

0 °C 0.8 kW -20 °C 0.5 kW

Heating capacity 2 kW
Band limit without
Flow pressure 0.4 bar

Bath fluid JULABO Thermal H10

Reactor 10 liters glass reactor (Normag)

filled with 10 liter JULABO Thermal H10

Jacket volume 5.0 l

Control External (ICC)

Environment

Room temperature 20 °C Humidity 45 %

Voltage 230 V / 50 Hz



Test Results

See chart on back page: The FP50-HL heating process from -20 $^{\circ}$ C to +20 $^{\circ}$ C in 70 min without overshoot.

Tip


You can also use the robust Pt100 with PTFE coating.


More tips on back page >>




JULABO GmbH Eisenbahnstraße 45 77960 Seelbach / Germany Tel. +49 (0) 7823 51-0







TipUse the free of charge *EasyTEMP* software to control the units with the PC and to show the temperature curves graphically.



JULABO GmbH Eisenbahnstraße 45 77960 Seelbach / Germany Tel. +49 (0) 7823 51-0

