Juliaho Case Study

JULABO FP50-HL

Cooling a 10 liters reactor from +120°C to -20 °C



Objective

This case study tests the cooling power of **JULABO FP50-HL** with a 10 liters glass reactor. The FP50-HL is connected to the reactor via two 2 m metal tubings. The FP50-HL is programmed to cool down from +120 °C to -20 °C.

Test Conditions

JULABO unit JULABO FP50-HL Cooling power +20 °C 0.9 kW 0 °C 0.8 kW

-20 °C 0.5 kW

Heating capacity 2 kW
Band limit without
Flow pressure 0,4 bar

Bath fluid JULABOThermal H10

Reactor 10 liters glass reactor (Normag)

filled with 10 liter JULABO Thermal H10

Jacket volume 5

Control External (ICC)

Environment

Room temperature $20 \, ^{\circ}\text{C}$ Humidity $45 \, \%$

Voltage 230V / 50 Hz





Test Results

See chart on back page: The FP50-HL cooling process from $+120^{\circ}$ C to -20° C in 3 h 25 min without overshoot.

JULABO GmbH Eisenbahnstraße 45 77960 Seelbach / Germany Tel. +49 (0) 7823 51-0









JULABO GmbH Eisenbahnstraße 45 77960 Seelbach / Germany Tel. +49 (0) 7823 51-0

