

JULABO FP50-HL

Temperature stability with a 10 l reactor at +20 °C

Objective

This case study tests the temperature stability of **JULABO FP50-HL** with a 10 liters glass reactor. The FP50-HL is connected to the reactor via two 2 m metal tubings. The FP50-HL was set to setpoint +20 °C.

JULABO FP50-HL +20 °C 0.9 kW

Test Conditions

JULABO unit Cooling power

Heating capacity Band limit Flow pressure Bath fluid Reactor Jacket Volume

Control

0 °C 0.8 kW -20 °C 0.5 kW 2 kW without 0.4 bar JULABO Thermal H10 10 liters glass reactor (Normag) filled with 10 liter JULABO Thermal H10 5.0 l External (ICC)

Test Results

See chart on back page: The FP50-HL cooled down the reactor to +20 °C. After reaching the temperature of +20 °C, the temperature within the reactor fluctuated for 10 min about \pm 0.01 K.

Environment

Room temperature	20 °C
Humidity	45 %
Voltage	230 V / 50 Hz



Tip You can also use the robust Pt100 with PTFE coating.




More tips on back page >>

> JULABO GmbH Eisenbahnstraße 45 77960 Seelbach / Germany Tel. +49 (0) 7823 51-0




www.julabo.de





Tip

Use the free of charge *EasyTEMP* software to control the units with the PC and to show the temperature curves graphically.



JULABO GmbH Eisenbahnstraße 45 77960 Seelbach / Germany Tel. +49 (0) 7823 51-0



www.julabo.de