Juliubio Case Study

JULABO FP50-HL

Cool-down a 10liters reactor from +20 °C to maximum low temperature



Objective

This case study tests the maximum low temperature of the **JULABO FP50-HL** with a 10 liters glass reactor. The FP50-HL is connected to the reactor via 2 m metal tubings. The FP50-HL cools down from +20 °C to maximum low temperature.

Test Conditions

JULABO unit JULABO FP50-HL Cooling power $+20~^{\circ}\text{C}$ 0.9 kW $0~^{\circ}\text{C}$ 0.8 kW

-20 °C 0.5kW

Heating capacity 2 kW
Band limit without
Flow pressure 0.4 bar

Bath fluid JULABO Thermal H10

Reactor 10 liters glass reactor (Normag)

filled with 10 liter JULABO Thermal H10


Jacket volume 5.0 l

Control External (ICC)

Environment

Room temperature $20 \, ^{\circ}\text{C}$ Humidity $45 \, \%$

Voltage 230 V / 50 Hz



Test Results

See chart on back page: The FP50-HL cooled the reactor from $+20^{\circ}$ C down to maximum low temperature of -25 °C in 5 h.

Tip


You can also use the robust Pt100 with PTFE coating.

More tips on back page >>




JULABO GmbH Eisenbahnstraße 45 77960 Seelbach / Germany Tel. +49 (0) 7823 51-0





Tip Use the free of charge EasyTEMP software to control the units with the PC and to show the temperature curves graphically.



JULABO GmbH Eisenbahnstraße 45 77960 Seelbach / Germany Tel. +49 (0) 7823 51-0

