

## **JULABO Presto A30**

Cooling and heating a 5.5 liters reactor between +20 °C and -20 °C

### **Objective**

This case study tests the heating and cooling power of JULABO Presto A30 with a 5.5 liters glass reactor. The A30 is connected to the reactor with two 1.0 m metal tubings. The A30 is programmed to cycle between +20 °C and -20 °C.

JULABO Presto A30

0 °C 0.4 kW

+20 °C 0.5 kW

## **Test Conditions**

JULABO unit Cooling power

Heating capacity
Band limit
Flow pressure
Bath fluid
Reactor

Control

-20 °C 0.2 kW 2.7 kW no 0.31 bar JULABO Thermal HL45 5.5 liters glass reactor (Bruno Kummer) filled with 5 l Thermal HL45 external (ICC)



#### Environment

Room temperature	+20 °C
Humidity	45 %
Voltage	230 V / 50 Hz



#### **Test Results**

See chart on back page: The A30 heating process from -20 °C to +20 °C in 1 h. Hitting exactly +20 °C without overshoot. The cooling process from +20 °C to -20 °C in 1 h 15 min. Hitting exactly -20 °C without overshoot.

#### Tip

Elbow fittings 90° helps relieving the connectors of the glass reactor.



More tips on back page >>

JULABO GmbH Eisenbahnstraße 45 77960 Seelbach / Germany Tel. +49 (0) 7823 51-0



#### www.julabo.de



## Tip

You can also use the robust Pt100 sensor with PTFE coating.



# Tip

The Ethernet interface permits full access to all operational functions of the PRESTO<sup>®</sup>.



JULABO GmbH Eisenbahnstraße 45 77960 Seelbach / Germany Tel. +49 (0) 7823 51-0



## www.julabo.de