Juliaho Case Study

JULABO Presto A30

Cooling and heating a 2 liters reactor between +20 °C and +100 °C



Objective

This case study tests the heating and cooling power of JULABO Presto A30 with a 2 litre glass reactor. The A30 is connected to the reactor with two 1.0 m metal tubings. The A30 is programmed to cycle between +20 °C and +100 °C.

Test Conditions

JULABO unit JULABO Presto A30 Cooling power +20 °C 0.5 kW 0 °C 0.4 kW

-20 °C 0.2 kW

Heating capacity 2.7 kW
Band limit no
Flow pressure 0.35 bar

Bath fluid JULABO Thermal HL45

Reactor 2 liters glass reactor (Schott Duran)

filled with 1.8 liters Thermal HL45

Control External (ICC)

Environment

Room temperature $20 \, ^{\circ}\text{C}$ Humidity $45 \, \%$

Voltage 230 V / 50 Hz




Test Results

See chart on back page: The A30 heating process from +20 °C to +100 °C in 27 min. Hitting exactly +100 °C without overshoot. The cooling process from +100 °C to +20 °C in 32 min. Hitting exactly +20 °C without overshoot.



JULABO GmbH Eisenbahnstraße 45 77960 Seelbach / Germany Tel. +49 (0) 7823 51-0







JULABO GmbH Eisenbahnstraße 45 77960 Seelbach / Germany Tel. +49 (0) 7823 51-0

