Juliaho Case Study

JULABO PRESTO® A40

Heating a 10 liters reactor from -25 °C to 0 °C



Objective

This case study tests the heating power of JULABO PRESTO® A40 with a 10 liters glass reactor. The A40 is connected to the reactor via two 2.0 m metal tubings. The A40 is programmed to heat up from -25 °C to 0 °C.

Environment

Room temperature +20 °C Humidity 45 %

Voltage 230 V / 50 Hz

Test Conditions

JULABO unit

Cooling power

+20 °C 1.2 kW

0 °C 0.9 kW

-20 °C 0.6 kW

Heating capacity 2.7 kW
Band limit No
Flow pressure 0.40 bar

Bath fluid JULABO Thermal HL40

Reactor 10 liters glass reactor (Normag)

filled with 10 liter JULABO Thermal HL40

Control External (ICC)


Test Results


See chart on back page: The A40 heating process from -25 $^{\circ}$ C to 0 $^{\circ}$ C in 30 min without overshoot.



JULABO GmbH Eisenbahnstraße 45 77960 Seelbach / Germany Tel. +49 (0) 7823 51-0









JULABO GmbH Eisenbahnstraße 45 77960 Seelbach / Germany Tel. +49 (0) 7823 51-0

