

JULABO PRESTO® A40

Heating a 5 liters reactor from +25 °C to +100 °C

Objective

This case study tests the heating power of JULABO PRESTO[®] A40 with a 5 liters glass reactor. The A40 is connected to the reactor via two 2.0 m metal tubings. The A40 is programmed to heat up from +25 °C to +100 °C.

JULABO PRESTO® A40

Test Conditions

JULABO unit Cooling power

Heating capacity Band limit Flow pressure Bath fluid Reactor

Test Results

See chart on back page: The A40 heating process from

+25 °C to +100 °C in 40 min without overshoot.

Control

+20 °C 1.2 kW 0 °C 0.9 kW -20 °C 0.6 kW 2.7 kW No 0.40 bar JULABO Thermal HL40 5 liters glass reactor (Rettberg) filled with 5 liter JULABO Thermal HL40 External (ICC)



Environment

Room temperature	+20 °C
Humidity	45 %
Voltage	230 V / 50 Hz



Tip You can also use the robust Pt100 with PTFE coating.


More tips on back page >>




JULABO GmbH Eisenbahnstraße 45 77960 Seelbach / Germany Tel. +49 (0) 7823 51-0



www.julabo.de





Tip

The Ethernet interface permits full access to all operational functions of the PRESTO[®].



JULABO GmbH Eisenbahnstraße 45 77960 Seelbach / Germany Tel. +49 (0) 7823 51-0



www.julabo.de