

JULABO PRESTO® A40

Temperature stability with a 5 l reactor at +100 °C

Objective

This case study tests the temperature stability of JULABO PRESTO[®] A40 with a 5 liters glass reactor. The A40 is connected to the reactor via two 2.0 m metal tubings. The A40 was set to a set point of +100 °C.

JULABO PRESTO® A40

Test Conditions

JULABO unit Cooling power

Heating capacity Band limit Flow pressure Bath fluid Reactor

Control

+20 °C 1.2 kW 0 °C 0.9 kW -20 °C 0.6 kW 2.7 kW No 0.40 bar JULABO Thermal HL40 5 liters glass reactor (Rettberg) filled with 5 liter JULABO Thermal HL40 External (ICC)

Test Results

See chart on back page: The A40 heats up the reactor to +100 °C. After reaching the temperature of +100 °C, the temperature within the reactor fluctuated for 10 min about ± 0.01 K.



Environment

Room temperature	+20 °C
Humidity	45 %
Voltage	230 V / 50 Hz



Tip You can also use the robust Pt100 with PTFE coating.


More tips on back page >>



JULABO GmbH Eisenbahnstraße 45 77960 Seelbach / Germany Tel. +49 (0) 7823 51-0



www.julabo.de



Tip

Make use of the option to regulate the pump pressure. You can define the desired pressure in the PRESTO® settings.



Tip

The Ethernet interface permits full access to all operational functions of the PRESTO[®].



JULABO GmbH Eisenbahnstraße 45 77960 Seelbach / Germany Tel. +49 (0) 7823 51-0



www.julabo.de