Juliubo Case Study

JULABO PRESTO® A40

Heating a 5 liters reactor from -25 °C to +200 °C



Objective

This case study tests the heating power of JULABO PRESTO® A40 with a 5 liters glass reactor. The A40 is connected to the reactor via two 2.0 m metal tubings. The A40 is programmed to heat up from -25 °C to \pm 200 °C.

Test Conditions

JULABO unit

Cooling power

JULABO PRESTO® A40

+20 °C 1.2 kW

+20 °C 1.2 kW 0 °C 0.9 kW

-20 °C 0.6 kW

Heating capacity 2.7 kW
Band limit No
Flow pressure 0.40 bar

Bath fluid JULABO Thermal HL40

Reactor 5 liters glass reactor (Rettberg)

filled with 5 liter JULABO Thermal HL40

Control External (ICC)

Environment

Room temperature +20 °C Humidity 45 %

Voltage 230 V / 50 Hz




Test Results


See chart on back page: The A40 heating process from -25 °C to +200 °C in 1 h 10 min without overshoot.



JULABO GmbH Eisenbahnstraße 45 77960 Seelbach / Germany Tel. +49 (0) 7823 51-0









JULABO GmbH Eisenbahnstraße 45 77960 Seelbach / Germany Tel. +49 (0) 7823 51-0

