Juliaho Case Study

JULABO PRESTO® A40

Heating a 5 liters reactor from 0 °C to +100 °C



Objective

This case study tests the heating power of JULABO PRESTO® A40 with a 5 liters glass reactor. The A40 is connected to the reactor via two 2.0 m metal tubings. The A40 is programmed to heat up from 0 °C to \pm 100 °C.

Test Conditions

JULABO unit JULABO PRESTO® A40

Cooling power +20 °C 1.2 kW

0 °C 0.9 kW

-20 °C 0.6 kW

Heating capacity 2.7 kW
Band limit No
Flow pressure 0.40 bar

Bath fluid JULABO Thermal HL40

Reactor 5 liters glass reactor (Rettberg)

filled with 5 liter JULABO Thermal HL40

Control External (ICC)

Environment

Room temperature +20 °C Humidity 45 %

Voltage 230 V / 50 Hz




Test Results

See chart on back page: The A40 heating process from $0 \,^{\circ}$ C to $+100 \,^{\circ}$ C in 35min without overshoot.




JULABO GmbH Eisenbahnstraße 45 77960 Seelbach / Germany Tel. +49 (0) 7823 51-0





SetpointTemperature in reactor's interiorTemperature in reactor's jacket





JULABO GmbH Eisenbahnstraße 45 77960 Seelbach / Germany Tel. +49 (0) 7823 51-0

