

# **JULABO PRESTO® A40**

Cooling a 20 liters reactor from +170 °C to +100 °C

# **Objective**

This case study tests the cooling power of JULABO PRESTO<sup>®</sup> A40 with a 20 liters vacuum insulated glass reactor. The A40 is connected to the reactor via two 2.0 m metal tubings. The A40 is programmed to cool down from +170 °C to +100 °C.

JULABO PRESTO® A40

# **Test Conditions**

JULABO unit Cooling power

Heating capacity
Band limit
Flow pressure
Bath fluid
Reactor
Jacket volume

Control

+20 °C 1.2 kW 0 °C 0.9 kW -20 °C 0.6 kW 2.7 kW No 0.40 bar JULABO Thermal HL40 Triple walled 20 liters glass reactor (Asahi) filled with 18 liter JULABO Thermal HL40 7.0 l External (ICC)


## **Test Results**

See chart on back page: The A40 cooling process from +170 °C to +100 °C in 1 h 5 min without overshoot.



#### Environment

Room temperature	+20 °C
Humidity	45 %
Voltage	230 V / 50 Hz



**Tip** You can also use the robust Pt100 with PTFE coating.


More tips on back page >>



JULABO GmbH Eisenbahnstraße 45 77960 Seelbach / Germany Tel. +49 (0) 7823 51-0



#### www.julabo.de



# Tip

Make use of the option to regulate the pump pressure. You can define the desired pressure in the PRESTO® settings.



## Tip

The Ethernet interface permits full access to all operational functions of the PRESTO<sup>®</sup>.



JULABO GmbH Eisenbahnstraße 45 77960 Seelbach / Germany Tel. +49 (0) 7823 51-0



## www.julabo.de