# Julabo Case Study

## **JULABO PRESTO® A80**

Heating a 20 liters reactor from +100 °C to +160 °C



## **Objective**

This case study tests the heating power of JULABO PRESTO® A80 with a 20 liters glass reactor. The A80 is connected to the reactor via two 2.0 m metal tubings. The A80 is programmed to heat up from +100 °C to +160 °C.

## **Environment**Room temperature

Room temperature +20 °C Humidity 45 %

Voltage 230 V / 50 Hz

## **Test Conditions**

JULABO unit

Cooling power

JULABO PRESTO® A80

+20 °C 1.2 kW

0 °C 1.2 kW

-20 °C 1.1 kW

Heating capacity 1.8 kW
Band limit No
Flow pressure 0.40 bar

Bath fluid JULABO Thermal HL80
Reactor 20 liters glass reactor (Asahi)

filled with 18 liter JULABO Thermal HL40

Jacket volume 7.0 l

Control External (ICC)


### **Test Results**

See chart on back page: The A80 heating process from +100 °C to +160 °C in 3 h 45 min without overshoot.




JULABO GmbH Eisenbahnstraße 45 77960 Seelbach / Germany Tel. +49 (0) 7823 51-0





SetpointTemperature in reactor's interiorTemperature in reactor's jacket





JULABO GmbH Eisenbahnstraße 45 77960 Seelbach / Germany Tel. +49 (0) 7823 51-0

