

PRESTO® A80t

Cool-down a 20 liters reactor from +20 °C to lowest possible temperature

Objective

This case study tests the lowest temperature of the PRESTO® A80t with a 20 liters glass reactor. The PRESTO® A80t is connected to the reactor via 1 m metal tubings. The PRESTO® A80t cools down from +20 °C to the lowest possible temperature.



Room temperature +20 °C Humidity 45 %

Voltage 208 V / 60 Hz



Test Conditions

JULABO unit PRESTO® A80t
Cooling power +20 °C 1.2 kW

0 °C 1.2 kW -20 °C 1.1 kW

Heating capacity 3.4 kW
Band limit with
Flow pressure 0.5 bar
Bath fluid Ethanol

Reactor 20 l glass reactor (Chemglass)

filled with 19 I Ethanol

Jacket volume 8 l

Control External (ICC)








Test Results

The PRESTO® A80t cooled the reactor from +20 °C down to the lowest possible temperature in 5 h. Within these test conditions the lowest possible temperature is -70 ° C.



Setpoint
Temperature in reactor's interior
Temperature in reactor's jacket

Tip Take advantage of our wide range of accessories. The M+R adapter enables you to display and record an additional temperature.



Tip Use the free of charge <i>EasyTEMP</i> software to control the units with the PC and to show the temperature curves graphically.	Edsy TEMP