

# PRESTO® A80t Cooling a 20 liters reactor from +20 °C to -60 °C

# Objective

0

This case study tests the cooling power of PRESTO<sup>®</sup> A80t with a 20 liters glass reactor. The PRESTO<sup>®</sup> A80t is connected to the reactor via two 1 m metal tubings. The PRESTO<sup>®</sup> A80t is programmed to cool down from +20 °C to -60 °C.

#### Environment

Room temperature	+20 °C
Humidity	45 %
Voltage	208 V / 60 Hz

# **Test Conditions**

JULABO unit Cooling power

Heating capacity Band limit Flow pressure Bath fluid Reactor

Jacket volume

Control

0 °C 1.2 kW -20 °C 1.1 kW 3.4 kW with 0.5 bar Thermal HL80 20 I glass reactor (Chemglass) filled with 19 I Ethanol 8 I External (ICC)

PRESTO® A80t

+20 °C 1.2 kW








#### **Test Results**

The PRESTO® A80t cooling process from +20 °C to -60 °C in 2 h 50 min without overshoot.



Measured with EasyTEMP Professional

Tip

E I

Ret C

You can also use the robust Pt100 with PTFE coating.



### Tip

Take advantage of our wide range of accessories. The M+R adapter enables you to display and record an additional temperature.



h