# **Operating Manual**

## SemiChill Recirculating Coolers

**Professional Series**  $(\mathcal{F})$ 9 SC 5000a air cooled Julabo SC 10000 SC 5000w water cooled . SC 10000w water cooled THE TEMPERATURE CONTROL COMPANY JULABO GmbH 77960 Seelbach / Germany Tel. +49 (0) 7823 / 51-0 Fax +49 (0) 7823 / 24 91 info.de@julabo.com

1.951.4576-V10 11/17

www.julabo.com

#### Congratulations!

You have made an excellent choice.

JULABO thanks you for the trust you have placed in us.

This operating manual has been designed to help you gain an understanding of the operation and possible applications of our circulators. For optimal utilization of all functions, we recommend that you thoroughly study this manual prior to beginning operation.

#### The JULABO Quality Management System



Temperature control devices for research and industry are developed, produced, and distributed according to the requirements of ISO 9001 and ISO 14001. Certificate Registration No. 01 100044846

#### Unpacking and inspecting

Unpack the circulator and accessories and inspect them for possible transport damage. Damage should be reported to the responsible carrier, railway, or postal authority, and a damage report should be requested. These instructions must be followed fully for us to guarantee our full support of your claim for protecting against loss from concealed damage. The form required for filing such a claim will be provided by the carrier.

Printed in Germany

Changes without prior notification reserved

Important: keep original operating manual for future use

### TABLE OF CONTENTS

Opera	ating manual	5
1. lı	ntended use	5
1.1.	Description	5
2. 0	Dperator responsibility – Safety instructions	6
2.1.	Disposal	
2.2.	EC Conformity	9
2.3.	Warranty conditions	11
2.4.	Technical specifications	12
2.5.	Cooling water connection	17
Opera	ating instructions	
3. C	Dperating controls and functional elements	
4. S	Safety notes for the user	
4.1.	Explanation of safety notes	
4.2.	Explanation of other notes	
4.3.	Safety instructions	
5. F	Preparations	
5.1.	Installation	
5.2.	Connecting the external system	
	5.2.1.Tubing	
	5.2.2.Bath fluids	
5.3.	<b>C</b>	
5.4.	5	
6. C	Operating procedures	
6.1.		
6.2.	Switching on / Start - Stop	
7. L	$\underline{T}$ / Setting the temperatures	
8.	Safety installations, warning functions	
8.1.	Excess temperature protection	
8	3.1.1.Early warning system, low level protection	
	Over and Sub temperature warning functions	
8 Г	3.2.1.Change-over of the warning function to shutdown function	
9.	MENU Menu functions	

9.1. MENU PUMP - Setting the pump pressure	35
9.2. MENU CONTROL – Control parameters	
9.2.1.CONTROL – internal / external control	
9.2.2.SELFTUNING	
9.2.3.DYN INT - Dynamic internal 9.2.4.Control parameters – XP, TN, TV internal	
9.2.5.COSPEED - external	
9.2.6.Control parameters – XPU, XP, TN, TV external	
9.3. MENU CONFIG - configuration	42
9.3.1.SETPOINT – Keypad control or remote control	43
9.3.2.A-START – Autostart	
9.3.3.OFF-MODE – Pump motor on / off	
9.3.4.RESET – Factory settings 9.3.5.ACTVAR - actuating variable	
9.3.6.TIME / DATE – setting time and date	
9.4. MENU SERIAL - BAUDRATE, HANDSHAKE, PARITY	
9.5. MENU LIMITS	49
9.6. MENU PROGRAM – Integrated programmer	51
9.7. MENU ADJUST – ATC Absolute Temperature Calibration	55
9.8. MENU ANALOG – Analog inputs/outputs	59
10. Troubleshooting guide / Error messages	65
11. Electrical connections	67
12. Remote control	69
12.1. Setup for remote control	69
12.2. Communication with a PC or a superordinated data system	69
12.3. List of commands	70
12.4. Status messages	74
12.5. Error messages	74
13. JULABO Service – Online remote diagnosis	75
14. Cleaning / repairing the unit	76

## **Operating manual**

## 1. Intended use

JULABO recirculating coolers have been designed to control the temperature of specific fluids in a bath tank. The units feature pump connections for temperature control of external systems (loop circuit).



JULABO recirculating coolers are not suitable for direct temperature control of foods, semi-luxury foods and tobacco, or pharmaceutical and medical products. Direct temperature control means unprotected contact of the object with the bath medium (bath fluid).

implemented microprocessor technology allows to set and to store different values that can be indicated on the VFD COMFORT-DISPLAY. Three menu keys facilitate

## 1.1. Description





☑ The integrated programmer allows storing and running temperature and timedependent processes.

"ICC - Intelligent Cascade Control" represents the supreme solution temperature control. ICC overs perfect temperature control with self-optimizing PID control

✓ The TCF - Temperature Control Features allow the user to have access to all

☑ The recirculating coolers are operated via the splash-proof keypad. The

adjusting setpoints, warning and safety functions and menu functions.





parameters.





PUMP

Absolute Temperature Calibration (ATC3) provides a high temperature stability in the bath. With the 3-point calibration an offset is adjusted at three temperatures to ensure an accurate temperature pattern at the selected spot in the bath over the full temperature range.

important temperature control parameters. This means: Full control on the control

mode and the chance to manually adjust or adapt control to the specific application.

- Electrical connections: The serial interface, switchable from RS232 to RS485, allows modern process **RS232** technology without additional interface. Connection for Pt100 external sensor for external temperature measurement and control. Pt100 The electronic module (option) provides 3 further analog connections (alarm input, standby input, recorder output, programmer input). ☑ The excess temperature protection conforming to IEC 61010-2-010 is a safety installation independent from the control circuit. This protection can be indicated and set on the VFD COMFORT-DISPLAY. ☑ The early warning system for low level signals that bath fluid needs to be refilled. before the low level protection conforming to IEC 61010-2-010 causes a complete shutdown of the main functional elements. SMART
  - ☑ Intelligent pump system: The pump capacity (electronically adjustable via the motor speed) enables to adapt to varying conditions for internal and external temperature applications.

## 2. Operator responsibility – Safety instructions

The products of JULABO ensure safe operation when installed, operated, and maintained according to common safety regulations. This section explains the potential dangers that may arise when operating the circulator and also specifies the most important safety precautions to preclude these dangers as far as possible.

The operator is responsible for the qualification of the personnel operating the units.

- The personnel operating the units should be regularly instructed about the dangers involved with their job activities as well as measures to avert these dangers.
- Make sure all persons tasked with operating, installing, and maintaining the unit have read and understand the safety information and operating instructions.
- When using hazardous materials or materials that could become hazardous, the circulator may be operated only by persons who are absolutely familiar with these materials and the circulator. These persons must be fully aware of possible risks.

If you have any questions concerning the operation of your unit or the information in this manual, please contact us!

Contact:	JULABO GmbH	Tel. +49 (0) 7823 / 51-0	info.de@julabo.com
	Gerhard-Juchheim-Strasse 1 77960 Seelbach / Germany	Fax +49 (0) 7823 / 24 91	www.julabo.com

### Safety instructions for the operator:

- Avoid strikes to the housing, vibrations, damage to the operating-element panel (keypad, display), and contamination.
- Make sure the product is checked for proper condition regularly (depending on the conditions of use). Regularly check (at least every 2 years) the proper condition of the mandatory, warning, prohibition and safety labels.
- Make sure that the mains power supply has low impedance to avoid any negative effects on the instruments being operated on the same mains.
- This unit is designed for operation in a controlled electromagnetic environment. This means that transmitting devices (e.g., cellular phones) should not be used in the immediate vicinity.
- Magnetic radiation may affect other devices with components sensitive to magnetic fields (e.g., monitors). We recommend maintaining a minimum distance of 1 m.
- > Permissible ambient temperature: max. 40 °C, min. 5 °C.
- Permissible relative humidity: 50% (40 °C).
- > Do not store the unit in an aggressive atmosphere. Protect the unit from contamination.
- > Do not expose the unit to sunlight.

#### Appropriate operation

Only qualified personnel is authorized to configure, install, maintain, or repair the circulator. Persons who operate the circulator must be trained in the particular tasks by qualified personnel. The summarized user guidance (short manual) and the specification table with information on individual parameters are sufficient for this.

#### Use

The bath can be filled with flammable materials. Fire hazard!

There might be chemical dangers depending on the bath medium used.

Observe all warnings for the used materials (bath fluids) and the respective instructions (safety data sheets).

Insufficient ventilation may result in the formation of explosive mixtures. Only use the unit in well ventilated areas.

Only use recommended materials (bath fluids). Only use non-acid and non corroding materials.

When using hazardous materials or materials that could become hazardous, **the operator must** affix the enclosed safety labels (1 + 2) to the front of the unit so they are highly visible:

	•	
1		Danger area. Attention! Observe instructions. (operating manual, safety data sheet)
2 or		Carefully read the user information prior to beginning operation. Scope: EU
2		Carefully read the user information prior to beginning operation. Scope: USA, NAFTA

Particular care and attention is necessary because of the wide operating range. There are thermal dangers: Burn, scald, hot steam, hot parts and surfaces that can be touched.



Hot surface warning. (The label is put on by JULABO)

Observe the instructions in the manuals for instruments of a different make that you connect to the circulator, particularly the corresponding safety instructions. Also observe the pin assignment of plugs and technical specifications of the products.

### 2.1. Disposal

The recirculating cooler contains a back-up battery that supplies voltage to the memory chips when the unit is switched off. Do not dispose of the battery with household waste!

Depending on battery regulations in your country, you may be obligated to return used or defective batteries to collection sites.

The product may be used with oil as bath fluid. These oils fully or partially consist of mineral oil or synthetic oil. For disposal, follow the instructions in the material safety data sheets.

This unit contains the refrigerant R404A or R452A, which at this time is not considered harmful to the ozone layer. However, over the long operating period of the unit, disposal rules may change. Therefore, only qualified personnel should handle the disposal.

#### Valid in EU countries



See the current official journal of the European Union – WEEE directive. Directive of the European Parliament and of the Council on waste electrical and electronic equipment (WEEE).

This directive requires electrical and electronic equipment marked with a crossedout trash can to be disposed of separately in an environmentally friendly manner. Contact an authorized waste management company in your country.

Disposal with household waste (unsorted waste) or similar collections of municipal waste is not permitted!

## 2.2. EC Conformity

EG-Konformitätserklärung nach EG Maschinenrichtlinie 2006/42/EG, Anhang II A EC-Declaration of Conformity to EC Machinery Directive 2006/42/EC, Annex II A

Hersteller / Manufacturer:		JULABO GmbH Gerhard-Juchheim-Strasse 1 77960 Seelbach / Germany Tel: +49(0)7823 / 51 - 0		
Hiermit erklären wir, We hereby declare, tha	dass das nachfolgend bezeichnet t the following product	e Produkt		
Produkt / Product:	Umlaufkühler / Recirculating Cod	oler		
Тур / <i>тур</i> е:	SC5000a; SC5000w	Serien-Nr. / Serial-No.: siehe Typenschild / see type label		
Sicherheits- und Ges due to the design and c	sundheitsanforderungen den nach	ns in Verkehr gebrachten Ausführung den grundlegenden folgend aufgeführten EG-Richtlinien entspricht. ted by our Company – complies with fundamental safety and health		
EMV-Richtlinie 20	nie 2006/42/EG; Machinery Direc 014/30/EU; EMC-Directive 2014/30 2011/65/EU; RoHS-Directive 2011	D/EU		
	nonisierte Normen und tech nduct is in compliance with the follo	n. Spezifikationen: wing harmonized standards and technical specifications:		
EN 50581 : 2012 Technische Dokumentation Technical documentation fr	) zur Beurteilung von Elektro- und Elektronikgeräte or the assessment of electrical and electronic prod	en hinsichtlich der Beschränkung gefährlicher Stoffe ucts with respect to the restriction of hazardous substances		
EN ISO 12100 : 2010 Sicherheit von Maschinen - Allgemeine Gestaltungsleitsätze - Risikobeurteilung und Risikominderung (ISO 12100:2010) Safety of machinery - General principles for design - Risk assessment and risk reduction (ISO 12100:2010)				
EN 61010-1 : 201 Sicherheitsbestimmungen Safety requirements for ele	0 für elektrische Mess-, Steuer-, Regel- und Laborge actrical equiment for measurement, control, and lat	eräte, Teil 1: Allgemeine Anforderungen boratory use, Part 1: General requirements		
EN 61010-2-010 : Sicherheitsbestimmungen Safety requirements for ele materials	für elektrische Mess-, Steuer-, Regel- und Laborge	eräte, Teil 2-010: Besondere Anforderungen an Laborgeräte für das Erhitzen von Stoffen boratory use, Part 2-010: Particular requirements for laboratory equipment for the heating of		
EN 61326-1 : 201 Elektrische Mess-, Steuer- Electrical equipment for m	3 , Regel- und Laborgeräte- EMV-Anforderungen- Te easurement, control, and laboratory use - EMC req	eil 1: Allgemeine Anforderungen juirements - Part 1: General requirements		
Auswahlkriterien		nte Anforderungen – Teil 1: Grundlegende Anforderungen, Begriffe, Klassifikationen und		
EN 378-2 : 2016		ents - Part 1: Basics requirements, definitions, classification and selection criteria nte Anforderungen – Teil 2: Konstruktion, Herstellung, Prüfung, Kennzeichnung und		
Refrigerating systems and EN 378-3 : 2016 Kälteanlagen und Wärmen	oumpen – Sicherheitstechnische und umweltreleva	ents - Part 2: Design, construction, testing, marking and documentation nte Anforderungen – Teil 3: Aufstellungsort und Schutz von Personen		
EN 378-4 : 2016 Kälteanlagen und Wärmer	oumpen – Sicherheitstechnische und umweltreleva	ents - Part 3: Installation site and personal protection nte Anforderungen – Teil 4: Betrieb, Instandhaltung, Instandsetzung und Rückgewinnung ents - Part 4: Operation, maintenance, repair and recovery		
Authorized represent	r für die Zusammenstellung ( ative in charge of administering tec ke, im Hause / on the manufacturer's	chnical documentation:		
	erklärung wurde ausgestellt			

2017\_080\_SC5000(a)(w)-Umlaufkühler\_d\_e.docx

#### EG-Konformitätserklärung nach EG Maschinenrichtlinie 2006/42/EG, Anhang II A EC-Declaration of Conformity to EC Machinery Directive 2006/42/EC, Annex II A

Hersteller /	Manufacturer:
--------------	---------------

JULABO GmbH Gerhard-Juchheim-Strasse 1 77960 Seelbach / Germany Tel: +49(0)7823 / 51 - 0

Hiermit erklären wir, dass das nachfolgend bezeichnete Produkt We hereby declare, that the following product

Produkt / Product: Umlaufkühler / Recirculating Cooler

Typ / Type: SC10000w Serien-Nr. / Serial-No.: siehe Typenschild / see type label

aufgrund seiner Konzipierung und Bauart in der von uns in Verkehr gebrachten Ausführung den grundlegenden Sicherheits- und Gesundheitsanforderungen den nachfolgend aufgeführten EG-Richtlinien entspricht. due to the design and construction, as assembled and marketed by our Company - complies with fundamental safety and health requirements according to the following EC-Directives.

#### Maschinenrichtlinie 2006/42/EG; Machinery Directive 2006/42/EC EMV-Richtlinie 2014/30/EU; EMC-Directive 2014/30/EU

RoHS-Richtlinie 2011/65/EU; RoHS-Directive 2011/65/EU

#### Angewandte harmonisierte Normen und techn. Spezifikationen: The above-named product is in compliance with the following harmonized standards and technical specifications:

#### EN 50581 : 2012

Technische Dokumentation zur Beurteilung von Elektro- und Elektronikgeräten hinsichtlich der Beschränkung gefährlicher Stoffe Technical documentation for the assessment of electrical and electronic products with respect to the restriction of hazardous substances

#### EN ISO 12100 : 2010

Sicherheit von Maschinen - Allgemeine Gestaltungsleitsätze - Risikobeurteilung und Risikominderung (ISO 12100:2010) Safety of machinery - General principles for design - Risk assessment and risk reduction (ISO 12100:2010)

#### EN 61010-1 : 2010

Sicherheitsbestimmungen für elektrische Mess-, Steuer-, Regel- und Laborgeräte, Teil 1: Allgemeine Anforderungen Safety requirements for electrical equiment for measurement, control, and laboratory use, Part 1: General requirements

#### EN 61010-2-010 : 2014

Sicherheitsbestimmungen für elektrische Mess-, Steuer-, Regel- und Laborgeräte, Teil 2-010: Besondere Anforderungen an Laborgeräte für das Erhitzen von Stoffen Safety requirements for eletrical equipment for measurement, control, and laboratory use, Part 2-010: Particular requirements for laboratory equipment for the heating of

#### EN 61326-1 : 2013

Elektrische Mess-, Steuer-, Regel- und Laborgeräte- EMV-Anforderungen- Teil 1: Allgemeine Anforderungen Electrical equipment for measurement, control, and laboratory use - EMC requirements - Part 1: General requirements

#### EN 378-1:2016

Kälteanlagen und Wärmepumpen – Sicherheitstechnische und umweltrelevante Anforderungen – Teil 1: Grundlegende Anforderungen, Begriffe, Klassifikationen und Auswahlkriterien Refrigerating systems and heat pumps - Safety and environmental requirements - Part 1: Basics requirements, definitions, classification and selection criteria

#### FN 378-2:2016

Kälteanlagen und Wärmepumpen – Sicherheitstechnische und umweltrelevante Anforderungen – Teil 2: Konstruktion, Herstellung, Prüfung, Kennzeichnung und

Dokumentation Refrigerating systems and heat pumps - Safety and environmental requirements - Part 2: Design, construction, testing, marking and documentation

#### EN 378-3 : 2016

Kälteanlagen und Wärmepumpen – Sicherheitstechnische und umweltrelevante Anforderungen – Teil 3: Aufstellungsort und Schutz von Personen Refrigerating systems and heat pumps - Safety and environmental requirements - Part 3: Installation site and personal protection

#### EN 378-4 : 2016

Kälteanlagen und Wärnepumpen – Sicherheitstechnische und umweltrelevante Anforderungen – Teil 4: Betrieb, Instandhaltung, Instandsetzung und Rückgewinnung Refrigerating systems and heat pumps - Safety and environmental requirements - Part 4: Operation, maintenance, repair and recovery

Bevollmächtigter für die Zusammenstellung der techn. Unterlagen: Authorized representative in charge of administering technical documentation:

Hr. Torsten Kauschke, im Hause / on the manufacturer's premises as defined above

Die Konformitätserklärung wurde ausgestellt The declaration of conformity was issued and valid of

Seelbach, 20.10.2017

M. Juchheim, Geschäftsführer / Managing Director

2017\_081\_SC10000(w)-Umlaufkühler\_d\_e.docx

### 2.3. Warranty conditions

JULABO GmbH warrants its products against defects in material or in workmanship, when used under appropriate conditions and in accordance with appropriate operating instructions

### for a period of ONE YEAR.

Extension of the warranty period - free of charge



With the '1PLUS warranty' the user receives a free of charge extension to the warranty of up to 24 months, limited to a maximum of 10 000 working hours.

To apply for this extended warranty the user must register the unit on the JULABO web site <u>www.julabo.com</u>, indicating the serial no. The extended warranty will apply from the date of JULABO GmbH's original invoice.

JULABO GmbH reserves the right to decide the validity of any warranty claim. In case of faults arising either due to faulty materials or workmanship, parts will be repaired or replaced free of charge, or a new replacement unit will be supplied.

Any other compensation claims are excluded from this guarantee.

### 2.4. Technical specifications

Professional Serie	25			SC 5000a	SC 5000w
Working temperatu				00 00000	00 0000
Standard unit	no rangee.		°C	+5 +35	+5 +35
Option Low Te	mn		°C	-20 +35	-20 +35
Option Low /Hi	•		°C		-20 +80
Option Low /Hi	•		°C	-20 +130	-20 +130
Temperature stabil	•		°C	0.1	0.1
Absolute Tempera	•	INT/E			±3 / ±9
Heater wattage Op	tion H5	ł	w	5.0	5.0
Heater wattage Op		ŀ	w	12.0	12.0
Cooling capacity			°C	<u>+20 0 -10</u>	<u>+20 0 -10</u>
Medium ethanol		ł	kW	5.0 2.5 1.2	5.0 2.5 1.2
Refrigerant				R404A, R452A*	R404A, R452A*
Pump capacity P3	(Standard) / P4 *	* (Option)			
see table 1 page 1	4				
Flow rate	P3 / P4	l/min at 0 k	bar	33 / 43	33 / 43
Pressure max.	P3 / P4	bar at 0 lite	ers	3.5 / 4.4	3.5 / 4.4
Overall dimensions	s (WxDxH)	(	cm	59x67/112	59x67/112
Filling volume		lite	ers	43 60	43 60
Weight			kg	153	153
Ambient temperatu	ire		°C	5 40	5 40
Mains power conne			Hz	400/50/3 Phases	400/50/3 Phases
	) V/3PNPE/50 Hz				
Current consumption			Α	7 (at 400 V)	7 (at 400 V)
Current consumption			Α	15/Phase	15/Phase
Current consumption	on 12 kW heater	/ P3 / P4	Α	30/Phase	30/Phase
Mains power conne		V/	Hz	208-230/60/3	208-230/60/3 Phases
	) V/3PPE/60 Hz			Phases	
Current consumpt.	•		A	9 (208 V) / 8 (230V)	
Current consumpt.			A	12 (208 V) /13(230V)	12 (208 V)/ 13 (230V)
Current consumpt.	-		A	23 (208 V)/22(230 V)	23 (208 V)/ 22 (230 V)
Current consumpt.			A	25 (208 V)/26(230 V)	
Current consumption		•		37 (208 V)/39(230 V)	37 (208 V)/39(230 V)
Current consumpti	on 12 kW heater	/phase / P4	A	41 (208 V)/43(230 V)	41 (208 V)/43(230 V)

All measurements have been carried out at: rated voltage and frequency ambient temperature: 20 °C Technical changes without prior notification reserved.

\* at 400 V / 50 Hz

\*\* Pump P4- reduces cooling capacity by 0.3 kW

Professional Seri	es		SC 10000w
Working temperati			
Standard unit		°C	+5 +35
Option Low Te	mp	°C	-20 +35
Option Low /H	•	°C	-20 +80
Option Low /H	•	°C	-20 +130
Temperature stabi	lity	°C	0.1
Absolute Tempera	ture <b>C</b> alibration	INT/EXT	±3 / ±9
Heater wattage Op	otion H5	kW	5.0
Heater wattage Op	otion H12	kW	12.0
Cooling capacity		°C	<u>+20 0 -10</u>
Medium ethanol		kW	10 5.0 2.5
Refrigerant			R404A, R452A*
Pump capacity P3	(Standard) / P4 ** (Option)		
see table 1 page 1			
	P3 / P4	lpm at 0 bar	
Pressure max.	P3 / P4	bar at 0 liters	3.5 / 4.4
-			
Overall dimension	s (WxDxH)	cm	59x67/112
Filling volume		liters	43 60
Weight		kg	155
Ambient temperati	ure	°C	5 40
<b></b>		· · · · · ·	
•	ection 365 V-440 V/3PNPE/50 H		400/50/3 phases
•	ion without heater /P3 /P4	A	11 (400 V)
•	ion (at 400 V) 5 kW heater /P3 /		17 (400 V)
Current consumpt	ion (at 400 V) 12 kW heater /P3 /	'P4 A	31 (400 V)
N.4. 1		\ <i>//</i> + +	
•	ection 197 V-254 V/3PPE/60 Hz		208-230/60/3 phases
	ion without heater / P3	A	15 (208 V) / 16 (230 V)
	ion without heater / P4	A	17 (208 V) / 17 (230 V)
	ion with 5 kW heater / P3	A	25 (208 V) / 26 (230 V)
Current consumpti	ion with 12 kW heater / P4	A	44 (208) / 46 (230 V)

All measurements have been carried out at: rated voltage and frequency ambient temperature: 20 °C Technical changes without prior notification reserved.

\* at 400 V / 50 Hz

\*\* Pump P4- reduces cooling capacity by 0.3 kW

Professional series				
Temperature selection		digital		
via keypad		indication on VFD COMFORT-DISPLAY		
remote control via personal computer		indication on monitor		
Temperature indication	°C	VFD COMFORT-DISPLAY		
Resolution	°C	0.1		
Temperature control		ICC - Intelligent Cascade Control		

Computer interface RS232

External	Pt100	sensor
LAGINA	1 1100	301301

<u>Optional</u>	
Programmer input	-100 °C to 400 °C = 0 - 10 V or 0 - 20 mA or 4 - 20 mA
Input for the signal of a flow meter	er or external manipulated variable
Temperature recorder outputs	0 - 10 V (0 V = -100 °C, 10 V = 400 °C)
	0 - 20 mA (0 mA = -100 °C, 20 mA = 400 °C)
	4 - 20 mA (4 mA = -100 °C, 20 mA = 400 °C)
Standby input	for external emergency switch-off
Alarm output	for external alarm signal

### Table 1

Pump capacity Bath fluid: Water, Silicone -oil		XX.X PUMP	<b>X</b> 1		KX.X PUMP	<b>X</b> 2		XX.XX PUMP 3
Circulating pump:	P0	P3		P0	P3		Р0	
Flow rate max. Lpm at 0 bar	31	30		42	33		48	
Pressure max. bar at 0 liter	0.75	1.8		1.2	3.5		1.78	
Bath fluid: Galden <sup>®</sup> eg. Fluorinert <sup>®</sup> 3283					•			
Pressure max. bar at 0 liter		3.5		1,78				



Notice: If Galden<sup>®</sup> or Fluorinert<sup>®</sup> is used the the charge of the motor increases. Wrong adjustment causes overheating and eventually destruction of the motor. With the circulation pump P0 a maximum pump pressure stage >PUMP 2< may be adjusted. With the circulation pump P3 a maximum pump pressure stage >PUMP 1< may be adjusted. The pump type can be recognized in the order no. on the name plate. 95 xx xxx \_ xx **PX** xx

Safety installations according to IEC 61010-2-010: Excess temperature protection Low liquid level protection Classification according to DIN 12876-1

adjustable from 0 °C ... 220 °C float switch class III

optical + audible (in intervals)

Supplementary safety installations	
Early warning system for low level	float switch
High temperature warning function	optical + audible (in intervals)
Low temperature warning function	optical + audible (in intervals)
Supervision of working sensor	plausibility control
Reciprocal sensor monitoring between	
working and safety sensors	difference >35 K
Alarm message	optical + audible (permanent)

Alarm message Warning message

Environmental conditions according to IEC 61 010-1:

	Use indoors only.		
	Altitude up to 2000 m - normal zero.		
	Ambient temperature: see Technical specifications		
	Humidity:		
	Max. relative humidity 80% for temperatures up to +31 °C,		
	linear decrease down to 50% relative humidity at a temperature of +40 °C		
	Max. mains voltage fluctuations of ±1	0% are permissible.	
Protection class according to IEC 60 529 The unit corresponds to Class I		IP21	
Overvoltage category		II	
Pollution degree		2	



Caution:

The unit is not for use in explosive environment.

#### **EMC** requirements

The device is an ISM device of group 1 per CISPR 11 (uses HF for internal purposes) and is classified in class A (industrial and commercial sector).

### Notice:

- Devices of class A are intended for the use in an industrial electromagnetic environment.
- When operating in other electromagnetic environments, their • electromagnetic compatibility may be impacted.

#### Information about the used refrigerants

The **Regulation (EU) No. 517/2014 on fluorinated greenhouse gases** applies to all systems which contain fluorinated refrigerants and replaces (EC) 842/2006.

The aim of the Regulation is to protect the environment by reducing emissions of fluorinated greenhouse gases.

Among other things it regulates the emission limits, use and recovery of these substances. It also contains requirements for operators of systems which require / contain these substances to function.

Under Regulation 517/2014, the operator of a system of this nature has the following duties:

- The operator must ensure that the equipment is checked at regular intervals for leaks.
- These intervals depend on the CO<sub>2</sub> equivalent of the system. This is calculated from the refrigerant fill volume and type of refrigerant. The CO<sub>2</sub> equivalent of your system is shown on the model plate.
- The operator undertakes to have maintenance, repair, service, recovery and recycling work carried out by certified personnel who have been authorized by JULABO.
- All such work must be documented. The operator must keep records and archive them for at least five years. The records must be submitted to the relevant authority on request.

Refer to the text of the Regulation for further information.

#### 2.5. **Cooling water connection**

Cooling water pressure (IN/OUT)	max.	6 bar
Pressure difference (IN - OUT)		3.5 to 6 bar
Cooling water consumption (IN with 15 °C),	SC 5000w	10 l/min
	SC 10000w	26 l/min
Cooling water temperature		< 20 °C



#### Notice: Cooling water circuit

Risk of oil leaking from the refrigeration system (compressor) of the recirculating cooler into the cooling water in case of a fault in the cooling water circuit! Observe the laws and regulations of the water distribution company valid in the location where the unit is operated.



#### Notice:

Danger of corrosion of heat exchanger due to unsuitable quality of cooling water.


- Due to its high content of lime, hard water is not suitable for cooling and causes • scale in the heat exchanger.
- Ferrous water or water containing ferrous particles will cause formation of rust even in heat exchangers made of stainless steel.
- Chlorinated water will cause pitting corrosion in heat exchangers made of • stainless steel.
- Due to their corrosive characteristics, distilled water and deionized water are • unsuitable and will cause corrosion of the bath.
- Due to its corrosive characteristics, sea water is not suitable.
- Due to its microbiological (bacterial) components, which settle in the heat exchanger, untreated and unpurified river water and water from cooling towers is unsuitable.
- Avoid particulate matter in cooling water.
- Avoid putrid water. •

#### Recommended quality of cooling water:

рН	7.5 to 9.0
Sulfate [SO4 2- ]	< 100 ppm
Hydrocarbonate [HCO 3-]/sulfate [SO4 2-]	> 1 ppm
Hardness [Ca 2+, Mg 2+]/[HCO 3-]	> 0.5 °dH
Alkalinity	60 ppm < [HCO 3-] < 300 ppm
Conductivity	< 500 µS/cm
Chloride (Cl -)	< 50 ppm
Phosphate (PO4 3-)	< 2 ppm
Ammonia (NH3)	< 0.5 ppm
Free chlorine	< 0.5 ppm
Trivalent iron ions (Fe 3+)	< 0.5 ppm
Manganese ions (Mn 2+)	< 0.05 ppm
Carbon dioxide (CO2)	< 10 ppm
Hydrogen sulfide (H2S)	< 50 ppm
Content of oxygen	< 0.1 ppm
Algae growth	impermissible
Suspended solids	impermissible

## **Operating instructions**

## 3. Operating controls and functional elements



	ESC	Escape key 1. Cancel entries 2. Return to a higher menu level
5		VFD COMFORT-DISPLAY
	S 15000 E 145.30	<ul> <li>Header: Control indicators see sections 2.1 and 2.2</li> <li>Line 1: Actual value internal or external The display is depending on the selected control mode in the menu &gt; Control &lt; (internal or external).</li> <li>Line 2: Working temp. setpoint, constantly S xxx.xx</li> <li>Line 3: Actual value (E = external or I = internal) Alternating with the display in line 1</li> </ul>
	XXXXX S xxxxx FL GOOD XXXXX S xxxxx R 005	Use the keys to indicate further values in line 3. However, the functions of these keys are different with the programmer started. PI Capacity in % - with manipulated variable set to >control<* or PS Capacity in % - with manipulated variable set to >SERIAL<* or >EPROG<* F Flow rate in liters/minute (providing EPROG input set to >FLOWRATE<) *see 9.3.5. ACTVAR - page 46 FL Status indication for flow >GOOD< = Pump switched on R Resistifity measurement and actual value display in the range from 0.5 5 MOhm/cm Press ESC to return to actual value (E = external or I = internal)
5.1	<u>₩</u>	Control indicators in the header: Heating / Cooling / Alarm / <b>R</b> emote control
5.2	°C Int °F Ext	<u>Control indicators in the header:</u> Temperature indication <b>Int</b> ernal or <b>Ext</b> ernal actual value Temperature indication in °C (°F not possible on this unit)
5.3	*	Display for the adjusted pump pressure stage. Adjustable via the MENU button, in the menu >PUMP<.
6		Filling level indication
7		Manometer (feed pressure)
8	Selo-	Handle
9		Filling opening
10		Venting grid, removable (only air cooled recirculating cooler)

Optio	on: Electronic mod	dule
13	ALARM	Alarm output (for external alarm signal)
14	STAND-BY	Standby input (for external emergency switch-off)
15	REG+E-PROG	Programmer input and temperature recorder output
16	ext Pt100	Socket for external measurement and control sensor or external setpoint programming
19		2 Connectors for solenoid valve. 230 V / max. 0.1 A No control voltage in the <b>-OFF-</b> condition
20		4 Safety cutouts: Mains fuses 16 A (with option H5)
		2 Safety cutouts: Mains fuses 10 A
		Motor protection circuit breaker for compressor motor
21		Mains power cable with plug
22	OUT IN	
23	A T	Pump connectors: 3/4" NPT male
	00	OUT / Feed IN / Return
24		Overflow connector, M10x1 female
		Order-No. 8 970 460 Barbed fitting for tubing 8 mm inner dia.
		$\textcircled{I}$ Closable when using e.g. 3M Fluorinert $^{ extsf{ iny B}}$ as temperature liquid.
25		Discharge nozzle with cap nut, Connection: $\frac{1}{2}$ " male
	Ţ	<b>Recommendation:</b> Before filling please install a drain cock at the discharge nozzle. (not included in delivery)
		Order-No. 8 920 100 Drain cock, stainless steel
26		Only for water cooled models: Cooling water OUTLET and INLET

## 4. Safety notes for the user

### 4.1. Explanation of safety notes

In addition to the safety warnings listed above, warnings are posted throughout the manual. These warnings are designated by an exclamation mark inside an equilateral triangle. "Warning of a dangerous situation (Attention! Please follow the documentation)." The danger is classified using a signal word. Read and follow these important instructions.
Warning: Describes a possibly highly dangerous situation. If these instructions are not followed, serious injury and danger to life could result.



#### Caution:

Describes a possibly dangerous situation. If this is not avoided, slight or minor injuries could result. A warning of possible property damage may also be contained in the text.



#### Notice:

Describes a possibly harmful situation. If this is not avoided, the product or anything in its surroundings can be damaged.

## 4.2. Explanation of other notes



Draws attention to something special.

### Important!

Note!

Indicates usage tips and other useful information.

## 4.3. Safety instructions

Follow the safety instructions to avoid personal injury and property damage. Also, the valid safety instructions for workplaces must be followed.



- Only connect the unit to a power socket with an earthing contact (PE protective earth)!
- The power supply plug serves as a safe disconnecting device from the line and must always be easily accessible.
- Place the unit on an even surface on a base made of nonflammable material.
- Do not stay in the area below the unit.
- Make sure you read and understand all instructions and safety precautions listed in this manual before installing or operating your unit.
- Set the excess temperature safety installation at least 25 °C below the fire point of the bath fluid.
- Never operate the unit without bath fluid in the bath.
- Pay attention to the thermal expansion of bath oil during heating to avoid overflowing of the fluid.
- Prevent water from entering the hot bath oil.

afety notes fo	or the user
	• Do not drain the bath fluid while it is hot! Check the temperature of the bath fluid prior to draining (e.g., by switching the unit on for a short moment).
	Use suitable connecting tubing.
	<ul> <li>Avoid sharp bends in the tubing, and maintain a sufficient distance from surrounding walls.</li> </ul>
	Make sure that the tubing is securely attached.
	Regularly check the tubing for material defects (e.g., for cracks).
	Never operate damaged or leaking units.
	<ul> <li>Always turn off the unit and disconnect the mains cable from the power source before performing any service or maintenance procedures, or before moving the unit.</li> </ul>
	<ul> <li>Always turn off the unit and disconnect the mains cable from the power source before cleaning the unit.</li> </ul>
	Always empty the bath before moving the unit.
	Transport the unit with care.
	Sudden jolts or drops may cause damage in the interior of the unit.
	Observe all warning labels.
	Never remove warning labels.
	Never operate units with damaged mains power cables.
	Repairs are to be carried out only by qualified service personnel.
	<ul> <li>Some parts of the bath tank and the pump connections may become extremely hot during continuous operation. Therefore, exercise particular caution when touching these parts.</li> </ul>
<b>A</b>	Caution:
	The unit may be used, for example, to control the temperature of fluids in a reactor.
	We do not know what substances are contained in these vessels. Many substances are:
	inflammable, easily ignited, or explosive
	hazardous to health
	environmentally hazardous
	i.e.: dangerous
	The user alone is responsible for the handling of these substances!
	The following questions should help to recognize possible dangers and to reduce the risks to a minimum
	<ul> <li>risks to a minimum.</li> <li>Are all tubes and electrical cables connected and laved?</li> </ul>

Are all tubes and electrical cables connected and layed? Note:

sharp edges, hot surfaces in operation, moving machine parts, etc.

- Do dangerous vapors or gases develop during heating? Must the work be done in a fume hood?
- What to do when a dangerous substance was spilled on or in the unit? • Before starting to work, obtain information concerning the substance and determine the method of decontamination.



#### Notice:

Check the safety installations at least twice a year!

- Excess temperature protection according to IEC 61010-2-010 With a screwdriver, turn back the adjustable excess temperature protection until the shutdown point (actual temperature).
- Low level protection according to IEC 61010-2-010 To check the function of the float, it can be manually lowered with a screwdriver, for example.


## 5. Preparations

### 5.1. Installation

- Place the unit in an upright position. For better stability, apply the holding brakes on the front casters.
- The place of installation should be large enough and provide sufficient air ventilation to ensure the room does not warm up excessively because of the heat the instrument rejects to the environment. (Max. permissible ambient temperature: 35 °C). For a fault (leakage) in the refrigeration system, the standard EN 378 prescribes a certain room space to be available for each kg of refrigerant.

The refrigerant quantity is specified on the type plate.

- > For 0.52 kg of refrigerant R404A, 1  $m^3$  of space is required.
- > For 0.423 kg of refrigerant R452A, 1  $m^3$  of space is required.
- Keep at least 20 cm of open space on the front and rear venting grids.
- Do not set up the unit in the immediate vicinity of heat sources and do not expose to sun light
- Before operating the unit after transport, <u>wait about one hour after</u> <u>installation</u>. This will allow any oil that has accumulated laterally during transport to flow back down, thus ensuring that the compressor can develop its maximum capacity.



#### **Recommendation:**

Before filling please install a drain cock at the discharge nozzle. (25) Connection:  $\frac{1}{2}$  male (not included in delivery)


Order-No. 8 920 100 Drain cock, stainless steel

#### **Cooling water connection (26)**



Only for water cooled models:				
Ensure circulation of cooling water by connecting the tubing to cooling				
water inlet and outlet on the rear of the recirculating cooler.				
Cooling water connectors	G3/4" male			
Cooling water temperature	< 20 °C			
Cooling water see page 17.				

### 5.2. Connecting the external system





Caution: Securely attach all tubing to prevent slipping.

### 5.2.1. Tubing

#### **Recommended tubing:**

	Maximum pressure
Textile-reinforced tubing	> 4.5 bar



#### Warning: Tubing:

At high working temperatures, the tubing used for temperature control and for the cooling water supply represents a danger source.

A damaged tubing line may allow a large amount of hot bath fluid to be pumped out within a short time.

#### This may result in:

- Burning of skin
- Breathing difficulties due to hot atmosphere

#### Safety instructions

- Use suitable connecting tubing.
- Make sure that the tubing is securely attached.
- Avoid sharp bends in the tubing and maintain a sufficient distance from surrounding walls.
- Regularly check the tubing for material defects (e.g., for cracks), at least once a year.
- Preventive maintenance: replace the tubing from time to time.

#### 5.2.2. Bath fluids



### Caution:

Carefully read the material safety data sheet of the bath fluid used, particularly with regard to the fire point!

If a bath fluid with a fire point of  $\leq$  65  $^\circ C$  is used, only supervised operation is possible.

#### Water:

The quality of water depends on local conditions.

- Due to the high concentration of lime, hard water is not suitable for temperature control because it leads to scale in the bath
- Ferrous water can cause corrosion, even on stainless steel.
- Chlorinated water can cause pitting corrosion.
- Distilled water and deionized water are unsuitable. Their special properties cause corrosion in the bath, even on stainless steel.

#### **Recommended bath fluids:**

Bath fluid	Temperature range
soft/decalcified water	5 °C to 80 °C



See website for list of recommended bath fluids. **Contact:** see page 6

	Caution:
<u>.</u>	<ul> <li>Fire or other dangers when using bath fluids that are not recommended:</li> <li>Please contact JULABO before using other than recommended bath liquids.</li> <li>JULABO assumes no liability for damage caused by the selection of an unsuitable bath fluid.</li> <li>Unsuitable bath fluids are fluids which, e.g.,</li> <li>are highly viscous (much higher than 30 mm<sup>2</sup> x s<sup>-1</sup> at the respective working temperature)</li> </ul>
	have a low viscosity and have creep characteristics
	<ul> <li>have corrosive characteristics or</li> <li>tend to crack.</li> </ul> No liability for use of other bath fluids!
ATTENTION:	The maximum permissible viscosity is 30 mm <sup>2</sup> /s

### 5.3. Filling




#### Notice:

Pay attention to the thermal expansion of bath oil during heating to avoid overflowing of the liquid.

Check to make sure that the drain port (25) is closed.

Remove the cap from the filling opening (9).



Filling of the external system.
Turn the mains switch (1) on

•

•

•


• Press the key  $\frac{1}{6}$  for filling the cooling loop for the external system. <u>Make sure that air can evacuate from the system.</u>

Connect the tubing from the external system to the pump connectors

Fill the bath tank using a funnel while monitoring the filling level (6).

- Check the filling level (6) and keep on filling the bath liquid using the funnel.
- After having finished the filling process, the liquid level should be below "MAX".
- Close the filling opening.

and check for leaks.




### 5.4. Draining



### Notice:

- Do not drain the bath fluid while it is hot or cold! Check the temperature of the bath fluid prior to draining (by switching the unit on for a short moment, for example).
- Store and dispose of the used bath fluid according to the environmental protection laws.



- Turn the mains switch (1) off.
- Place a suitable vessel for accepting the used bath liquid underneath the drain.
- Unscrew the cup nut (A) from the drain port (25) and empty the unit completely.
- Close the drain port.

## 6. Operating procedures

### 6.1. Power connection

#### Caution:

- Only connect the unit to a power socket with an earthing contact (PE protective earth)!
- The power supply plug serves as a safe disconnecting device from the line and must always be easily accessible.
- Never operate the unit with a damaged mains power cable.
- Regularly check the mains power cables for damage.
- We disclaim all liability for damage caused by incorrect line voltages!

Make sure that the line voltage and frequency match the supply voltage specified on the type plate. Deviations of  $\pm 10$  % are permissible.

### 6.2. Switching on / Start - Stop

### Switching on:


- Turn on the mains power switch (1).
- During the self-test all segments of the VFD-Info-Display light up. Then the software version number (example: V 4x.1x) and the order number of the recirculating cooler appears. (Example: [ 95 20 025 03 P0H1 ] ).

The display "OFF" or "R OFF" indicates the unit is ready to operate.

 The recirculating cooler enters the operating mode activated before switching the recirculating cooler off: keypad control mode (manual operation) or

remote control mode (operation via personal computer).






- **Start:** Press the start/stop key  $\frac{1}{6}$ 
  - The actual bath temperature is displayed on the VFD COMFORT-DISPLAY.


Stop: Press the start/stop key

The VFD COMFORT-DISPLAY indicates the message "OFF".

Autostart: see chapter 9.3.2. A-START – Autostart

The Autostart function enables the start of the recirculating cooler directly by pressing the mains switch or using a timer.





## 7. **T** Setting the temperatures



Factory settings: SETP 1 25 °C SETP 2 37 °C SETP 3 70 °C

3 different working temperatures are adjustable. Their values are freely
selectable within the operating temperature range.

Press the  $\lfloor \mathbf{T} \rfloor$  key to call up the menu for temperature selection.

This setting may be carried out with the recirculating cooler being in the Start or Stop condition!

### Example: Setting working temperature "SETPoint 3"

- **1.** Press the **T** / key until the desired menu window is indicated on the VFD COMFORT-DISPLAY Example: SETP 3 / 70.0 °C (last digit blinks)
- 70.0

3

XXX

SETP

XXXXX	
SETP	3
85	.0

<u>SSS</u>	
7	70.7
S	85.0
E	



SETP 1 SETP 2 SETP 3 2. Change the value to 85 °C.

Use the cursor keys **t** is move left or right on the display until the numeral you wish to change blinks.

Use the increase/decrease arrows **V** to change the selected numeral (-, 0, 1, 2, 3, ... 9). see example left: SETP 3 / 85.0.

- **3.** Press enter **C** to store the value.
- In the >Start< condition this value is immediately used for controlling the working temperature.</p>

The indication on the VFD COMFORT-DISPLAY is updated.

The heater control indicator blinks.

Notice: See SETMAX and SETMIN in chapter 9.5. MENU LIMITS

### Example: Selecting the working temperature

**1.** Press the  $\Box$  / key until the desired menu item is indicated on the VFD COMFORT-DISPLAY.

**2.** Press enter **C**.

(1) The recirculating cooler uses the new working temperature value for temperature control.

## 8. Safety installations, warning functions

	Check the safety installation at least twice a year! See page 23
40,, 40,, 230 °C	Settings for the excess temperature protection > SAFETMP< according to IEC 61010-2-010 and for the high > OVERTMP< and Iow> SUBTMP< temperature warning functions are made in a menu that is called up with the key .
SAFETMP	
> OVERTMP	Menu item > LIMITSR <: "Warning" or "Alarm"
SUBTMP	For the two menu items > OVERTMP< and >SUBTMP< choose between
➢ LIMITSR	a warning message being signalled or a complete shutdown of the main functional elements such as heater and circulating pump being effected.

### 8.1. Excess temperature protection





This safety installation is independent of the control circuit. When the temperature of the bath fluid has reached the safety temperature, a complete shutdown of the heater and pump is effected.

The alarm is indicated by optical and audible signals (continuous tone) and on the VFD COMFORT-DISPLAY appears the error message "ALARM-CODE 1".

Setting range: 20 °C ... 220 °C

- Press the button until the menu item > SAFETMP < is displayed.</li>
- Set the new cut-out value using a screwdriver via the VFD COMFORT-DISPLAY (Example: 100 °C)
- Press ESC to update the display immediately, or the unit automatically returns to the effective display after about 30 seconds ①.

#### **Recommendation:**

Set the excess temperature protection at 5 to 10 °C above the working temperature setpoint.



#### Warning:


ing: 🗢 📭

The excess temperature protection >SafeTemp< should be set at least 25 °C below the fire point of the bath fluid used. There is a risk of fire in the event of a wrong setting! We disclaim all liability for damage caused by wrong settings!

### 8.1.1. Early warning system, low level protection



(patented)



This low level protection is independent of the control circuit and is divided in two sections.

 Switch in stage 1 recognizes a critical fluid level . An audible warning (interval tone) sounds and a message appears on the VFD COMFORT-DISPLAY.



Switch in stage 2 recognizes a low fluid level O.
 If stage 2 of the low level protection device (according to IEC 61010-2-010) is triggered, a complete shutdown of the heater and circulating pump is effected.

A continuous alarm tone sounds and a message >ALARM< >CODE 01< appears on the VFD COMFORT-DISPLAY.

XXXXX		
ALARM		
CODE	1	

Turn off the unit with the mains switch, refill bath fluid and turn the unit on again!

**Important:** Check the safety installation at least twice a year! See page 23.



### Warning:

When adding bath fluid, always use the same bath fluid type that is already in the bath.


Bath oils must not contain any water and should be pre-heated approximately to the current bath temperature! Explosion hazard at high temperatures!

### 8.2. Over and Sub temperature warning functions

If for a sensitive temperature application task adherence to a working temperature value > SETP < is to be supervised, then set over and sub temperature warning values.

In the example below, the > SETP < of 85 °C is surrounded by the values > OVERTEMP < 87 °C and > SUBTEMP < 83 °C. The electronics immediately registers when the actual temperature attains a temperature out of the limits and it follows a reaction according to what is set in the menu item >LIMITSR<.

(see chapter 8.2.1. Change-over of the warning function to shutdown function)



- **1.** Press the → button until the menu item > OVERTMP < or >SUBTMP< is displayed.
- **2.** Set value:

Use the cursor keys **T** to move left or right on the display until the numeral you wish to change blinks.

Use the increase/decrease arrows  $\checkmark$  to change the selected numeral (-, 0, 1, 2, 3, ... 9).

- **3.** Press enter **C** to store the value.
  - The warning functions are only triggered when the actual bath temperature, after start from the "OFF" or "rOFF" mode, lies within the set limits for 3 seconds.


#### **Recommendation:**

Set the high temperature warning value > OVERTMP < at 5 °C to 10 °C above the working temperature setpoint.

Set the low temperature warning value > SUBTMP < at 5  $^{\circ}$ C to 10  $^{\circ}$ C below the working temperature setpoint.



Over temperature




### 8.2.1. Change-over of the warning function to shutdown function

For the two menu items > OVERTMP< and >SUBTMP< choose between a warning message being signalled or a complete shutdown of the main functional elements such as heater and circulating pump being effected (see page 32).

#### • Setting >WARNING<

An audible warning (**interval tone**) sounds and a meassage appears on the VFD COMFORT-DISPLAY.




Factory setting: >WARNING<

XXX . X

LIMITSR WARNING


XXX.X LIMITSR

ALARM




Setting >ALARM<

A complete shutdown of heater and circulating pump is effected. An audible alarm (**continuous tone**) sounds and a message appears on the VFD COMFORT-DISPLAY.




- 2. Select the parameter with the keys ▼▲. (>WARNING< or >ALARM<)
- **3.** Press enter **C** to store the selected parameter.

9	. Mer	nu functions	
M	ENU	The term "menu functions" refers to adjustments such as	
>	PUMP	Electronically adjustable pump capacity	page 35
	CONTROL	Intelligent Cascade Control, control parameters CONTROL - internal or external control ELFTUNING DYNAMIC - internal COSPEED - external Control parameters - XP, TN, TV internal Control parameters - XP, TN, TV, XPU external	page 36
	CONFIG	Configurationen of the unit SET (Setpoint) – keypad control or remote control A-START – Autostart OFF MODE – Motor on / off RESET – Factory settings ACTVAR – actuating variable TIME / DATE – Setting time and date	page 42
~	SERIAL	Adjustable interface parameters BAUDRAT, H-SHAKE, PARITY (Baud rate, Handshake, Parity)	page 48
>	LIMITS	Limits to temperature or capacity SET MAX / MIN - Maximum and minimum setpoint HEAT MAX - Adjusted maximum heating COOL MAX - Adjusted maximum cooling INTERN MAX / MIN – Limitation of the working temperature r BAND HIGH / LOW – Band limit	page 49 ange
	PROGRAM	Integrated programmer	page 51
>	ADJUST	ATC - Absolute Temperature Calibration, Sensor calibration, 3-point calibration	page 55
the	ANALOG dication only when e Analog Interface odul is mounted.		page 59




#### **MENU PUMP - Setting the pump pressure** 9.1.





Pump capacity: See table 1 page 14

### 9.2. MENU CONTROL – Control parameters



#### 9.2.1. CONTROL – internal / external control

XXX.X
CONTROL
INT

The recirculating cooler is conceived for internal and external temperature control. Switching is carried out in this submenu. Depending on what is set, only the respective set of parameters is indicated.

Factory setting: INT

Possible parameters:

- **INT** internal temperature control
- **EXT** external temperature control with external Pt100 sensor

(1) The control type can only be adjusted in the -OFF- condition

- **1.** Press the button until the submenu >CONTROL< is displayed.
- **2.** Select the parameter with the keys  $\checkmark$  (INT / EXT).
- **3.** Press enter **Continue:** Press **ESC** or **C**.



## Notice:

Place the external sensor into the bath medium and securely fix the sensor.



Pt100

## IMPORTANT:

Additional measures for external temperature control.

- For external control and temperature measurement an external Pt100 sensor must be connected to the socket (17) on the rear of the recirculating cooler.
- Suggested adjustments for external temperature control: BAND HIGH / LOW and INTERN MAX / MIN see chapter > LIMITS < on page 49.</p>
- Sensor calibration of the external Pt100 sensor is carried out in the >MENU / ADJUST<, in the submenu >ATC SEN / EXT< (see page 55).

#### Accessory: Pt100 external sensor

Order No.	Description	Material	Cable
8981003	200x6 mm Ø,	stainless steel	1.5 m
8981005	200x6 mm Ø,	glass	1.5 m
8981006	20x2 mm Ø,	stainless steel	1.5 m
8981010	300x6 mm Ø,	stainless steel	1.5 m
8981015	300x6 mm Ø,	stainless steel / PTFE coated	3 m
8981013	600x6 mm Ø,	stainless steel / PTFE coated	3 m
8981016	900x6 mm Ø,	stainless steel / PTFE coated	3 m
8981014	1200x6 mm Ø,	stainless steel / PTFE coated	3 m
8981103	Extension cable for	Pt100 sensor	3.5 m
8981020	M+R in-line Pt100	sensor	

Pt100

M+R

The M+R in-line Pt100 sensor is a flow sensor and can be installed loop circuit

#### 9.2.2. SELFTUNING

XXXXXX SELFTUN ALWAYS

Factory setting: ONCE

#### Selftuning:

When performing a selftuning for the controlled system (temperature application system), the control parameters Xp, Tn and Tv are automatically determined and stored.

Possible parameters:

OFF - no selftuning

The control parameters ascertained during the last identification are used for control purposes.

**ONCE** - <u>single</u> selftuning (factory setting)

The instrument performs a single selftuning of the controlled system after each start with the start/stop key  $\boxed{\frac{1}{6}}$  or after receiving a start command via the interface.

ALWAYS - continual selftuning


The instrument performs a selftuning of the controlled system whenever a new setpoint is to be reached.

Use this setting only when the temperature application system changes permanently.

- **1.** Press the  $\frac{MENU}{D}$  button until the submenu > SELFTUN < is displayed.
- Select the parameter with the keys ▼▲. (ALWAYS / OFF / ONCE)
- **3.** Press enter **C** to store the selected parameter.

Continue: Press ESC or

#### 9.2.3. DYN INT - Dynamic internal



The control parameters preset in factory are in most cases adequate for achieving an optimum temperature pattern for the samples requiring temperature application.

Each parameter may be manually set via the keypad if necessary, to allow optimum control performance.

#### Proportional range >Xp<

The proportional range is the range below the selected temperature value in which the control circuit reduces the heating power from 100 % to 0 %.

Setting	range:0.1	99.9
---------	-----------	------

1.0

XX	XXX
TN	INT
	100

XXXXX

XP INT

Setting range:1 ...9999

#### Resetting time >Tn< (Integral component)

Compensation of the remaining control deviation due to proportional regulation. An insufficient resetting time may cause instabilities to occur. Excessive resetting time will unnecessarily prolong compensation of the control difference.

XXXXX	
TV	INT
	5

Setting range: 0 ... 999

#### Lead time >Tv< (Differential component)

The differential component reduces the control settling time. An insufficient lead time will prolong the time required to compensate for disturbance effects and cause high overshooting during run-up. An excessive lead time could cause instabilities (oscillations) to occur.


- **1.** Press the  $\frac{MENU}{D}$  button until the desired submenu is displayed XP INT, TN INT, TV INT.
- 2. Set value:

Use the cursor keys **t** to move left or right on the display until the numeral you wish to change blinks.

Use the increase/decrease arrows  $\checkmark$  to change the selected numeral (-, 0, 1, 2, 3, ... 9).


**3.** Press enter **Continue**: Press MENU / **ESC** or **Continue**: Press MENU / **P** = Pre

#### Optimization instructions for the PID control parameters:



The heat-up curve reveals inappropriate control settings.

Inappropriate settings may produce the following heat-up curves:



#### 9.2.5. COSPEED - external

XXXXX
COSPEED
0.1

This parameter affects the temperature pattern only in case of **external** control.

Possible parameters: 0.0 ... 5.0


- **1.** Press the  $\frac{\text{MENU}}{\text{button until the submenu}}$  button until the submenu > COSPEED < is displayed.
- 2. Set value:

Use the cursor keys **t** to move left or right on the display until the numeral you wish to change blinks.

Use the increase/decrease arrows  $\checkmark$  to change the selected numeral (-, 0, 1, 2, 3, ... 9).

**3.** Press enter **Continue**: Press MENU / **ESC** or **Continue**: Press MENU / **P** = Pre

During selftuning, the control parameters Xp, Tn and Tv of a controlled system are automatically determined and stored. Depending on the controlled system, time for tuning can be unequally longer. This controller layout allows protection of sensitive objects requiring temperature application.



As soon as a co-speed factor is set, it is considered for calculating the control parameters. As shown in the diagram, tuning times become shorter the higher the co-speed factor is, but overshooting can happen in the internal system.

#### 9.2.6. Control parameters – XPU, XP, TN, TV external

XXXXX		
XP	EXT	
	0.7	

Setting range: 0.1 ... 99.9

XXXXX	
TN EXT	
720	
Setting range: 1	9999


XXXXX		
TV	EXT	
	55	


Setting range: 0 ... 999

XXXXX	
XPU	EXT
	5.0

Setting range: 0.1 ... 99.9

#### 9.3. **MENU CONFIG - configuration**





XX	XXX	
TV	EXT	
	55	

XXXXX		
XPU	EXT	
	5.0	

# **1**. Press the button until the desired submenu is displayed -XP EXT, TN EXT, TV EXT, XPU EXT.

Each parameter may be manually set via the keypad if necessary, to

The control parameters preset in factory are in most cases adequate for achieving an optimum temperature pattern for the samples requiring

2. Set value:

temperature application.

allow optimum control performance.

- to move left or right on the display Use the cursor kevs until the numeral you wish to change blinks.
- Use the increase/decrease arrows  $\checkmark$   $\checkmark$  to change the selected numeral (-, 0, 1, 2, 3, ... 9).
- **3.** Press enter **C** to store the value. Continue: Press MENU / ESC or 🤅

#### Proportional range >Xpu<

The proportional range Xpu of the cascaded controller is only needed for external control.

# 9.3.1. SETPOINT – Keypad control or remote control

XXXXXX SETP KEY Factory setting: KEY The selected mode is indicated on the VFD	<ol> <li>The recirculating cooler provides four possibilities for setpoint setting.</li> <li>Press the <sup>MENU</sup> button until the submenu &gt; SETP &lt; is displayed.</li> <li>Select the parameter with the keys ▼▲ (KEY / SERIAL / PT100 / EPROG)</li> <li>Press enter To store the selected parameter. Continue: Press <sup>MENU</sup> / ESC or <sup>(1)</sup></li> </ol>
COMFORT-DISPLAY	<b>KEY</b> – Setpoint setting with the keys and <b>a</b> and <b>b</b>
	<b>SERIAL</b> - Setpoint setting via the serial RS232 interface through a PC or superordinated data system.
R -OFF- s xx.xx z xx.xx serial	Important: Connect the recirculating cooler to a PC using an interface cable. Check the interface parameters of both interfaces (on recirculating cooler and PC) and make sure they match. (see chapter 12.1. Setup for remote control page 69)
-OFF- ST XX.XX T PT100	<b>PT100 -</b> Setpoint setting via the analog socket "ext. Pt100" using an external temperature sensor or an appropriate voltage/current source.
	<b>EPROG</b> - Can only be adjusted when an electronic module with analog connections is used (option). Setpoint setting via the analog interface REG+E-PROG connection with an external voltage or current source or a programmer.
I XX EPROG	Important: Connect the external voltage or current source or a programmer to the circulator via the socket REG+E-PROG (see page 62). In the menu >MENU ANALOG< set the parameter >EPROG< and the input variables (see page 63).
	The E-Prog input can only be used either under menu item > <b>SETP</b> < or under menu item > <b>ACTVAR</b> < (see page 46).

#### 9.3.2. A-START – Autostart

XXXXX
A-START
OFF

The AUTOSTART function (automatic start mode) is allowing the start of the instrument directly by pressing the mains power switch or using a timer.

- **1.** Press the  $\stackrel{\text{MENU}}{\longrightarrow}$  button until the submenu > A-START < is displayed.
- **2.** Select the parameter with the keys  $\checkmark$  (OFF / ON).
- **3.** Press enter **C** to store the selected parameter.

Continue: Press MENU / ESC or

Possible parameters:

**ON** - AUTOSTART on

**OFF** - AUTOSTART off

#### Note:

The temperature system has been configured and supplied by JULABO according to N.A.M.U.R. recommendations. This means for the start mode, that the unit must enter a safe operating state after a power failure (non-automatic start mode). This safe operating state is indicated by **OFF**, resp. **R OFF** on the VFD COMFORT-DISPLAY. A complete shutdown of the main functional elements such as heater and circulating pump is effected simultaneously.

The values set on the recirculating cooler remain stored, and the unit is returned to operation by pressing the start/stop key (in manual control mode).

In remote control mode, the values need to be resent by the PC via the interface.

Should such a safety standard not be required, the AUTOSTART function (automatic start mode) may be activated, thus allowing the start of the instrument directly by pressing the mains power switch or using a timer.

The AUTOSTART function can only be used, if setpoint setting is carried out via > KEY <, > EPROG < or >PT100<.



#### Warning:

For supervised or unsupervised operation with the AUTOSTART function, avoid any hazardous situation to persons or property. The instrument does no longer conform to N.A.M.U.R. recommendations.

The safety and warning functions of the instrument should always be used to their fullest capacity.

Factory setting: OFF

#### 9.3.3. OFF-MODE – Pump motor on / off

XXXXX	
OFFMODE	
PMP OFF	

Normally the circulating pump is switched via the start/stop signal. However, if circulation should be maintained also for the -OFF- condition, the parameter **>PMP ON<** needs to be set.

Possible parameters:

Factory setting: PMP OFF **PMP ON** Pump motor on

PMP OFF Pump motor off

- **1.** Press the  $\stackrel{\text{MENU}}{\longrightarrow}$  button until the submenu > OFFMOD < is displayed.
- 2. Select the parameter with the keys ▼▲. (PMP ON / PMP OFF)
- **3.** Press enter **C** to store the selected parameter.


Continue: Press MENU / ESC or

(1) In case of an alarm state, a shutdown of the pump motor is still effected.

#### 9.3.4. RESET – Factory settings

XXXXX
RESET
NO

Factory setting: NO



Possible parameters:

NO / YES

- **1.** Press the  $\frac{MENU}{D}$  button until the submenu > RESET < is displayed.
- 2. Select the parameter with the keys  $\checkmark$  (NO / YES).
- **3.** Press enter **C** to store the selected parameter.



(1) As long as the message -RUN- appears all parameters a reset to factory settings.

#### 9.3.5. ACTVAR - actuating variable



Factory setting: CONTROL The variable (*ACTuating VARiable*) corresponds to the extent to which the heater or cooling machine of the recirculating cooler is controlled. Heat or cold is applied to the bath according to this variable. If this happens with the control electronics of the recirculating cooler, called > CONTROL < in this particular case, the bath temperature is exactly heated and maintained constant at the adjusted setpoint.

- Programming of variables for the parameters > SERIAL < or > EPROG < is only accepted, if the unit is in Start mode.</p>
- **1.** Press the button until the submenu > ACTVAR < is displayed.
- 2. Select the parameter with the keys (CONTROL / SERIAL / EPROG )
- **3.** Press enter **Continue:** Press MENU / **ESC** or **Continue:** Press MENU

Possible parameters:

- **CONTROL** The internal control electronics of the recirculating cooler controls the heater and cooling machine. Self-tuning is possible.
- **SERIAL** The heater or cooling machine receives the control signal via the serial interface. Self-tuning is not possible.
- **EPROG -** The heater or cooling machine receives the control signal via the E-Prog input. Self-tuning is not possible.

#### Important:

Under >MENU ANALOG< set the input variable to >EPROG x / ACTVAR< (see page 63).



ACTVAR Example: EPROG U / ACTVAR

#### Note:

The E-Prog input can only be used either under menu item > **SETP** < (page 43) or under menu item > **ACTVAR** <.

XXXXXX ACTVAR SERIAL

XXXXXX ACTVAR EPROG



#### Warning:

The working temperature range of the recirculating cooler is determined during configuration. If set to >CONTROL<, this range cannot be exceeded.

If set to > **SERIAL** < and > **EPROG** <, heat or cold is applied to the bath without control. The permissible maximum temperature can be exceeded. The user has to take adequate precautions for temperature control.

Materials, such as gaskets or insulations for example, may be damaged or destroyed, if the permissible maximum temperature is exceeded.

The safety and warning functions of the instrument should always be used to their fullest capacity.\_\_\_\_

See chapter 8. Safety installations, warning functions, page 30.

#### 9.3.6. TIME / DATE – setting time and date


WWWW	The integrated clock allows starting a profile at any date and time. The clock is preset in the factory.
XXXXX <i>TIME</i> 16h43,17 hh mm	<ol> <li>Press the MENU button until the submenu &gt; TIME &lt; or &gt; DATE &lt; is displayed.</li> </ol>
	2. Setting time / date:
XXXXXX DATE	Use the cursor keys to move left or right on the display until the numeral you wish to change blinks. Use the increase/decrease arrows to change the selected
10/12 <u>0</u> 3	numeral (-, 0, 1, 2, 3, 9).
TT/MM.JJ	<b>3.</b> Press enter <b>Continue:</b> Press MENU / <b>ESC</b> or <b>Continue:</b> Press MENU / <b>Conti</b>
	Clock: Only hours and minutes are set

Clock: Only hours and minutes are set. Settings are checked for plausibility.

# 9.4. MENU SERIAL - BAUDRATE, HANDSHAKE, PARITY


XXXXXX MENU SERIAL Factory settings: 4800 Bauds even hardware handshake	MENU >BAUD MENU >PARIT MENU >H-SHA	<sup>-</sup> Y< ∖KE< between recirculating cooler and a PC or a
	<ul> <li>must be identical.</li> <li>1. Press the MENU</li> <li>2. Select the paran</li> <li>3. Press enter </li> </ul>	ess system the interface parameters of bath units button until the desired menu item is displayed. neter with the keys $\checkmark$ $\bigstar$ . to store the selected parameter. $\boxed{MENU}$ / $\boxed{ESC}$ or $\bigcirc$ .
XXXXXX BAUDRAT 4800	Adjustable interfac BAUDRATE	e parameters 4800 bauds 9600 bauds
XXXXXX PARITY EVEN	PARITY	no odd even
XXXXXX H-SHAKE HARD	HANDSHAKE SOFT = HARD = Data bits = 7	software handshake hardware handshake ; Stop bits = 1





#### Menu functions

Factory settings



#### Adjusted maximum heating / cooling.

Heating and cooling powers of the recirculating cooler are adjustable. 100 % corresponds to the values in the technical specifications of the equipment.

Setting range:

HEAT MAX – 0 to 100 % in steps of 1 % COOLING MAX – 0 to 100 % in steps of 1 %

#### INTERNAL MAX / MIN Limit setting

Setting range: -94.9 °C ... +300.0 °C

The limits INT MAX and INT MIN are only valid for external control. INT MAX and INT MIN are used to limit the expected internal bath temperatures to any upper and lower values. The temperature controller cannot exceed these limits even if it would be required for reaching the temperature in the external system. Consequently the external setpoint may thus not be reached.

Sense of a limit setting:

- $\square$  Protects the bath fluid from overheating.
- ✓ Prevents an undesired alarm shutdown by the excess temperature protection >ALARM CODE 14<.</p>

Set >INTMAX:< to a value at least 5 °C below the >SAFETEMP:< value.

Protects the pump motor from high viscosity of the bath fluid at low temperatures.

For recirculating coolers: Freezing protection when using water as bath fluid.

#### BAND HIGH / LOW - Band limit

For the heat-up and cool-down phases different settings to conform to the requirements of the particular application are possible.

Setting range: 0 ... 200 °C

Using **BAND HIGH** and **BAND LOW**, the difference between the temperatures in the internal bath and the external system can be limited to any maximum value for the heat-up or the cool-down phase. During the heat-up phase the difference value always adds to the actual external temperature. During the cool-down phase, the difference value is subtracted.

Factory settings

XXXXX
INT MAX
300.00
_







Factory settings:








# 9.6. MENU PROGRAM – Integrated programmer



The integrated programmer allows any desired temperature program sequences to be realized. Such a temperature sequence is called profile. A profile consists of individual sections defined by duration (t:) and target temperature. Target temperature is the setpoint (T:), that is achieved at the end of a section. The programmer uses time and temperature difference values within a section to calculate a temperature ramp.



Select the number of profile repetitions. Example: RUNS 10

The profile can be started immediately with the parameter (NOW),

14:25 hrs

1 XXXXX PS RUNS 10

>PS RUNS< (1 ... 99)

>PS GO<

The profile is repeated 10 times.

Example below: 19 December 2003

or at the set time with the parameter (TIME).

XXXXX

PS STEP

XXXXX PS GO TIME

	The started programmer
$\begin{bmatrix} \mathbf{X}\mathbf{X}\mathbf{X}\mathbf{X}\mathbf{X}\\ s & xx. x\\ step & \mathbf{X} \end{bmatrix}$ $\begin{bmatrix} s & xx. x\\ xxhxxxx} \end{bmatrix}$ $\begin{bmatrix} s & xx. x\\ xxhxxx} \end{bmatrix}$ $\begin{bmatrix} s & xx. x\\ xxhxx} \end{bmatrix}$ $\begin{bmatrix} s & xx. x\\ xxhxx} \end{bmatrix}$ $\begin{bmatrix} s & xx. x\\ zxhxx} \end{bmatrix}$	The started programmer indicates the actually calculated setpoint S XX.X in line 2. This value increases within the time period >Px TIM< until the target temperature >Px SEP< of the section is reached. If time in a section is set to "0", the next section starts only as soon as the target temperature is reached. Use the edit keys $\checkmark$ $\checkmark$ to scroll in line 3. The display changes in intervals of approx. 4 seconds between the valid section (STEP X) and A Remaining time of the valid section or B Remaining total time: Profile x number of repetitions or C Actual bath temperature I xxx.x – internal actual value or E xxx.x – external actual value D RUN – the programmer is started or PAUSE – the profile was interrupted with the key $\checkmark$ . Time is stopped and the temperature is maintained constant at the last calculated setpoint value. Continue: Press $\checkmark$
-OFF- S xxx.x I xx.x	<ul> <li>Interrupting / Termination of a profile</li> <li>A profile can be terminated by pressing ESC.</li> <li>Power failure with the programmer started: <ul> <li>The reaction of the circulator is the same as when switched off and on again with the mains switch. The VFD COMFORT-DISPLAY indicates "OFF".</li> <li>If the AUTOSTART function is activated, the programmer starts again at a point approx. 20 seconds before the interruption took place. However, an uncontrolled change of the bath temperature happened.</li> </ul> </li> </ul>
	Setting time/date for the start
XXXXX TIME 14h25ss XXXXX DATE 19/1203	<ul> <li>&gt;TIME</li> <li>Enter the start time. Example: 14:25 hrs</li> <li>&gt;DATE</li> <li>Enter the start date. Example: 19 December 2003</li> <li>① Check the setting of the internal real time clock (see page 47).</li> </ul>
	Status at the end of the profile


	Status at the end of the prome
XXXXX PS END SETP	>PS END< (STBY / SETP) Set the status for the end of the profile. With the parameter STandBY the circulator enters the –OFF– state. With the parameter SETPoint the circulator maintains the temperature at the value of the last section.

	Compiling profiles, indicating sections
Example: Section 2          XXXXXX         PE STEP         2	<ol> <li>Press the MENU button until the submenu &gt;PE STEP x&lt; is displayed.</li> <li>1.1 Use the increase/decrease arrows ✓ ▲ to set the number of the desired section (1, 2, 3, 10).</li> <li>1.2 Press ▲ to enter menu level 3.</li> <li>Menu level 3:</li> <li>Submenu &gt;Px SEP&lt; (SETPOINT) Set a temperature value: Example: 180 °C</li> </ol>
<b>XXXXX</b> P2 SEP 180.0	<ul> <li>2.1 Use the cursor keys to move left or right on the display until the numeral you wish to change blinks.</li> <li>2.2 Use the increase/decrease arrows to change the selected numeral (-, 0, 1, 2, 3, 9).</li> </ul>
XXXXX P2 TIM 00h10	<ul> <li>2.3 Press enter to store new the value.</li> <li>2.4 Continue: Press MENU</li> <li>3 Submenu &gt;Px TIM&lt; (TIME) Set a time. Example: 10 minutes.</li> <li>3.1 Use the cursor keys to move left or right on the display until the numeral you wish to change blinks.</li> </ul>
XXXXX P2 DEL NO	<ul> <li>3.2 Use the increase/decrease arrows ▲ to change the selected numeral (-, 0, 1, 2, 3, 9).</li> <li>3.3 Press enter ▲ to store new the value.</li> <li>3.4 Continue: Press ▲</li> <li>4 Submenu &gt;Px DEL&lt; (DELETE) Standard setting (NO) Set the parameter to YES to delete the values in this section.</li> <li>4.1 If necessary, use the arrows ▲ to set the parameter to YES and press ▲.</li> </ul>
	<ul> <li>5 Press <sup>MENU</sup> to return to submenu &gt;Px SEP&lt; in menu level 3.</li> <li>6 Press ESC to return to submenu &gt;PE STEP x&lt; in menu level 2.</li> </ul>
<b>XXXXX</b> P1 xxx 	Sections without value or time indication (including the value 0) are skipped. But they can be integrated in the profile at a later time. Example: Section 1

	Deleting a profile
XXXXX PE DEL YES	>PE DEL< (YES/NO) Set the parameter to YES to delete all sections >PE STEP / 1 to 10<.

# 9.7. MENU ADJUST – ATC Absolute Temperature Calibration

ATC serves to compensate a temperature difference that might occur between recirculating cooler and a defined measuring point in the bath tank because of physical properties.



- **1.** Press the  $\stackrel{\text{MENU}}{=}$  button until the desired submenu is displayed.
- 2. Set value /parameter:


Use the cursor keys **t** to move left or right on the display until the numeral you wish to change blinks.

Use the increase/decrease arrows to change the selected numeral (-, 0, 1, 2, 3, ... 9).
or select the parameter with the keys 
3. Press enter to store the value / parameter Continue: Press MENU / ESC or .

#### Menu functions

	ATC SENSOR - INTERNAL / EXTERNAL
XXX.X ACT SEN INT	In the first submenu the ATC function is set for the >INT< internal or >EXT< external temperature sensor.
	Calibration can be carried out for the internal temperature sensor and for the external temperature sensor connected to the socket "ext. Pt100".
XXX.X ACT SEN EXT	The recirculating cooler is able to store both parameter sets. However, only the one set under this menu item is indicated.
	ATC STATUS - ON / OFF
XXX.X ATCSTAT	In the second submenu the ATC function for the temperature sensor selected above is activated >ON< or deactivated >OFF<.
OFF	>OFF< The controller of the recirculating cooler uses the original curve of the temperature sensor.
XXX.X ATCSTAT	<b>Important:</b> During the calibration process > <b>OFF</b> < needs to be set.
ON	>ON< The controller of the recirculating cooler uses the new calibration curve.
	In the ATC STATUS > ON<, the ATC calibration curve always affects the effective working temperature (also the one set via the interface).
Calibration type	CALIBRATION ART: 1, 2, 3-point calibration
XXX.X C ART x POINT	A >1-point<, >2-point < or >3-point < calibration can be carried out.
	First geometrically define the location for calibration (measuring point CT), then determine the temperature values of the calibration points. The calibration type also determines the number of pairs of values indicated on the VFD COMFORT-DISPLAY.

XXX.X TTEMP 1 x 80.00	XXX.X CTEMP 1 x 79.73	<ul> <li>2 values per calibration point</li> <li>TTEMP X: Temperature on circulator (actual value TT) 1 or 2 or 3 Defined temperature value of the calibration point. This value is simultaneously stored with &gt; CTEMP &lt; and can be</li> </ul>
XXX.X TTEMP 2 x120.00	XXX.X CTEMP 2 x119.51	<ul> <li>indicated for control purposes.</li> <li>CTEMP X: Calibration temperature (actual value CT) 1 or 2 or 3 The "Calibration value" is determined with a temperature measuring</li> </ul>
XXX.X TTEMP 3 x160.00	XXX.X CTEMP 3 x159.34	device and stored under menu item > CTEMP <.



## Example:

3-point calibration for internal control.

In the temperature range of 80 °C to 160 °C the calibration curve of the temperature sensor (TT) should be assimilated to the actual temperatures on the measuring point (CT).

#### Set controller to internal control:

(1) The control type can only be adjusted in the -OFF- condition

- 1. Press the button until the menu item > MENU / CONTROL< is displayed and press enter
  - 1.1. Under the menu >CONTROL< set the parameter to > INT < and

press enter (see page 36).

Continue: Press ESC

**2.** Press the start/stop key  $\sqrt{\frac{1}{6}}$ 





#### Menu functions

SE	<b>XX.X</b> TP 1 80.00
ME	XXXX NU DJUST
	XXX.X ACT SEN INT
	XXX . X ATCSTAT OFF
	XXX.X C ART 3 POINT
	XXX.X CTEMP 1 I 79.73
	XXX.X TTEMP 1 I 80.00
	↑ "I" for interr

↑ "I" for internal control

#### Setting working temperature SETP:

- **3.** Press **⊥** / and set the first temperature value under > SETPoint 1 < for example (example 1<sup>st</sup> value = 80 °C).
  - 3.1. Wait until this temperature is maintained constant in the bath for about 5 minutes.

#### Calibration procedure:

- 4. Press the button until the menu item > MENU / ADJUST < is displayed and press enter
  - 4.1. Set menu item >ATC SEN< to >INT< ,
  - 4.2. Set menu item >ATCSTAT< to >OFF<,
  - 4.3. Set menu item >C ART< to >3 POINT <.
- These 3 settings are maintained for the length of the 3-point calibration procedure.
- 5. Read the value of CT on the temperature measuring device and enter the respective value under menu item > CTEMP 1<
  - 5.1. Set value: (79.73 °C)

Use the cursor keys  $\checkmark$  to move left or right on the display until the numeral you wish to change blinks. Use the increase/decrease arrows  $\checkmark$  to change the selected numeral (-, 0, 1, 2, 3, ... 9).

- 5.2. Press enter and the recirculating cooler also stores the value of TT as value for >TTMP 1< (80.00 °C). The first of the 3 points is now calibrated. Continue: Press ESC.</li>
- 6. Repeat the calibration procedure for 120 °C and 160 °C. (point **3.** to **5.**).

#### Examples:

TT = 80.00 °C	TT = 120.00 °C	TT = 160.00 °C
CT = 79.73 °C	CT = 119.51 °C	CT = 159.34 °C

# 9.8. MENU ANALOG – Analog inputs/outputs






(1) In order to use the analog inputs and outputs, the recirculating cooler must be equipped with the Analog Interface Modul available as option.

This submenu enables setting of the input and output values for the programmer input and the temperature recorder outputs of socket REG+E-PROG.

The >STAND-BY< input and the >ALARM< output are configurable.

	MENU ANALOG: Press en	ter 📕 to switch to menu level 2
F-ALARM Function at alarm A-ALARM Type of alarm	MENU >F-ALARM< MENU >A-ALARM<	( <b>STANDBY</b> / ALARM / AL-STBY) ( <b>NORMAL</b> / INVERS )
EX-STBY STAND-BY input	MENU >EX-STBY <	(INACTIV / ACTIV)
ACTINT [°C] internal actual temperature value (bath temperature) ACTEXT [°C] external actual temperature value (external sensor) S-POINT [°C] active setpoint temperature POWER [%]	MENU       >CHANNEL <	(1 / 2 / 3) (ACTINT/ ACTEXT/ POWER/ S-POINT) (-99.90) (300.00) (000) (POWER/ S-POINT/ ACTINT/ ACTEXT) (0.00) (100.00) (100.00) (S-POINT/ ACTINT/ ACTEXT/ POWER) (0-20MA / 4-20MA) (-99.90) (300.00)
SETP [°C] Setpoint FLOWRAT [L/M] ACTVAR [%]	MENU >EPROG < MENU > EPROG U< MENU > UW °C<	( <b>VOLTAGE</b> (U) / <sub>CURRENT</sub> (I)) ( <b>SETP</b> / FLOWRAT / ACTVAR) (-99.90 )
Actuating variable	MENU > OW °C<	(300.00)



- **1.** Press the  $\stackrel{\text{MENU}}{\longrightarrow}$  button until the desired submenu is displayed.
- 2. Set value/parameter:

Use the cursor keys **t** to move left or right on the display until the numeral you wish to change blinks.

Use the increase/decrease arrows  $\checkmark$  to change the selected numeral (-, 0, 1, 2, 3, ... 9). or select the parameter with the keys  $\checkmark$   $\checkmark$ .

3. Press enter to store the value / parameter Continue: Press MENU / ESC or

## Alarm output

(for external alarm signal)

Possible parameters:

>F-ALARM< (STANDBY or ALARM or AL-STBY)

>A-ALARM< (NORMAL or INVERS)


This socket is a potential-free change-over contact. With the adjustments in the menu item > **F**unktion-**ALARM** < all operating conditions can be signaled without having to change the pin assignments.

Signification of the terms under menu item >Function< :

The recirculating cooler is in condition >STANDBY< or >ALARM<

• • • • • • • • •	•••	
-OFF-		XXXXX
S xx.xx		ALARM
E		CODE XX


or >AL-STBY< - both conditions are signaled.








Type of alarm



STAND-BY

XXXXX
EX-STBY
ACTIV



For > **A-ALARM** / **NORMAL** < pins 2 and 3 are connected in any case according to the selected > Function-ALARM <.

For > **A**-**ALARM** / **INVERS** < pins 2 and 1 are connected in any case according to the selected >Function<.

Switching capacitymax.30 W / 40 VASwitching voltagemax.125 V~/-Switching currentmax.1 A

#### EX-STBY: External Stand-by input

(for external emergency switch-off) (Connector see page 68)

Possible parameters: >EX-STBY <

**INACTIV** - standby input is ignored

ACTIV - standby input is active

Activate the standby input:

- 1. Under menu item > EXT-STBY <, set the parameter to >ACTIV<.
- 2. Connect an external contact ,AK' (e.g. for emergency switch-off) or an alarm contact of the superordinated system.

In case the connection between Pin 2 and Pin 3 is interrupted by opening the contact ,AK', a complete shutdown of the circulating pump and heater is effected, and the unit enters the condition **"E OFF**".

If the contact is reclosed, the instrument returns to the standby state and **"E OFF**" is displayed. Press  $\boxed{\frac{1}{6}}$  to start.

## (i) Additional tips for using the STANDBY input:

The standby function can be used in conjunction with the AUTOSTART feature (see page 44)

- 1. If the Autostart function is NOT turned ON, the standby input is used as described above.
- 2. If the Autostart function is enabled, the instrument will revert back to the original method of entering the setpoint (i.e. keypad, RS232, Analog input, etc.).


#### Entering the setpoint with the keypad. e.g.

As described above, a bipolar shutdown is accompanied by displaying the **"E OFF**" state. The recirculating cooler starts again when the contact is reclosed. The temperature of the bath fluid changed during the **"E OFF**" state.


Entering the setpoint with the programmer (see pages 51). The display **E OFF** appears. The setpoint value and the time are both held at the current value. The temperature of the bath fluid will be held constant at this temperature. The programmer continues once the contact is reclosed.

Notice: This is not an actual shutoff feature.

#### Menu functions







# REG+E-PROG :

#### Outputs of the connector REG+E-PROG

- **1.** Select CHANNEL 1, 2 or 3
- First define the desired output value for CHANNELs 1 to 3:
   ACTINT internal actual temperature value (bath temperature))
   ACTEXT external actual temperature value (external sensor)
   POWER periodic or intermittent heating or cooling
  - S-POINT active setpoint temperature

(SETPoint 1, 2, 3,/ integr. programmer /external programmer)

3. Then select the display size for CHANNELs 1 to 3:

#### Voltage outputs CHANNELs 1 and 2

Assign the voltage values of 0 V to the lowest and 10 V to the highest necessary temperature (°C) or power rating (%) required as an output value.


#### **Current output CHANNEL 3**

The current output (channel 3) offers 2 ranges for selection: **0 mA to 20 mA** or

#### 4 mA to 20 mA


Assign the current values 0 mA or 4 mA to the lowest and 20 mA to the highest temperature (°C) or power rating (%) required as an output value.

XXXXX	ХХХХХ	XXXXX	XXXXXX	or
CH3	СНЗ	CH3 OMA	CH3 20A	
ACTINT	0-20МА	xx.xx	xx.xx	
XXXXX CH3 ACTINT	ХХХХХ СНЗ 4-20МА	XXXXX CH3 4MA xx.xx	XXXXXX CH3 20A xx.xx	 (20A ≙ 20 N



Example 1:

	lowest temperature value:	10 °C
٢	highest temperature value Fig. shows 200 °C scaled to pa rise: 50 mV/°C	210 °C aper width
Example 2:	lowest temperature value: highest temperature value: Fig. shows 5 °C scaled to pape rise: 2000 mV/°C	197 °C 202 °C er width



I xx.xx

#### E-PROG - input

Setting needs to be carried out, if

>ACTVAR< to >EPROG<.

1. setpoint programming is to be made via an external voltage or current source or programmer.

For this, in the menu > MENU / CONFIG < first set the menu item > SETP < to >EPROG<.

 the heater variable should be controlled via an external control pulse.
 For this, in the menu > MENU / CONFIG < set the menu item</li>

3. the signal of an external flow meter should be registered.

The E-Prog input can only be used either under menu item > SETP
 or under menu item > ACTVAR <.</li>
 If the input is neither occupied by > SETP < or > ACTVAR <, the signal of a flow meter can be connected.</li>

(i) First set > MENU / CONFIG <, then > MENU / ANALOG <.

- 1. Connect the external voltage or current source or programmer to socket (12) REG+E-PROG of the circulator.
- **2.** Selecting the signal:

The programmer (E-PROG) input of the circulator can be matched to the output signal of the external voltage or current source. **VOLTAGE** voltage input = U **CURRENT** current input = I

**3.** Define the input variable:

SETP	Setpoint programmed by external voltage or current
FLOWRAT	Signal of an external flow sensor
ACTVAR	Manipulated (ACTuated) variable for the heater
	with an external control pulse

## 4. UW – Setting the LOW value: (See below $\bigcirc$ )

First adjust and set the lowest voltage or current on the external voltage or current source (e.g. 0 V or 0 mA).

Then after approx. 30 secs enter the corresponding temperature value (e.g. 20.00°C).on the circulator by pressing the appropriate buttons

 $(\checkmark \checkmark \checkmark)$  on the keypad and press enter  $\checkmark$  to set.

5. OW – Setting the HIGH value: (See below ⊃)
First adjust and set the highest voltage or current on the external voltage or current source (e.g. 10 V or 20 mA).
Then after approx. 30 secs enter the corresponding temperature value (e.g. 300 °C).on the circulator by pressing the appropriate buttons (▲ ▲) on the keypad and press enter to set.

Example: SP 50.00 °C

• Set the external voltage or current source output for the equivalent of 50 °C temperature setpoint.

°C

300

250

200

150

100

80 60 40

20

Ż

2 4 6 8 10

3 4

5678

14

#### SP = Setpoint Programmer

R

9 10

18 20 mA



The temperature value adjusted and set on the external voltage or current source is displayed in line 2 of the VFD COMFORT-DISPLAY for control purposes.

This EPROG input enables the use of different voltage and current values as program parameters.

#### UW – Setting the LOW value:: (See below $\bigcirc$ )

- Adjust and set the lowest desired value on the voltage or current source resp. (Example A: 1 V). Wait appr. 30 seconds.
- 2. Assign a lower temperature threshold value to this adjusted voltage/current value by pressing the appropriate buttons on the keypad of the instrument (Example A: 20 °C ) and set by pressing enter

#### **OW – Setting the HIGH value:** (See below **그**)

- Adjust and set the highest desired value on the voltage or current source resp. (Example A: 10 V). Wait appr. 30 seconds.
- Assign an upper temperature threshold value to this adjusted voltage/current value by pressing the appropriate buttons on the keypad of the instrument (Example A: 200°C) and set by pressing enter .



(i) Example B in the diagram serves to illustrate that the end point values are freely selectable (Ex: 8 mA and 16 mA).

#### Example out of diagram A:

• Adjusting the voltage source for an output of 7.6 V!

The instrument calculates this value from the rise angle of the two predecided end points (in example A: 7.6 V correspond to an external setpoint temperature of 152.0 °C).

After returning the VFD display to standard display by pressing **ESC** this value is displayed in line 2 (Example: SP 152.00 °C).



# Notice:

If this adjustment is not correctly performed at two different points, the setpoint setting will be incorrect.



#### Important:

The usable temperature range between > **UW** < and > **OW** < is limited to the configured working temperature range of the circulator resp. of the unit combination (working temperature range see technical specifications on page 12)




# 10. Troubleshooting guide / Error messages

XXXXX ALARM CODE 01	Alarm with a complete shutdown of the unit Whenever the microprocessor electronics registers a failure, a complete shutdown of the heater and circulating pump is performed. The alarm light "A" illuminates and a continuous signal tone sounds. The VFD COMFORT-DISPLAY indicates the cause for the alarm in form of a code.
XXXXX WARNING CODE 40	<b>Warning</b> without a complete shutdown of the unit The VFD COMFORT-DISPLAY indicates the cause for the warning in form of a code and an acoustic signal sounds in regular intervals. These messages appear every 10 seconds.
	Press enter to quit the audible signal.
ALARM CODE 01	• The recirculating cooler is operated without bath fluid, or the liquid level is insufficient. Replenish the bath tank with the bath fluid.
	<ul> <li>Tube breakage has occured (insufficient filling level due to excessive bath fluid pumped out). Replace the tubing and replenish the bath tank with the bath fluid.</li> </ul>
	<ul> <li>The float is defect (e. g., because damaged in transit). Repair by authorized JULABO service personnel.</li> </ul>
	Excess temperature sensors defect.
	• The excess temperature value lies below the working temperature setpoint. Set the excess temperature to a higher value.
WARNING	Excess temperature warning     or
CODE 03	Excess temperature alarm
ALARM CODE 03	Warning type: Set to>Warning< or >Alarm< (see page 33)
WARNING	Low temperature warning
CODE 04	or
ALARM CODE 04	Low temperature alarm. Warning type: Set to>Warning< or >Alarm< (see page 33)
ALARM CODE 05	Cable of the working temperature sensor interrupted or short-circuited.
ALARM CODE 06	<ul> <li>Defect of the working or excess temperature sensor. Working temperature and excess temperature sensors report a temperature difference of more than 35 K.</li> </ul>
ALARM CODE 07	Other errors (I <sup>2</sup> C-BUS errors)

Troubleshooting guide	ror messages
ALARM CODE 12	Error in A/D converter
ALARM	Excess temperature sensors defect.
CODE 14	The excess temperature value lies below the working temperature setpoint. Set the excess temperature to a higher value.
ALARM CODE 15	External control selected, but external Pt100 sensor not connected or defect.
ALARM CODE 33	Cable of the excess temperature sensor interrupted or short-circuited.
ALARM CODE 38	Ext. Pt100 sensor input without signal, but setpoint programming set to external Pt100.
WARNING CODE 40	The early warning system for low level signals a critical fluid level. Replenish the bath tank with the bath fluid.



After eliminating the malfunction, press the mains power switch off and on again to cancel the alarm state. If the unit cannot be returned to operation, contact an authorized JULABO service station.



Mains circuit breakers on the rear of the unit

4 Mains circuit breakers (resettable)16 A (with option H5)

2 Mains circuit breakers (resettable)10 A

Motor protection circuit breaker for compressor motor

# 11. Electrical connections




# Notice:

Use shielded cables only.

The shield of the connecting cable is electrically connected to the plug housing.


The unit ensures safe operation if connecting cables with a maximum length of 3 m are used. The use of longer cables does not affect proper performance of the unit, however external interferences may have a negative impact on safe operation.



#### Socket for external Pt100 sensor

Pin assignment:		
<u>Pin</u>	Signal	
1	+	
2	U+	
3	U-	
4	I-	

The shield of the connecting cable is electrically connected to the plug housing and the sensor tube.



#### **RS232 serial interface**

This port can be used to connect a computer with an RS232 cable for remote control of the circulator.

#### Pin assignments RS232:


Pin 2	RxD
Pin 3	TxD
Pin 5	0 V
Pin 7	RTS
Pin 8	CTS

**Receive Data** Transmit Data Signal GND Request to send Clear to send

Pin 1; 4; 6, 9 Reserved - do not use!

Accessories:	Order No.	Description	
	8 980 073	RS232 interface cable 9-pol./9-pol., 2,5 m	
	8 900 110	USB interface adapter cable	

#### **Electrical connections**



#### Alarm output

(for external alarm signal)

This potential-free change-over contact is activated in case of an alarm when pins 2 and 3 are connected.

Switching capacitymax.30 W / 40 VASwitching voltagemax.125 V~/-Switching currentmax.1 A

Functional description see page 60

# STAND-BY input

(for external emergency switch-off)


Pin assignment:	Pin	Signal
-	1	not connected
	2	5 V / DC
	3	0 V

Functional description see page 6

#### Programmer input / temperature recorder output

Signal	
Channel 1	0 10 V
Channel 2	0 10 V
	0 V
EPROG	0 to 10 V / 0 to 20 mA
Channel 3	0 to 20 mA / 4 to 20 mA
	0 V
	Channel 1 Channel 2 EPROG

Functional description see page 62




#### Control connector

#### Return flow safety device.

If the liquid levels in the recirculating cooler and the external system are at different heights, overflowing must be prevented after the power has been turned off. For this reason, solenoid valves for loop circuit can be integrated in the loop circuit.

The control output is not powered in the OFF condition. Output voltage: 230 V~ / max. 0.1 A




STAND-BY

REG+E-PROG

# 12. Remote control

# 12.1. Setup for remote control



- Check the interface parameters for both interfaces (on recirculating cooler and PC) and make sure they match. (Serial interface see page 48)
- In the menu > MENU / CONFIG < set the menu item > SETPoint < to > SERIAL < .</li>

(see 9.3.1.SETPOINT – Keypad control or remote control on page 43)

3. Connect both units with an interface cable..

(a)	Like all parameters which can be entered through the keypad, interface parameters are stored in memory even after the recirculating cooler is turned
	off.

# 12.2. Communication with a PC or a superordinated data system

If the recirculating cooler is put into remote control mode via the configuration level, the VFD COMFORT-DISPLAY will read "R -OFF-" = REMOTE STOP. The recirculating cooler is now operated via the computer. In general, the computer (master) sends commands to the recirculating cooler (slave). The recirculating cooler sends data (including error messages) only when the computer sends a query.

In remote control mode: After a power interruption the order to start and all values which have to be adjusted must be resent from the personal computer via the interface.

AUTOSTART is not possible.

A transfer sequence consists of:

- command
- space (⇔; Hex: 20)
- parameter (the character separating decimals in a group is the period)
- end of file (,,; Hex: 0D)
- The response (data string) after an **in** command is always followed by a line feed (LF, Hex: 0A).

The commands are divided into **in** or **out** commands. **in** commands: asking for parameters to be displayed **out** commands: setting parameters

(P

Important times for a command transmission:

To ensure a safe data transfer, the time gap between two commands should be at least 250 ms.

The circulator automatically responds to an **in** command with a data string followed by a LF (Line Feed). The next command should only be sent after 10 ms.

The out commands are valid only in remote control mode.

Examples: Command to set the working temperature >Setpoint1< to 55.5 °C out\_sp\_00 ⇔ 55.5↓ Command to ask for the working temperature >Setpoint1< in\_sp\_00↓

Response from the circulator: **55.5**⊣ **LF** 

# 12.3. List of commands

out commands: Setting temperature values or parameters. Command Parameter **Response of recirculating cooler** 0 out\_mode\_01 Use working temperature > SETP 1< 1 Use working temperature > SETP 2< out\_mode\_01 2 out\_mode\_01 Use working temperature > SETP 3< 0 Selftuning "OFF". out\_mode\_02 Temperature control by using the stored parameters. out\_mode 02 1 Selftuning "ONCE" Single selftuning of controlled system after the next start. out\_mode\_02 2 Selftuning "ALWAYS" Continual selftuning of controlled system whenever a new setpoint is to be reached. 0 out mode 03 Set external programmer input to voltage. Voltage 0 V ... 10 V 1 out\_mode\_03 Set external programmer input to current. Current 0 mA ... 20 mA 0 Temperature control of internal bath. out mode 04 out\_mode\_04 1 External control with Pt100 sensor. 0 Stop the unit = R - OFF-. out mode 05 out mode 05 1 Start the unit.

Command	Parameter	Response of recirculating cooler		
out_mode_08	0	Set the control dynamics - aperiodic	Set the control dynamics - aperiodic	
out_mode_08	1	Set the control dynamics - standard		
out_sp_00	xxx.xx	Set working temperature. "SETP 1"		
out_sp_01	xxx.xx	Set working temperature. "SETP 2"		
out_sp_02	XXX.XX	Set working temperature. "SETP 3"	Set working temperature. "SETP 3"	
out_sp_03	xxx.xx	Set high temperature warning limit "OVERTEMP"		
out_sp_04	XXX.XX	Set low temperature warning limit "SUBTEMP"	Set low temperature warning limit "SUBTEMP"	
out_sp_06	XXX.XX	Set manipulated variable for the heater via serial interface -99.99 +100 [%]		
out_sp_07	x	Set the pump pressure stage. (1 4)		
out_par_04	X.X	CoSpeed 0 5.0 Band limit during external control. Setting the maximum difference between the temperatures in the internal bath and external system.		
out_par_06	xxx	Xp control parameter of the internal controller. 0.1 9	99.9	
out_par_07	ххх	Tn control parameter of the internal controller. 0 99	99	
out_par_08	xxx	Tv control parameter of the internal controller. 0 99	9	
out_par_09	xxx	Xp control parameter of the cascade controller. 0.1 9	99.9	
out_par_10	xxx	Proportional portion of the cascade controller. 1 99	.9	
out_par_11	xxx	Tn control parameter of the cascade controller. 0 99	99	
out_par_12	xxx	Tv control parameter of the cascade controller. 0 99	9	
out_par_13	xxx	Maximum internal temperature of the cascade controller.		
out_par_14	xxx	Minimum internal temperature of the cascade controller.		
out_par_15	xxx	Band limit (upper) 0 200 °C		
out_par_16	xxx	Band limit (lower) 0 200 °C		

Command	Parameter	Response of recirculating cooler	
version	none	Number of software version (V X.xx)	
status	none	Status message, error message (see page 74)	
in_pv_00	none	Actual bath temperature.	
in_pv_01	none	Heating power being used (%).	
in_pv_02	none	Temperature value registered by the external Pt100 sensor.	
in_pv_03	none	Temperature value registered by the safety sensor.	
in_pv_04	none	Setpoint temperature of the excess temperature protection	
in_sp_00	none	Working temperature "SETP 1"	
in_sp_01	none	Working temperature "SETP 2"	
in_sp_02	none	Working temperature "SETP 3"	
in_sp_03	none	High temperature warning limit "OVERTEMP"	
in_sp_04	none	Low temperature warning limit "SUBTEMP"	
in_sp_05	none	Setpoint temperature of the external programmer (REG+E-PROG).	
in_sp_07	none	Adjusted pump stage	
in_sp_08	none	Value of a flowrate sensor connected to the E-Prog input	
in_par_00	none	Temperature difference between working sensor and safety sensor	
in_par_01	none	Te - Time constant of the external bath.	
in_par_02	none	Si - Internal slope	
in_par_03	none	Ti - Time constant of the internal bath.	
in_par_04	none	CoSpeed - Band limit (max. difference between the temperatures in the internal bath and external system).	
in_par_05	none	Factor pk/ph0: Ratio of max. cooling capacity versus max. heating capacity	
in_par_06	none	Xp control parameter of the internal controller.	
in_par_07	none	Tn control parameter of the internal controller.	
in_par_08	none	Tv control parameter of the internal controller.	
in_par_09	none	Xp control parameter of the cascade controller.	
in_par_10	none	Proportional portion of the cascade controller.	

Command	Parameter	Response of recirculating cooler
in_par_11	none	Tn control parameter of the cascade controller.
in_par_12	none	Tv control parameter of the cascade controller.
in_par_13	none	Adjusted maximum internal temperature of the cascade controller.
in_par_14	none	Adjusted minimum internal temperature of the cascade controller.
in_par_15	none	Band limit (upper)
in_par_16	none	Band limit (lower)
in_mode_01	none	Selected setpoint:
		0 = SETP 1
		1 = SETP 2
		2 = SETP 3
		3 = Last setpoint setting was carried out
n mode 02		through an external programmer
n_mode_02	none	Selftuning type: 0 = Selftuning "OFF"
		1 = Selftuning "ONCE"
		2 = Selftuning "ALWAYS"
n_mode_03	none	Type of the external programmer input:
m_mode_00	none	0 = Voltage  0  V to  10  V
		1 = Current 0 mA to 20 mA
n mode 04	none	Internal/external temperature control:
II_IIIOUC_04	none	0 = Temperature control with internal sensor.
n mode 05		1 = Temperature control with external Pt100 sensor.Recirculating cooler in Stop/Start condition:
II_III0de_05	none	
		0 = Stop
		1 = Start
n_mode_08	none	Adjusted control dynamics
		0 = aperiodic
		1 = standard
in_hil_00	none	Max. cooling power (%).
in_hil_01	none	Max. heating power (%).

# 12.4. Status messages

Status messages	Description	
00 MANUAL STOP	Recirculating cooler in "OFF" state.	
01 MANUAL START	Recirculating cooler in keypad control mode.	
02 REMOTE STOP	Recirculating cooler in "r OFF" state.	
03 REMOTE START	Recirculating cooler in remote control mode.	

# 12.5. Error messages

Error messages	Description
-01 LOW LEVEL ALARM	Low liquid level alarm or Excess temperature protector alarm
-03 EXCESS TEMPERATURE WARNING	High temperature warning.
-04 LOW TEMPERATURE WARNING	Low temperature warning.
-05 WORKING SENSOR ALARM	Working temperature sensor short-circuited or interrupted.
-06 SENSOR DIFFERENCE ALARM	Sensor difference alarm. Working temperature and safety sensors report a temperature difference of more than 35 K.
-07 I <sup>2</sup> C-BUS ERROR	Internal error when reading or writing the I <sup>2</sup> C bus.
-08 INVALID COMMAND	Invalid command.
-09 COMMAND NOT ALLOWED IN CURRENT OPERATING MODE	Invalid command in current operating mode.
-10 VALUE TOO SMALL	Entered value too small.
-11 VALUE TOO LARGE	Entered value too large.
-12 TEMPERATURE MEASUREMENT ALARM	Error in A/D converter.
-13 WARNING : VALUE EXCEEDS TEMPERATURE LIMITS	Value lies outside the adjusted range for the high and low temperature warning limits. But value is stored.
-14 EXCESS TEMPERATURE PROTECTOR ALARM	Excess temperature protector alarm
-15 EXTERNAL SENSOR ALARM	External control selected, but external Pt100 sensor not connected.
-33 SAFETY SENSOR ALARM	Excess temperature sensor short-circuited or interrupted.
-38 EXTERNAL SENSOR SETPOINT PROGRAMMING ALARM	Ext. Pt100 sensor input without signal and setpoint programming set to external Pt100.
-40 NIVEAU LEVEL WARNUNG	Low liquid level warning in the internal reservoir.

# 13. JULABO Service – Online remote diagnosis

×

 $\nabla$ 

JULABO circulators of the HighTech series are equipped with a black box. This box is implemented in the controller and records all significant data for the last 30 minutes.

In case of a failure, this data can be read out from the unit by using special software. This software is available as a **free** download from <u>www.julabo.com</u>\ EasyBlackBox.

- Installation is easy and is performed step by step. Please observe the instructions.
- Data read-out is possible in the conditions "OFF", "R OFF" or "ALARM".
- Connect the circulator to the computer using an interface cable.
- Start the EasyBlackBox program. The program asks for the port used (COM1, .....) and the baud rate of the unit. You do not have this information on hand? Simply try it out!

The program continues to send the request until the correct settings are made.

EasyBlackBox.vi		_ <u> </u>
Julabo	Ea	syBlackBox Version 1.0
Einstellungen/Settings	Alarmspeicher/Alarms stored	Blackbox
*** TEMPERATURE LIMITS Working Temperature Rang Select Temperature Rangu SavePoti Temperature: 90	: *** : 61.00 C soint T2: 37.00 C Setpoint T3: 61.00 : *** ge: -94.90 C to 200.00 C : 94.90 C to 200.00 C : 0 : 10: -99.90 C , High Limit: 105.00 C	c
Bandlimit Up: 200.00 Bandlimit Down: 200.00 Bad-InternMax: 200.00 Bad-InternMin: -94 90	-	<b>•</b>
Speichern/Save	Hilfe/Help	Beenden/Quit

4800 Baud

NortDef.vi

Mit OK bestätigen!

the used baud rate! Confirm with OK!

 $\nabla$ 

COM1

Bitte den verwendeten COM Port und die Übertragungsrate auswählen!

Please choose desired COM port and

OK


- Data is read out and shown on the monitor divided into the sections
   >Einstellungen/Settings<,</li>
   >Alarmspeicher/Alarms stored<,</li>
   >Blackbox
  - ← see example
- After pressing >Speichern/Save<, a text file is created. The program suggests a filename ->C:\model description and barcode no.<. Modifications are possible.
- E-mail this file to <u>service.de@julabo.com</u>, JULABO's service department. JULABO is thus able to provide rapid support.

# 14. Cleaning / repairing the unit



# Caution:

- Always turn off the unit and disconnect the mains cable from the power source before cleaning the unit.
- Prevent humidity from entering into the circulator.
- Service and repair work may be performed only by authorized electricians.



Air cooled recirculating coolers:

To maintain the full cooling performance, clean the condenser from time to time.

- Switch off the unit, disconnect mains power cable.
- Remove the venting grid.
- Clean the ribbed condenser with a vacuum cleaner.
- Replace the venting grid.
- Switch on the unit.

#### **Cleaning:**

For cleaning the bath tank and the immersed parts of the recirculating cooler, use low surface tension water (e.g., soap suds). Clean the outside of the unit using a wet cloth and low surface tension water.

The recirculating cooler is designed for continuous operation under normal conditions. Periodic maintenance is not required.

The tank should be filled only with a bath fluid recommended by JULABO. To avoid contamination, it is essential to change the bath fluid from time to time.

#### Repairs

# Before asking for a service technician or returning a JULABO instrument for repair, please contact an authorized JULABO service station.

When returning the unit:

- Clean the unit in order to avoid any harm to the service personnel.
- Attach a short fault description.
- During transport the unit has to stand upright. Mark the packing correspondingly.
- When returning a unit, take care of careful and adequate packing.
- JULABO is not responsible for damages that might occur from insufficient packing.

JULABO reserves the right to carry out technical modifications with repairs for providing improved performance of a unit.

