PURA 4  
PURA 10  
PURA 14  
PURA 22  
PURA 30

**Wasserbäder – Water Baths - Bain-maries - Baños de Agua**  
**Betriebsanleitung - Operating Manual -**  
**Mode d'emploi - Manual de instrucciones**

1.950.6070ML-V1

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## **Water Baths**

### **Operating Manual**

Important: keep original operating manual for future use.

## Table of contents

<b>1</b>	<b>Use according to intended purpose.....</b>	<b>27</b>
1.1	Description.....	27
<b>2</b>	<b>Explanation of safety information .....</b>	<b>28</b>
2.1	Explanation of other information.....	28
<b>3</b>	<b>Operator responsibility – Safety recommendations.....</b>	<b>29</b>
3.1	Safety recommendations for the operator .....	29
3.2	Appropriate Operation.....	30
3.3	Use.....	30
<b>4</b>	<b>Control and function elements .....</b>	<b>31</b>
<b>5</b>	<b>Safety instructions .....</b>	<b>33</b>
<b>6</b>	<b>Preparations .....</b>	<b>34</b>
6.1	Installation.....	34
6.2	Bath fluids.....	35
6.3	Filling / Draining.....	35
<b>7</b>	<b>Commissioning.....</b>	<b>37</b>
7.1	Adjusting of temperature setpoint .....	37
<b>8</b>	<b>Error messages / Possible causes of faults .....</b>	<b>40</b>
<b>9</b>	<b>Technical specifications .....</b>	<b>41</b>
<b>10</b>	<b>Cleaning / repairing the unit .....</b>	<b>44</b>
10.1	Cleaning.....	45
10.2	Repairs .....	45
10.3	Returning a unit: .....	45
<b>11</b>	<b>Warranty .....</b>	<b>46</b>
<b>12</b>	<b>Waste disposal .....</b>	<b>46</b>
12.1	Packaging.....	46
12.2	Unit.....	46

## **Congratulations.**

You have made an excellent choice.

JULABO thanks you for the trust.

This operating manual is designed to familiarize you with the operation of our units and their possible applications. Please read the operating manual carefully.

Please call us if you have any questions about the operation of the unit or about the operating manual.

---

### **Contact:**

JULABO GmbH  
Gerhard-Juchheim-Strasse 1  
77960 Seelbach / Germany  
Tel. +49 (0) 7823 / 51-0  
Fax +49 (0) 07823 / 2491  
info.de@julabo.com  
www.julabo.com

---

## **The JULABO quality management system**

The standards for the development, production and distribution of temperature control devices for laboratory and industry use satisfy the requirements of ISO 9001 and ISO 14001.  
Registration certificate No. 01 100044846.

## **Unpacking and inspection**

If the packaging is damaged or if you discover any concealed transport damage when you have unpacked the devices and the accessories, please notify the supplier in the form of a statement of damage.

### **NOTICE**

The operating manual

- should be kept for future use.
- must be available to operating personnel at all times.

## 1 Use according to intended purpose

JULABO water baths have been designed for temperature application to specific fluids in a bath tank.

JULABO water baths are not conceived for direct temperature application to food and luxury articles or pharmaceutical and medico-technical products. Direct temperature application means: Unprotected contact of the object with the bath medium (bath fluid).

The unit is not for use in explosive atmosphere.

english

### 1.1 Description

**PID1**

- The water baths are operated via the splash-proof keypad.
- The implemented microprocessor technology allows to set and to store the setpoint that can be indicated on the temperature display.
- The PID temperature control adapts the heat supplied to the thermal requirements of the bath.
- The water baths conforms to the relevant requirements specified by European guidelines.
- The models PURA-4, -10, -14, -22, -30 provide a drain tap and handles for easy relocation.

## 2 Explanation of safety information

The operating manual contains warnings to increase safety when using the unit. The general warning sign, consisting of an equilateral triangle surrounding an exclamation sign and reproduced in various signal colors, is preceded by the signal words.

"Warning of a dangerous situation".

The significance of the danger is classified with a signal word. Read the instructions carefully and follow them.

### **⚠ DANGER**

This signal word designates a danger with a high level of risk which, if it not prevented, will result in death or serious injury.

### **⚠ WARNING**

This signal word designates a danger with a medium level of risk which, if it not prevented, may result in death or serious injury.

### **⚠ CAUTION**

This signal word designates a danger with a low level of risk which, if it not prevented, may result in minor or moderate injury.

### **NOTICE**

Designates a possibly harmful situation. If it is not prevented, the system or something near it may be damaged.

## 2.1 Explanation of other information

### **TIP**

Your attention is drawn to something special by this.

Designates user tips and other useful information.

### **Dangers at second glance**

Designates states which only occur after the start of an action and could have been prevented if the warning had been heeded.

### **Informative note**

Provides additional information.

### 3 Operator responsibility – Safety recommendations

The products of JULABO GmbH warrant a safe operation if installation, operation and maintenance are carried out according to common safety regulations. This section informs you about potential dangers that may arise from operating the circulator and also mentions the most important safety precautions.

The operator is responsible for the qualification of the personnel operating the units.

Make sure that the person who operates the water baths are instructed in the respect of the work.

The operator should be constantly informed about the dangers involved with their job activities as well as preventive actions.

Make sure all persons expected to carry out operation, installation and maintenance of the unit read and understand the safety information and operating instructions.

If you have any questions concerning the operation of your unit or the information in this manual, please contact us!

#### 3.1 Safety recommendations for the operator

- You received a product conceived for industrial use. Nevertheless, avoid strikes to the housing, vibrations, damages to the keypad foil (keys, display) or contamination.
- Make sure the product is regularly checked for proper condition. Regularly check (at least every 2 years) the proper condition of the mandatory, warning, prohibition and safety labels.
- Take care that the mains supply features a low impedance to avoid any negative affects on the instrument being operated in the same mains.
- This unit is designed for operation in a controlled electromagnetic environment. This means that transmitting devices (e.g. cellular phones) should not be used in the immediate vicinity. Magnetic radiation may influence other units with components susceptible to magnetic fields (e.g. a monitor). We recommend to keep a minimum distance of 1 m.
- Permissible ambient temperature: max. 40 °C, min. 5 °C.
- Permissible relative air humidity: 50 % (40 °C).
- Do not store in an aggressive atmosphere. Protect from contaminations.
- Do not expose to sunlight.

## 3.2 Appropriate Operation

Only qualified personnel is authorized to perform configuration, installation, maintenance and repairs of the water bath. Routine operation can also be carried out by untrained personnel who should however be instructed by trained personnel.

## 3.3 Use

Insufficient ventilation may result in the formation of explosive mixtures. Only use the unit in well ventilated areas. The unit is not for use in explosive atmosphere.

JULABO water baths have been designed for temperature application to water in a bath tank.

The heating zone approx. 4 cm from the rear panel of the device must remain unobstructed. No temperature control samples may be located in this area.

The bath may **not** be filled with flammable materials. Fire hazard! Only use non-acid and non corroding bath fluids.

When using hazardous materials or materials that could become hazardous, **the operator must** affix the enclosed safety labels (1 + 2) to the front of the unit so they are highly visible:

1		<p>Danger area. Attention! Observe instructions. (operating manual, safety data sheet)</p>
2a		<p>Carefully read the user information prior to beginning operation. <b>Scope: EU</b></p>
or		
2b		<p>Carefully read the user information prior to beginning operation. <b>Scope: USA, NAFTA</b></p>

Particular care and attention is necessary because of the wide operating range.

There are thermal dangers: Burn, scald, hot steam, hot parts and surfaces that can be touched.

	<p>Hot surface warning.</p>
--	-----------------------------

## 4 Control and function elements

1

Mains power switch

2

Display and control panel

2.1

LED temperature display,  
Actual value / Setpoint value (the decimal point flashes)  
Alarm display  
Menu indication (Parameter)

2.2

Control indicator – Heating

2.3

Control indicator – Alarm

2.4

Temperature units

2.5

Timer

2.6 **h.m.** Hours, minutes (Indication of the timer)

---

2.7  Edit keys  
Temperature setpoint increase or decrease  
Press the key briefly for step-by-step changes.  
Press and hold the key for fast change of setpoint.

---

2.8  OK key, with two functions  
1. Switch the device on/off.  
2. Store value / parameter.

---

2.9  MENU key

---

3  Right side of the device:  
Socket: USB-Host interface (Type B)  for remote control (JULABO EasyEcho) or for firmware update via PC..

---

4  Left side of the device:  
Drain tap and drain stud (inner tube diameter: 12 mm)

---

5  Hood cover

---

6  Rear side of the device:  
Mains connection: Integrated connector for voltage supply

---

Mains cable, country-specific.  
Only use original JULABO mains cables.

## 5 Safety instructions

It is important to follow these safety instructions to prevent personal injury and property damage. These instructions apply in addition to standard safety practices for working places.

- It is essential that you read the user information before starting the unit.
- Use PPE (safety gloves, safety shoes, safety goggles).
- Transport the unit carefully. The interior of the unit can also be damaged by impacts or if it is dropped.
- Do not stay under the unit when transporting and operating.
- Operation is permitted with **non-flammable** liquids only.
- Only operate the unit in rooms that are well-ventilated, dry and free of frost.
- Place the unit on a flat surface of non-flammable material.
- Do not start the unit if it is damaged or leaking.
- Compare the mains voltage and frequency with the specifications on the type plate.
- Only connect the unit to a fused mains connection via a FI circuit breaker ( $I_a=30$  mA).
- The power supply plug serves as safe disconnecting device from the power supply network and must be freely accessible at all times.
- Check the mains cable regularly for signs of damage.
- Never operate equipment with damaged mains power cables.
- Condensation that could appear in and on other units near the water bath may result in reduced operating safety.  
Be careful when setting up and operating the water bath!
- Some parts of the bath cover may become extremely warm during continuous operation. When lifting the bath cover, pay attention to hot steam!
- Never use the unit without bath fluid.
- Do not reach into the thermal bath fluid.
- Do not drain the bath fluid while it is hot!
- Check the temperature of the bath fluid prior to draining (by switching the unit on for a short moment for example).
- Always turn off the unit and disconnect the mains cable from the power source before performing any service or maintenance procedures, or before moving the unit.
- Have all service and repair work carried out by authorized specialists only.

- Switch off the unit and disconnect it from the power supply before cleaning it.
- Always empty the bath before moving the unit.

## 6 Preparations

### 6.1 Installation

---

Place the water bath in an upright position.

#### **WARNING**

##### **Explosion hazard**

The unit is not for use in explosive environment.

##### **Potential hazards from the samples**

Proper use of water baths includes immersion of samples contained in test tubes, Erlenmeyer flasks, or other containers for the purpose of controlling their temperature.

We do not know which substances are contained within these vessels. Many substances are:

- inflammable, easily ignited or explosive
- hazardous to health
- environmentally unsafe

##### **i.e.: dangerous**

The user is responsible for the handling of these substances!

- Always properly seal all sample containers.

#### **NOTICE**

**There is a danger of electrochemical oxidation or corrosion when using test-tube racks or samples made of non-ferrous metal.**

- Avoid using these types of racks or samples.
- Use only original JULABO test-tube racks.

## 6.2 Bath fluids

**Recommended bath fluids:** soft/decalcified water.

### NOTICE

#### Poor water quality may result in corrosion in the bath.

The quality of water (tap water) depends on local conditions.

- Ferrous water can cause corrosion - even on stainless steel.
- Chloric water can cause pitting corrosion.
- Distilled and deionized water is unsuitable. Their special properties cause corrosion in the bath, even in stainless steel.
- Due to the high concentration of lime, hard water is not suitable for temperature control because it leads to calcification in the bath.

JULABO takes no responsibility for damages caused by the selection of an unsuitable bath fluid.

Please contact JULABO before using other than recommended bath fluids.

#### Do not use flammable bath fluids!

#### Water bath protection products

The water bath protection product "Aqua-Stabil" is recommended to combat algae, bacteria and fungus formation.

Order No.	Designation
8 940 006	6x 100 ml bottles
8 940 012	12x 100 ml bottles

## 6.3 Filling / Draining

### ! WARNING

#### Thermal dangers when opening the bath cover:

Burn, scald, hot steam, hot parts and surfaces that can be touched.

- Do not drain the bath fluid while it is hot!
- Check the temperature of the bath fluid prior to draining, by switching the unit on for a short moment, for example.

#### Filling:

- Recommended **maximum filling level** is 2.5 cm below the tank rim.
- **Minimum filling level:** approx. 2 cm above the bath bottom.

## NOTICE


The working filling level depends on size and number of the items (fixtures) to be placed inside.

The recommended procedure is to fill the water bath only partially, place the items (fixtures) inside and then correct the filling level (adding or removing liquid) as required.

### Recommendation:

Use the water bath cover to keep temperature losses to a minimum. This is especially important for working temperatures above 60 °C.

### Draining:

- Press the mains switch to turn the water bath off
- Pull of the mains cable from the water bath device
- Place a suitable collecting bucket or tub underneath the unit for draining the used bath liquid.
- Push a hose (inner diameter 12 mm) onto the drain stud.
- To drain, turn the drain tap on the side of the water bath from position "A" to "B".
- When the water bath is completely drained, turn the drain tap back to position "A".
- Remove the hose.

## 7 Commissioning

### **⚠ WARNING**

#### **Danger from mains voltage.**

Risk of injury from electric power.

- Compare the mains voltage and frequency with the details on the model plate.
- Connect the device only to a safe power supply via FI-circuit breaker (IA = 30 mA).
- The device may only be connected to power outlets with a ground contact (PE – protective earth).
- The mains plug serves as a safe disconnecting device from the power supply network and must be freely accessible at all times.
- Do not start the device if it has a damaged mains cable.
- Check the mains cable regularly for signs of damage.
- We disclaim all liability for damage caused by incorrect line voltages!

english

### 7.1 Adjusting of temperature setpoint

Device settings

Switch-on the unit -> display *OFF*

Press **OK** -> indication of the actual value *28.5*

The temperature can be set in the start/stop condition.

The set value is saved so that it will be retained even after a power outage.

#### **Setpoint settings**

Example:


1. Press the edit key **▼** or **▲** briefly to switch from displaying the actual value to the desired value. The decimal point will flash.
2. Change the value:  
Press the **▲** key to set a higher value.  
Press the **▼** key to set a lower value.  
Press the key briefly for single steps, press and hold the key to adjust the values quickly.
3. Save the set value as the desired value by pressing the **OK** key.
4. The new desired value will flash three times.
5. A signal sounds once the setpoint is reached.

## Timer settings


### Set Timer

1. Press  "T" flashes, *OFF*
2. Press  or  *0n*
3. Press  *run*
4. Press  adjust h.m with  or 
5. Confirm with , Timer runs („T" lights up).

After the desired value is reached, the timer starts and the display alternates between the actual temperature (for 3 sec.) and the remaining time (for 2 seconds).

After the set time has elapsed, a two-time beep tone sounds.


### Stopp Timer

1. Press  "T" flashes, *0n*
2. Press  or  *OFF*
3. Confirm with , Timer stops *25.5* actual temperature

## Unit settings

### Set temperature unit

1. Press  2 twice *°C*
2. Press  or  *°C* or *°F*
3. Confirm with , temperature values are indicated with selected unit.

## 7.2 Operation of the water baths via the USB interface

- PC and PURA are switched on.
- Connect the water bath and PC with the USB cable.

For comfortable operation, you can use the JULABO software EasyTEMP.

There are also interface commands to the direct query available:

Input	Function
STATUS	current operating status
IN_SP_00	Setpoint query
IN_PV_00	Current value query
IN_PV_01	Heating capacity query
IN_SP_11	Temperature display °C (=0); °F (=1)
IN_MODE_05	Start/Stop query
OUT_SP_00	Send setpoint
OUT_MODE_05	Start/Stop device 1/0

### Remote operation

The unit can be remote-controlled via USB interface. For this purpose, connect the unit via the interface with a PC. Remote operation is activated and/or deactivated by pressing the  and  key combination (press and hold at the same time for approx. 3 seconds). During remote operation, the display alternates between the current bath temperature and " - r - ". The unit cannot be operated using the (membrane) keypad. The unit can be merely deactivated. The set remote operation is continued after the mains switch is switched on or off.

## 8 Error messages / Possible causes of faults

The following faults which trigger alarms result in the units's heater and circulating pump being shut down permanently.

The alarm indicator  lights up and a continuous signal tone will sound. The reason for the alarm or warning will be shown on the LED temperature display in coded form.

Warnings are displayed alternately with the actual value on the display.

The signal tone can be muted by pressing the OK key.

E 05

The cable for the working temperature sensor has been interrupted or short-circuited.

E 14

The cut-out value of the excessive temperature protector is below the defined working temperature.  
Set the safety temperature to a higher value.

### To cancel the alarm state

- Switch off the device at the main switch.
- Eliminate the cause of the alarm.
- Eliminate the cause of the alarm or wait for approx. 2 seconds, depending on error type.
- Switch on the device again at the mains power switch.
- If the error occurs again, a remote diagnostic is to be made.

### Faults which are not displayed:

The heater is protected from overloads. After a cooling phase the unit will restart automatically.

If necessary the unit should be inspected by a JULABO service technician.

### JULABO Technical Service

Phone: +49 (0) 07823 / 5166

Fax: +49 (0) 07823 / 5199

Email: [service.de@julabo.com](mailto:service.de@julabo.com)

## 9 Technical specifications

Water bath		PURA 4	PURA 10
Working temperature range (with cover)	°C	18 ... 99.9	18 ... 99.9
Temperature stability		±0.2	±0.2
Temperature selection		Digital	digital
Temperature indication		LED	LED
Resolution	°C	0.1	0.1
Temperature control		PID1	PID1
Heater wattage, (at 230 V)	kW	0.8	1.3
Heater wattage, (at 110 V / 115 V)	kW	0.6 / 0.8	0,9 / 1,1
Bath opening (WxL)	cm	12x27	22x27
Bath depth	cm	17	17
Filling volume	Liter	0.4 ... 4.8	1 ... 10
Overall dimensions (WxDxH) without Makrolon® -cover	cm	21x38x30	31x38x30
Bath, material		Coated aluminum, high-quality plastic	Coated aluminum, high-quality plastic
Weight	kg	5	6.2
Ambient temperature	°C	5 ... 40	5 ... 40
Mains power connection	V/ Hz	230 ±10 % / 50 / 60	230 ±10 % / 50 / 60
Current input, (at 230 V)	A	4	6
Mains power connection	V/ Hz	100-115 ±10 % / 50/60	100-115 ±10 % / 50/60
Current input (100 V / 115 V)	A	6 / 7	9 / 10

english

Water bath		PURA 14	PURA 22
Working temperature range (with cover)	°C	18 ... 99.9	18 ... 99.9*
Temperature stability		±0.2	±0.2
Temperature selection		digital	digital
Temperature indication		LED	LED
Resolution	°C	0.1	0.1
Temperature control		PID1	PID1
Heater wattage, (at 230 V)	kW	1.3	2.0
Heater wattage, (at 110 V / 115 V)	kW	0.9 / 1.1	0.9 / 1.1
Bath opening (WxL)	cm	30x27	54x27
Bath depth	cm	17	17
Filling volume	Liter	1 ... 14	2 ... 22
Overall dimensions (WxDxH) without Makrolon®-cover	cm	42x38x30	63x38x30
Bath, material		Coated aluminum, high-quality plastic	Coated aluminum, high-quality plastic
Weight	kg	9.5	13.0
Ambient temperature	°C	5 ... 40	5 ... 40
Mains power connection	V / Hz	230 ±10 % / 50 / 60	230 ±10 % / 50 / 60
Current input (at 230 V)	A	6	9
Mains power connection*	V / Hz	100-115 ±10 % / 50/60	115 ±10 % / 50/60
Current input (100 V / 115 V)	A	9 / 10	9 / 10

\* max. 92 °C at mains power connection 115 V

<b>Water bath</b>	<b>PURA 30</b>	
Working temperature range (with cover)	°C	18 ... 99.9*
Temperature stability		±0.2
Temperature selection		digital
Temperature indication		LED
Resolution	°C	0.1
Temperature control		PID1
Heater wattage, (at 230 V)	kW	2.0
Heater wattage, (at 110 V / 115 V)	kW	0.9 / 1.1
Bath opening (WxL)	cm	76x27
Bath depth	cm	17
Filling volume	Liter	2 ... 30
Overall dimensions (WxDxH) without Makrolon®-cover	cm	85x38x30
Bath, material		Coated aluminum, high-quality plastic
Weight	kg	16.5
Ambient temperature	°C	5 ... 40
Mains power connection	V/ Hz	230 ±10 % / 50 / 60
Current input, (at 230 V)	A	9
Mains power connection*	V/ Hz	100-115 ±10 % / 50/60
Current input (100 V / 115 V)	A	9 / 10

\* max. 92 °C at mains power connection 115 V

All measurements have been carried out at:  
 rated voltage and frequency  
 operating temperature: 70 °C    ambient temperature: 20 °C  
 bath fluid: water  
 Technical changes without prior notification reserved.

## Safety installations according to IEC 61010-2-010:

Excess temperature protection	130 °C - fixed value
Classification according to DIN 12876-1	class I
Alarm message	optical + audible (permanent)

## Environmental conditions according to EN 61 010, part 1:

Use only indoor.

Altitude up to 2000 m - normal zero.

Ambient temperature: +5 ... +40 °C (for storage and transportation)

### Air humidity:

Max. rel. humidity 80 % for temperatures up to +31 °C,  
linear decrease down to 50 % relative humidity at +40 °C

Max. mains fluctuation of  $\pm 10$  % are permissible.

Protection class according to EN 60 529 IP21

The unit corresponds to Class I

Overvoltage category II

Pollution degree 2

## Standards for interference resistance according to EN 61326-1

This unit is an ISM device classified in Group 1 (using high frequency for internal purposes) Class A (industrial and commercial range).

## 10 Cleaning / repairing the unit

### CAUTION

**Danger of injury during maintenance, repair and transport  
Danger from mains voltage.**

- Have all service and repair work carried out by authorized specialists only.
- Always turn off the unit and disconnect the mains cable from the power source before performing any service or maintenance procedures, or before moving the unit.
- Prevent humidity from entering into the water bath.
- Do not use alcohol-based or solvent-based cleaning agents.

## 10.1 Cleaning

---

For cleaning the bath tank and the immersed parts of the water bath, use low surface tension water (e.g. soap suds).

Clean the outside of the unit using a wet cloth and low surface tension water.

The JULABO Water Baths are designed for continuous operation under normal conditions. Periodic maintenance is not required.

The tank should be filled only with a bath fluid recommended by JULABO. To avoid contamination, it is essential to change the bath fluid from time to time.

If Aqua Stabil is used, blue residues on the inside of the bath can be removed with a cleaning agent based on anionic surfactants (e.g., dish detergent). For this purpose, fill the bath with a diluted cleaning agent solution and operate it for several hours at 60 °C.

## 10.2 Repairs

---

Before asking for a service technician or returning a JULABO instrument for repair, please contact an authorized JULABO service station.

### **JULABO Technical Service**

Phone: +49 (0) 07823 / 5166

Fax: +49 (0) 07823 / 5199

E-mail: service.de@julabo.com

## 10.3 Returning a unit:

---

When returning the unit:

- Clean the unit and, if necessary, decontaminate the unit in order to avoid endangering service personnel.
- Attach a short fault description.  
Before returning the device, please complete an online return form at <http://www.julabo.com/com/support/rma>.
- During transport the unit has to stand upright. Mark the packing correspondingly.
- When returning a unit, take care of careful and adequate packing.
- JULABO is not responsible for damages that might occur from insufficient packing.
- JULABO reserves the right to carry out technical modifications with repairs for providing improved performance of a unit.

## 11 Warranty

JULABO provides a warranty that the device will function perfectly as long as it is connected and handled correctly and as described in the operating manual.

**The warranty period is one year.**

**Free extended warranty**

2 Years Warranty  
**1Plus Warranty**  
Registration free of charge on [www.julabo.com](http://www.julabo.com)

The 1PLUS guarantee gives the user a free extended warranty to 24 months, limit to a maximum of 10,000 hours of service.

The requirement for this is that the user registers the device on the JULABO website at [www.julabo.com](http://www.julabo.com), quoting its serial number. The JULABO GmbH invoice date is used to define the warranty period.

In the event of a justified complaint the warranty is limited to refinishing or free repair or replacement at our discretion. Defective parts will be repaired or replaced free of charge if, in the event of a fault or a defect, there is evidence that it was due to a material or workmanship flaw.

All other compensation claims will be rejected.

## 12 Waste disposal

### 12.1 Packaging

---

Packaging materials must be disposed of as prescribed by the current local regulations.

### 12.2 Unit

---

In the European Economic Area (EEA) the disposal of waste equipment is regulated in the "**Directive of the European Parliament and of the Council on Waste Electrical and Electronic Equipment (WEEE)**". The current official journal on this matter is available on the European Parliament's homepage.

The symbol for the separate collection of electrical and electronic equipment is a crossed-out trash can.

Disposal with household waste (unsorted waste) or similar collections of municipal waste is not permitted!

Contact an authorized waste disposal contractor in your country.

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[info@julabo.us](mailto:info@julabo.us) [www.julabo.com](http://www.julabo.com)

The following Warranty Provisions shall apply to products sold in North America by Julabo (“**Seller**”) to the entity shown as buyer (“**Buyer**”) on Seller’s invoice.

1. **Initial Warranty.**

Upon Seller’s receipt of payment in full for the products and subject to Buyer’s compliance with the terms of sale and any other agreement with Seller relating to the products, Seller warrants to the Buyer that the products manufactured by the Seller are free from defects in material and workmanship for a period not to exceed two (2) years or ten thousand (10,000) hours of operation, whichever comes first, from the date the product is shipped by Seller to Buyer (the “**Initial Warranty**”).

2. **EXCLUSION OF ALL OTHER EXPRESS WARRANTIES; EXCLUSION OF ALL IMPLIED WARRANTIES. OTHER THAN THE INITIAL WARRANTY, NO OTHER EXPRESS WARRANTIES ARE MADE. ALL IMPLIED WARRANTIES OF EVERY TYPE AND KIND, INCLUDING BUT NOT LIMITED TO ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OR USE, ARE EXCLUDED IN ALL RESPECTS AND FOR ALL PURPOSES. SELLER DISCLAIMS AND MAKES NO IMPLIED WARRANTIES WHATSOEVER.**

3. **Exclusions.**

The Initial Warranty does not include damage to the product resulting from accident, misuse, improper installation or operation, unauthorized or improper repair, replacement or alteration (including but not limited to repairs, replacements, or alterations made or performed by persons other than Seller’s employees or authorized representatives), failure to provide or use of improper maintenance, unreasonable use or abuse of the product, or failure to follow written installation or operating instructions. Buyer must return the product’s record of purchase to the Seller or one of Seller’s authorized representatives within thirty (30) days of the date the product is shipped by Seller to Buyer in order to make a

claim under the Initial Warranty. Notwithstanding anything contained herein to the contrary, all glassware, including but not limited to reference thermometers, are expressly excluded from the Initial Warranty.

4. **Buyer's sole remedies; Limitations on Seller's Liability.** Buyer's sole and exclusive remedy under the Initial Warranty is strictly limited, in Seller's sole discretion, to either: (i) repairing defective parts; or (ii) replacing defective parts. In either case, the warranty period for the product receiving a repaired or replaced part pursuant to the terms of the Initial Warranty shall not be extended. All repairs or replacements performed by Seller pursuant to these Warranty Provisions shall be performed at Seller's facility in Allentown, Pennsylvania, U.S.A. or at the facility of an authorized representative of Seller, which location shall be determined by Seller in its sole discretion; provided, however, that Seller may, in its sole discretion perform such repairs or replacements at Buyer's facility in which case Buyer shall pay Seller's travel, living and related expenses incurred by Seller in performing the repairs or replacements at Buyer's facility. As a condition precedent to Seller's obligation to repair or replace a product part under the Initial Warranty, Buyer shall (i) promptly notify Seller in writing of any such defect; (ii) shall have returned the product's record of purchase to Seller or to one of Seller's authorized representatives within thirty (30) days of the date the product is delivered to Buyer; and (iii) assist Seller in all respects in its attempts to determine the legitimacy and basis of any claims made by or on behalf of Buyer including but not limited to providing Seller with access to the product to check operating conditions. If Buyer does not provide such written notice to Seller within the Initial Warranty period or fails to return the product's record of purchase as set forth above, Seller shall have no further liability or obligation to Buyer therefore. In no event shall Seller's liability under the Initial Warranty exceed the original purchase price of the product which is the subject of the alleged defect.

**5. THE REMEDIES PROVIDED IN THE INITIAL WARRANTY ARE THE SOLE AND EXCLUSIVE REMEDIES AVAILABLE TO THE BUYER. NOTWITHSTANDING ANYTHING TO THE CONTRARY CONTAINED HEREIN, AND EVEN IF THE SOLE AND EXCLUSIVE REMEDIES FAIL OF THEIR ESSENTIAL PURPOSE FOR ANY REASON WHATSOEVER, IN NO EVENT SHALL SELLER BE LIABLE FOR BUYER'S MANUFACTURING COSTS, LOST PROFITS, GOODWILL, OR ANY OTHER SPECIAL, INDIRECT, PUNITIVE, INCIDENTAL OR CONSEQUENTIAL DAMAGES TO**

**BUYER OR ANY THIRD PARTY AND ALL SUCH DAMAGES ARE HEREBY DISCLAIMED.**

6. **Assignment.** Buyer shall not assign any of its rights or obligations hereunder without the prior written approval of Seller; provided, however, that if Buyer is a distributor of Seller, the rights and obligations of Buyer under these Warranty Provisions shall inure to the benefit of and be binding upon Buyer's customers who provide the product's proof of purchase to Seller pursuant to the terms set forth herein. Seller may assign any or all of its rights or obligations hereunder without Buyer's prior consent.

7. **Governing Law.**

The Warranty Provisions and all questions relating to their validity, interpretation, performance, and enforcement shall be construed in accordance with, and shall be governed by, the substantive laws of the Commonwealth of Pennsylvania without regard to its principles of conflicts of law.

8. **Waiver.**

Any failure of the part of Seller to insist on strict compliance with the Warranty Provisions shall no way constitute a waiver of such right. No claim or rights arising out of a breach of the Warranty Provisions by Buyer may be discharged in whole or in part by a waiver of the claim or right, unless the waiver is in writing signed by an authorized representative of Seller. Seller's waiver or acceptance of any breach by Buyer of any provisions of the Warranty Provisions shall not constitute a waiver of or an excuse for nonperformance as to any other provision of the Warranty Provisions nor as to any prior or subsequent breach of the same provision.

9. **Freight.**

Buyer will arrange and pay for shipping and handling charges for the unit to be returned to the Seller. Seller will arrange and pay for shipping and handling for the return of the unit to the Buyer.

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*The above-named product is in compliance with the following harmonized standards and technical specifications:*

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Changes without prior notification reserved.

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**JULABO USA, Inc.**