

Determination of acidity/acid degree in milk and milk products



### Use

This method is used for the quantitative determination of acidity (acid degree) in milk and milk products. The acidity is calculated either as SH-, Dornic- or Therner degree. The different degree uses different concentrations of NaOH titrant.

### Appliances

- Titrator: TL 6000/7000 (TL 6000/7000 M2/20) consists of
- Basic device
- Magnetic stirrer TM 235
- 20 mL Exchange unit WA 20, with nrown glass bottle for titrant complete
- And pH combination electrode A 162 DIN ID

#### Electrodes

- Electrode: A 162 DIN ID
- Calibration: DIN buffer pH= 4.00 and pH= 7.00



#### Reagents

- Titrant: sodium hydroxide solution 0.1mol/l (Therner), 0.25 (SH) or 0.11 (Dornic) mol/l
- Soda lime for carbon dioxide uptake of the reagent.
- Titer: potassium hydrogen phthalate (reference material)

#### Description

#### Calibration

The pH combination electrode is calibrated in technical buffer pH=4.00 and pH= 7.00 or in DIN buffer pH= 4.01 and pH= 6.87. Example of the calibration documentation:

### Calibration

#### Buffers used

pH buffer 1:	TEC_4.000
pH buffer 2:	TEC_7.000

Measured values

pH buffer 1:	TEC_4.000	165.6 mV / 23.4 °C
pH buffer 2:	TEC_7.000	-11.2 mV / 23.0 °C

Calibration data

Slope:	99.4 % / -58.8 mV/pH	
Zero point:	pH 6.81 / -11.2 mV	
Temperature:	23.4 °C (a)	
Date and time:	07.03.13 / 15:04	

### Determination of the exact concentration of the standard solution

By carbon dioxide absorption from the air occurs in the sodium hydroxide solution of sodium bicarbonate, which changes the pH of the titrant. To prevent this, a drying tube filled with soda lime is placed on the reagent bottle. The exact concentration of the sodium hydroxide solution is determined using the standard potassium hydrogen phthalate. The potassium hydrogen phthalate is dried in the oven before the titer determination for 2 hours at 120°C and cooled in a desiccator.

#### Implementation

In a 50 mL beaker, 0.2 to 0.3g potassium hydrogen phthalate were weighed accurately and dissolved in 30 mL of dist. water with stirring. It is titrated with 0.1 mol/l sodium hydroxide solution. For 0.25 mol/l NaOH you need 0.5 - 0.75 g potassium hydrogen phthalate.



### Standardisation titration (result)



**GLP** documentation

Weight (W):	0.2049 g (m)	Factor 2 (F2):	1000.0000
Blank value (B):	0.0000 ml	Factor 1 (F1):	1.0000
Statistics:	3		



Method data

Method name:	Titer NaOH	Created at:	02/15/12 15:32:03
Method type:	Automatic titration	Last modification:	02/16/12 10:14:55
Measured value:	pH	Documentation:	GLP
Titration mode:	Dynamic		
Dynamic:	average		
Measuring speed / drift:	Normal:	minimum holding time:	02 5
reasoning opeca / and		maximum holding time:	15 s
		measuring time:	02 s
		drift:	20 mV/min
Initial waiting time:	0 s		
Titration direction:	Increase		
Pretitration:	Off		
End value:	10.500 pH		
EQ:	On		
slope value:	Steep	Value:	700
and a second		1000	

Dosing parameter			
Dosing speed:	100 %	Filling speed:	30 s
Maximum dosing volume:	30.00 ml		
Calculation formula			
Titer NaOH 0,1mol/l:	(W*F2)/((EQ1-B)*M*F1)		
Mol (M):	20.42230		
Unit:		Decimal places:	4
w - 15 040		5- h- 2 (52)	1000 0000
Weight (W):	man	Factor 2 (F2):	1000.0000
blank value (b);	0.0000 mi	Factor 1 (F1):	1.0000
Device information			
Device: TitroLine 6000			
Software version: 07 12	outh Titler NaOH 29 02 1	2-10 51 09 odf	
Dortware versions 0/_12	mur_mer_NaOH_23_02_12-10_51_03.pdf		



1/1



### Titration of the sample

Pipette accurately in 100 ml beaker either:

- 25 ml milk for SH degree
- 11 ml milk for Dornic degree
- 9 ml for Therner degree

Add app. 25 ml - 40 ml CO<sub>2</sub> free water into and mix the sample on the magnetic stirrer for a few seconds. Titrate with the suitable NaOH titrant to a fixed pH endpoint (8.2-8.7, depending on national norms).

### Example





### Method parameters: (optimized for SH degree)

Method data overall view			
Method name: Method type: Measured value: Titration mode: Linear steps:	acidity in milk Automatic titration pH End pt. 0.050 ml	Created at: Last modification: Damping settings: Documentation:	04/29/13 13:24:34 04/29/13 13:44:18 None GLP
Measuring speed / drift:	Normal:	minimum holding time: maximum holding time: Measuring time: Drift:	02 s 15 s 02 s 20 mV/min
Initial waiting time: Titration direction: Pretitration:	0 s Increase Off		
Endpoint 1:	pH 8.200	delta endpoint 1: Endpoint delay 1:	рН 1.000 5 s
Dosing parameter			
Dosing speed: Maximum dosing volume:	20.00 % 50.00 ml	Filling speed:	30 s
Unit values			
Unit size: Unit ID: Reagent: Batch ID: Concentration [mol/l]: Determined at: Expire date: Opened/compounded: Test according ISO 8655:	20ml 00360005 NaOH no Charge 1.00000 04/29/13 20:26:19 01/01/13 02/01/12 01/01/00		
Last modification:	04/29/13 13:26:21		

Device information

Device: TitroLine 7750 Serial number: 10018602 Software version: 1316

Hard Disk\result\acidity\_in\_milk\_29\_04\_13-13\_45\_40.pdf

2/2



### Notes

The acidity can also calculate as % lactic acid instead of acid degree. In this case the calculation is:

% lactic acid: (EP1-B) \* M\*T\*F1 /(W\*F2)

- EP1: ml consumption to pH endpoint
- B: blank value in ml. Here = 0
- M: molecular weight of lactic acid: 90.08
- T: exact concentration of the NaOH in mol/l (e.g. 0.2510 mol/l
- F1: 0.1 (conversion factor to %)
- W: sample amount in g 1
- F2:

If you have any questions on the application, you can feel free to contact us.

优莱博技术(北京)有限公司

地址:北京市朝阳区酒仙桥东路1号院M8号楼C厅3层301室 ■E-Mail: info@julabo.cn Web: www.julabo.cn

