

Determination of lodine number



Use

The method describes the procedure for the determination of the iodine number with the use of a catalyst. The application is suitable for fats, fatty acids, oils and different organic solvents.

Appliances

- Titrator: TL 6000/7000 (TL 6000/7000 M1/50) consists of
- Basic device
- Magnetic stirrer TM 235
- 50 mL exchange unit WA 50, with amber glass bottle for titrant complete

Electrodes

• Electrode: Pt 62 or Pt 61, or Pt 62 RG with cable L 1 A

Reagents

- Titration agent: Sodium thiosulphate solution $(Na_2S_2O_3)$ 0.1 mol/L
- Reaction solution: Iodine monochloride 0.1 mol/L in glacial acetic acid
- Solvent: Glacial acetic acid, dist. water
- Catalyst solution: Magnesium acetate solution 3 %
- Other reagents: Potassium iodine solution 10 %

Description

Preparation of Sodium thiosulfate solution 0.1 mol/L

The 0.1 mol/L titrant are prepared form ampules. Ready to use titrants are also available and recommended. The solution should be stored in a dark place/bottle.

Preparation of the lodine monochloride reaction solution

16.2 g iodine monochloride is dissolved in glacial acetic acid, and filled up to 1 liter with it.

Preparation of the potassium iodine solution

25 g KI are dissolved in dist. water filled up to 250 ml. It is recommended to prepare the solution freshly before use.

Preparation of the Magnesium acetate solution

45 g Magnesium acetate are dissolved in glacial acetic acid and filled up to 1 liter with it.



Titration

The sample weight depends to the expected iodine number (please refer to the table).

lodine no. g/100g	smpl. amount [g]
< 5 g	3
from 5 to 20	1
from 20 to 50	0.5
from 50 to 100	0.25
from 100 to 150	0.15
>150	0.1

Transfer the sample (weigh exactly to 0.0001 g) in an Erlenmeyer volumetric flask with stopper (if possible take a 200 ml one). Add 20 ml glacial acetic acid, 25.00 ml reaction solution and 10 ml catalyst solution. Close the Erlenmeyer flask with a stopper, shake/mix the solution for some seconds and let the sample stay in darkness for 5 minutes. Add 15 ml of the KI-solution and 50 ml dist water. Place the elektrode and burette tip in the sample and start the method (stirr very well). Carry out a blank titration without the sample in the same matter.

Electrode handling

If not in use, the electrolyte should be stored in the electrolyte solution. For further details, please refer to the electrode's operating instructions.



Methods

Blank value (page 1):



Calculation formula			
Blank:	EQ1 -> M01		
Statistics:	Off	Statistics	Off
		Julianes.	011

Blank value (page 2):



Method	data	overall view

Method name:	Iodine Blank	Created at:	05/08/13 13:31:57	
Method type:	Automatic titration	Last modification:	05/08/13 14:15:22	
Measured value:	mV	Damping settings:	None	
Titration mode:	Dynamic	Documentation:	GLP	
Dynamic:	Average			
Measuring speed / drift:	User-defined:	minimum holding time:	03 s	
		maximum holding time:	15 s	
		Measuring time:	03 s	
		Drift:	10 mV/min	
Initial waiting time:	0 s			
Titration direction:	Decrease			
Pretitration:	Off			
End value:	Off			
EQ:	On (1)			
Slope value:	Steep	Value:	700	

Dosing	parameter

Dosing speed:	100.00 %	Filling speed:	30 s
Maximum dosing volume:	50.00 ml		
Unit values			
Unit size:	50ml		
Unit ID:	10045135		
Reagent:	Natriumthiosulfat		
Batch ID:	no entry		
Concentration [mol/l]:	0.10000		
Determined at:	03/28/13 17:23:03		
Expire date:			
Opened/compounded:			
Test according ISO 8655:	03/14/13		
Last modification:	03/28/13 10:23:23		
C			



sample titration (page 1):

GLP documentation Titration graph mV Olive oil 375,0 350,0 325,0 300,0 275,0 250,0 29.35 ml 262.8 mV 225,0 ml 200,0 0 5,0 10,0 15,0 20,0 25,0 • mV/ml dmV/dml Method data Method name: Iodine Number Titration duration: 3 m 10 s End date: 08.05.13 End time: 15:16:40 Titration data Sample ID: Olive oil Weight: 0.25540 g End mV: Start mV: 387.7 mV 232.3 mV EQ: 29.348 ml / 262.8 mV Iodine Number: 82.14

Calculation formula

Iodine Number:	(B-EQ1)*T*M*F1/(W*F2)	Mol (M):	12.69000
Blank value (B):	45.8800 ml (M01)	Titre (T):	0.10000000 (a)
Factor 1 (F1):	1.0000	Weight (W):	0.25540 g (m)
Factor 2 (F2):	1,0000	Statistics:	Off



sample titration (page 2):

Method data overall view

Method name:	Iodine Number	Created at:	05/08/13 13:33:02
Method type:	Automatic titration	Last modification:	05/08/13 15:01:16
Measured value:	mV	Damping settings:	None
Titration mode:	Dynamic	Documentation:	GLP
Dynamic:	Average		
Measuring speed / drift:	User-defined:	minimum holding time:	03 s
		maximum holding time:	15 s
		Measuring time:	03 s
		Drift:	10 mV/min
Initial waiting time:	0 s		
Titration direction:	Decrease		
Pretitration:	Off		
End value:	Off		
EQ:	On (1)		
Slope value:	User-defined	Value:	350

Dosing parameter			
Dosing speed:	100.00 %	Filling speed:	30 s
Maximum dosing volume:	50.00 ml		
Unit values			
Unit size:	50ml		
Unit ID:	10045135		
Reagent:	Natriumthiosulfat		
Batch ID:	no entry		
Concentration [mol/l]:	0.10000		
Determined at:	03/28/13 17:23:03		
Expire date:	÷		
Opened/compounded:	-		
Test according ISO 8655:	03/14/13		
Last modification:	03/28/13 10:23:23		





Hints

If you have any questions concerning the application, you are welcome to contact us.

Literature

优莱博技术 (北京)有限公司

