

Quantitative analysis of ascorbic acid with specific titrant DCPIP

Date: 23.01.2014 Page 1 of 8



Use

This method is used to determine the content of ascorbic acid (Vitamin C) with the specific titrant 2,6-Dichlorphenolindophenol (DCPIP).

Appliances

Titrator: TL 6000/7000 (TL 6000/7000 M2/20) consists of

Basic device

Magnetic stirrer TM 235

20 mL Exchange unit WA 20, with brown glass bottle for titrant complete

Electrodes

Electrode: Pt 62 or Pt 6280 with cable L 1 A

Reagents

• Titration agent: 2,6-Dichlorphenolindophenol solution $(C_{12}H_7CI_2NO_2)_{0,001}$ mol/L

• Standardisation ascorbic acid solution $(C_6H_8O_6)$ (prepared freshly)

• Others: oxalic acid solution $(H_2C_2O_4)$ 10% and sodium acetate solution (CH_3COONa) 10%

Description

Preparation of the 2,6-Dichlorphenolindophenol solution (0,001 mol/L)

163,1 mg 2,6-Dichlorphenolindophenol are weighed into a 400 mL beaker and filled up to 250 mL with distilled water. This solution is than stirred for 20 minutes at 50°C. After this it is filltrated into a 500 mL volumetric flask. For stabilisation 50 mg potassium hydrogen carbonate may be added. After that the flask is filled up with distilled water.

chemical structure of 2,6-Dichlorphenolindophenol:

Preparation of the sodium acetate solution

10 g sodium acetate are weighed into a 100 mL volumetric flask and filled up with distilled water.

Preparation of the oxalic acid solution

10 g oxalic acid are weighed into a 100 mL volumetric flask and filled up with distilled water.

Date: 23.01.2014 Page 2 of 8



Preparation of the ascorbic acid solution

50 mg ascorbic acid (p.a.) are weighed into a 100 mL volumetric flask, dilluted in 10-20 mL oxalic acid and filled up with distilled water.

chemical structure of ascorbic acid:

Standardisation of the 2,6-Dichlorphenolindophenol

In a 50 mL beaker a mixture of 15 mL oxalic acid, 1 mL sodium acetate solution and 15 mL distilled water is made. After that 1 mL of the ascorbic acid solution is added. This solution is then titrated with the 2,6-Dichlorphenolindophenol solution.

calculation (t(DCPIP)):

$$t(DCPIP) = \frac{m(ascorbic acid)}{(V(DCPIP) \cdot 0,176)}$$

t(DCPIP): titer of the 2,6-Dichlorphenolindophenol solution m(ascorbic acid): amount of ascorbic acid [mg] (here 0,5 mg) V(DCPIP): consumption of the 2,6-Dichlorphenolindophenol solution [mL]

The factor 0,176 is the product of the concentration of the 2,6-Dichlorphenolindophenol solution and the molar mass of ascorbic acid. 1mL of the 2,6-Dichlorphenolindophenol solution complies with 0,176 mg ascorbic acid.

Titration

In a 150 mL beaker 0,5-5 g of the sample are weighed exactly. Then 40 mL oxalic acid and 1 mL sodium acetate solution are added. The whole solution is now stirred for 5 minutes. After adding 40 mL distilled water it is titrated with the 2,6-Dichlorphenolindophenol solution.

chemical equation:

 $C_6H_8O_6$

+ $C_{12}H_7CI_2NO_2$

 $C_6H_6O_6$

 $C_{12}H_9CI_2NO_2$

ascorbic acid

2,6-Dichlorphenolindophenol

dehydroascorbic acid

calculation:

1.) result in %:

$$\frac{V(DCPIP) \cdot t(DCPIP) \cdot 0,176 \cdot 100}{m(\text{sample}) \cdot 1000}$$

t(DCPIP): titre of the 2,6-Dichlorphenolindophenol solution

V(DCPIP): consumption of the 2,6-Dichlorphenolindophenol solution [mL]

m(sample): amount of the sample [g]

factor 0,176: see above

The factor 1000 (under the fraction stroke) is the result of the conversion from g into mg. Because the result should be indicated in %, the numerator has to be multiplied with 100.

2.) result in mg/100g:

$$\frac{V(DCPIP) \cdot t(DCPIP) \cdot 0,176 \cdot 100}{m(\text{sample})}$$

t(DCPIP): titer of the 2,6-Dichlorphenolindophenol solution

V(DCPIP): consumption of the 2,6-Dichlorphenolindophenol solution [mL]

m(sample): amount of the sample [g]

factor 0,176: see above

Date: 23.01.2014 Page 3 of 8




Because the result should be indicated in mg/100g, the numerator has to be multiplied with 100.

Methods

standardisation (page 1):

GLP-Dokumentation



Methodendaten

Methodenname: Titer 2,6 DCPIP Enddatum: 15.02.12 Titrationsdauer: 5 m 46 s Endzeit: 14:23:15

<u>Titrationsdaten</u>

Proben ID: ohne Start mV: 320.4 mV Einwaage: 0.5020 g End mV: 314.3 mV

EQ: 2.974 ml / 294.7 mV c; 0.958 mmol/l

Berechnungsformel

c: (W*F2)/((EQ1-B)*M*F1)

Mol (M): 0.17612

Einwaage (W): man Faktor 2 (F2): 1.0000 Blindwert (B): 0.0000 ml Faktor 1 (F1): 1.0000

Date: 23.01.2014 Page 4 of 8



standardisation (page 2):

Method data

Method name: Titer 2.6 DCPIP
Method type: Automatic titration
Measured value: mV

 Measured value;
 mV

 Titration mode:
 Linear

 Linear steps:
 0.050 ml

Created at: 04/26/13 12:09:46 Last modification: 04/29/13 12:45:51

Damping settings: None Documentation: GLP

Measuring speed / drift: User-defined: minimum holding time: 05 s

maximum holding time: 10 s
Measuring time: 02 s
Drift: 40 mV/min

10 s

Delay time:

Initial waiting time: 0 s
Titration direction: Increase
Pretitration: 1.200 ml
End value: Off

End value: Off
EQ: On (1)

Slope value: User-defined Value: 80

Dosing parameter

Dosing speed: 100.00 % Filling speed: 30 s

Maximum dosing volume: 6.00 ml

Calculation formula

Titer: (W*F2)/((EQ1-B)*M*F1) -> WA Mol (M): 0.17612

Unit: mmol/I Decimal places: 3

 Weight (W):
 man
 Factor 2 (F2):
 1.0000

 Blank value (B):
 M01
 Factor 1 (F1):
 1.0000

Blank value (B): M01 Factor 1 (F1):
Statistics: 3

Date: 23.01.2014 Page 5 of 8




7 m 22 s

14:11:50

sample titration (page 1):

GLP documentation

Titration graph



Method data

Method name: Vitamin C Titration duration: End date: 16.02.12 End time:

Titration data

 Sample ID:
 dosage
 Weight:
 1.1215 g

 Start mV:
 272.5 mV
 End mV:
 308.6 mV

EQ: 2.779 ml / 292.3 mV Vitamin C mg/100g: 41.5

Calculation formula

Vitamin C mg/100g: (EQ1-B)*T*M*F1/(W*F2)

Mol (M): 17.60000

Blank value (B): 0.0000 ml Titre (T): 0.95100000 Factor 1 (F1): 1.0000 Weight (W): man

Factor 2 (F2): 1.0000

Date: 23.01.2014 Page 6 of 8



sample titration (page 2):

Method data overall view

Method name: Vitamin C Created at: 02/16/12 13:54:17
Method type: Automatic titration Last modification: 02/16/12 13:58:25

Measured value: mV
Titration mode: Linear
Linear steps: 0.050 ml

Measuring speed / drift: 7 s

Initial waiting time: 0 s
Titration direction: Increase
Pretitration: Off
End value: Off
EO: On

slope value: user-defined Value: 80

Dosing parameter

Dosing speed: 100 % Filling speed: 30 s

Maximum dosing volume: 10.00 ml

Unit values

 Unit size:
 10ml

 Unit ID:
 10035409

 Reagent:
 La(NO3)3

 Batch ID:
 no comment

 Concentration [mol/l]:
 0.10000

Determined at: 09/21/11 23:15:06

Expire date: 12/31/21
Opened/compounded: 09/20/11
Test according ISO 8655: 01/01/00

Last modification: 09/21/11 16:15:13

Date: 23.01.2014 Page 7 of 8



Hints

Normally the solutions should be titrated under nitrogen. This is not necessary, if you work as quickly as possible. Additionally all solution (especially the ascorbic acid solution) should be prepared freshly.

If you have any questions concerning the application, you are welcome to contact us.

Literature

优莱博技术(北京)有限公司

地址:北京市朝阳区酒仙桥东路1号院M8号楼C厅3层301室

■E-Mail : info@julabo.cn
♦Web: www.julabo.cn



Date: 23.01.2014 Page 8 of 8