

Quantitative analysis of ascorbic acid with lodine titrant

Date: 21.10.2013 Page 1 of 5



Use

This method is used to determine the content of ascorbic acid (Vitamin C) with the lodine titrant in juices like orange or apple juice. The sulfite (SO₂) in the juice is masked before with glyoxal solution.

Appliances

- Titrator: TL 6000/7000 (TL 6000/7000 M1/10) consists of
- Basic device
- Magnetic stirrer TM 235
- 10 mL Exchange unit WA 10, with brown glass bottle for titrant complete

Electrodes

Electrode: Pt 1200 or with cable L 1 NN

Reagents

Titration agent: Iodine solution 0.01 mol/L

Other reagents: H₂SO₄ 25 % and Glyoxal solution 40 %

Description

Preparation of the 0.01 mol/L lodine solution

The iodine titrant is made out of 0.05 mol/l titrant solution which can made from ampuoles.

Preparation of the Glyoxal solution

A solution of 40 % Glyoxal is adjusted with NaOH to a pH value of 7.0

Preparation of the H₂SO₄ 25 %

To 75 ml distilled water are carefully added 14 ml H_2SO_4 96 % (conc.) This mixture is filled up to 100 ml with distilled water.

Titration

In a 100 or 150 mL beaker are added 50 ml sample and 2 ml of the Glyoxal solution. After 5 minutes 5 ml of the H_2SO_4 25 % are added. The mixture is directly titrated with the lodine solution to an μA endpoint.

Date: 21.10.2013 Page 2 of 5




1 m 12 s

Methods

GLP documentation

Titration graph



Method data

Method name: Vitamin C with Iodine Titration duration:

End date: 17.05.13 End time: 13:32:40

Titration data

Sample ID: Granini 5 Pattern: 50.000 ml Start μ A: -0.024 μ A End μ A: 3.046 μ A

EP: 0.846 ml/ 2.5 μA ASC: 101.31 mg/l

Calculation formula

ASC: (EP-B)*T*M*F1/(V*F2) Mol (M): 176.10000

 Blank value (B):
 0.0000 ml
 Titre (T):
 0.03400000 (m)

 Factor 1 (F1):
 1000.0000
 Pattern (V):
 50.000 ml (m)

Factor 2 (F2): 1.0000 Statistics: Off

Date: 21.10.2013 Page 3 of 5



sample titration (page 2):

Method data overall view

Linear steps:

Method name: Vitamin C with Iodine Created at: 05/17/13 12:43:14

Method type: Automatic titration Last modification: 05/17/13 13:29:30

Measured value: μA

Titration mode: d-stop Documentation: GLP

0.020 ml

Measuring speed / drift: 1 s

Initial waiting time: 0 s
Titration direction: Increase

Pretitration: 0.100 ml Delay time: 04 s

Endpoint: 2.5 μA delta endpoint: 2.0 μA

Endpoint delay: 5 s

Polarization voltage: 100 mV

Dosing parameter

Dosing speed: 20.00 % Filling speed: 30 s

Maximum dosing volume: 4.00 ml

Unit values

Unit size: 05ml Unit ID: 1296649042

Reagent: Iod
Batch ID: no Charge
Concentration [mol/l]: 0.03400

Determined at: 10/17/12 22:53:39

Expire date: 12/05/11

Opened/compounded: --Test according ISO 8655: --

Last modification: 10/17/12 15:53:47

Date: 21.10.2013 Page 4 of 5



Hints

If you have any questions concerning the application, you are welcome to contact us.

Literature

优莱博技术(北京)有限公司

地址:北京市朝阳区酒仙桥东路1号院M8号楼C厅3层301室

■E-Mail : info@julabo.cn
♦Web: www.julabo.cn



Date: 21.10.2013 Page 5 of 5