

Determination of Peroxide value

Date: 24.05.2013 Page 1 of 7



Use

The peroxide number gives the initial evidence of rancidity in unsaturated fats and oils. Other methods are available but peroxide value is the most widely used. It gives a measure of the extent to which an oil sample has undergone primary oxidation.

Appliances

- Titrator: TL 6000/7000 (TL 6000/7000 M1/10) consists of
- Basic device
- Magnetic stirrer TM 235
- 10 mL exchange unit WA 10, with brown glass bottle for titrant complete

Electrodes

Electrode: Pt 62 or Pt 61, or Pt 62 RG with cable L 1 A

Reagents

Titration agent: Sodium thiosulphate solution (Na₂S₂O₃) 0.01 or 0.001 mol/L

Solvent: Acetic acid/chloroform or acetic acid/decanol (3/2)

Other reagents: Potassium iodine solution sat.

Description

Preparation of Sodium thiosulfate solution 0.01 or 0.001 mol/L

The 0.01 or 0.001 mol/L titrant are prepared freshly from a 0.1 mol/L titrant solution. The solution should be stored in a dark place.

Preparation of the solvent mixture

For 1 L solvent mixture add 600 mL acetic acid to 400 mL chloroform or decanol (hexanole is also suitable)

Preparation of the potassium iodine solution

The solution should be freshly prepared before using. Mix 1.0 g potassium iodine and 1.3 g distilled water.

Titration

Weigh 1 g or more of the sample (exactly on 0,0001 g) in an Erlenmeyer volumetric flask. Add 30 mL of the solvent mixture to dissolve the sample. Add 0,5 mL of the potassium iodine solution, close the Erlenmeyer flask with a stopper and stir well for 60 sec. Add 30 mL of distilled water, place the electrode and titration tip in the sample and start the method (stir very well). Carry out a blank titration without the sample in the same matter before. The result of PV is calculated in milliequivalents O₂ per kg sample.

Electrode handling

If not in use, the electrolyte should be stored in the electrolyte solution. For further details, please refer to the electrode's operating instructions.


Date: 24.05.2013 Page 2 of 7



Methods

Blank value (page 1):

GLP documentation



Method data

Method name: POV Blank Titration duration: 4 m 20 s
End date: 03.05.13 End time: 16:52:27

Titration data

 Sample ID:
 Solvent
 Weight:
 1.00000 g

 Start mV:
 223.0 mV
 End mV:
 164.9 mV

EQ: 0.012 ml / 202.6 mV Blank: 0.012 ml

Calculation formula

Blank: EQ1 -> M01

Statistics: Off

Statistics: Off

Date: 24.05.2013 Page 3 of 7



Blank value (page 2):

Method data overall view

Method name: POV Blank Created at: 04/30/13 16:06:44
Method type: Automatic titration Last modification: 05/03/13 15:43:45

Measured value: mV Damping settings: None Titration mode: Linear Documentation: GLP Linear steps: 0.004 ml

Measuring speed / drift: User-defined: minimum holding time: 04 s maximum holding time: 15 s

Measuring time: 03 s

Drift: 10 mV/min
Initial waiting time: 5 s

Titration direction: Decrease
Pretitration: Off
End value: Off
EQ: Off

Dosing parameter

Dosing speed: 100.00 % Filling speed: 30 s

Maximum dosing volume: 0.20 ml

Unit values

 Unit size:
 10ml

 Unit ID:
 10035516

 Reagent:
 Na2S203

 Batch ID:
 no entry

 Concentration [mol/l]:
 0.01000

Determined at: 04/09/13 19:04:39


Expire date: -Opened/compounded: -Test according ISO 8655: --


Last modification: 04/09/13 12:04:42

Date: 24.05.2013 Page 4 of 7



sample titration (page 1):





Method data

 Method name:
 POV
 Titration duration:
 1 m 59 s

 End date:
 03.05.13
 End time:
 16:55:47

Titration data

 Sample ID:
 Olive oil
 Weight:
 1.00940 g

 Start mV:
 289.3 mV
 End mV:
 173.2 mV

EQ: 1.624 ml / 208.7 mV POV: 15.97

Calculation formula

POV: (EQ1-B)*T*M*F1/(W*F2) Mol (M): 1000.00000

Blank value (B): 0.0120 ml (M01) Titre (T): 0.01000000 (a) Factor 1 (F1): 1.0000 Weight (W): 1.00940 g (m)

Factor 2 (F2): 1,0000 Statistics: Off

Date: 24.05.2013 Page 5 of 7



sample titration (page 2):

Method data overall view

Method name: POV Created at: 04/30/13 16:11:15

Method type: Automatic titration Last modification: 05/03/13 15:14:08

Measured value: mV Damping settings: None Titration mode: Dynamic Documentation: GLP

Dynamic: Average

Measuring speed / drift: User-defined: minimum holding time: 03 s

maximum holding time: 15 s Measuring time: 03 s

Drift: 10 mV/min

Initial waiting time: 0 s
Titration direction: Decrease
Pretitration: Off
End value: Off

Slope value: Flat Value: 120

On (1)

Dosing parameter

Dosing speed: 100.00 % Filling speed: 30 s

Maximum dosing volume: 5.00 ml

Unit values

EQ:

 Unit size:
 10ml

 Unit ID:
 10035516

 Reagent:
 Na2S2O3

 Batch ID:
 no entry

 Concentration [mol/I]:
 0.01000

Determined at: 04/09/13 19:04:39

Expire date: -Opened/compounded: -Test according ISO 8655: --

Last modification: 04/09/13 12:04:42

Date: 24.05.2013 Page 6 of 7



Hints

If you have any questions concerning the application, you are welcome to contact us.

Literature

优莱博技术(北京)有限公司

地址:北京市朝阳区酒仙桥东路1号院M8号楼C厅3层301室

■E-Mail : info@julabo.cn
♦Web: www.julabo.cn



Date: 24.05.2013 Page 7 of 7