

## Quantitative analysis of hydrogen peroxide



# Application

## Content

Content.....	2
Use.....	2
Appliances.....	2
Electrodes .....	2
Reagents .....	2
Description.....	2
Method .....	4
Example .....	5
Hints.....	5

## Use

The concentration of hydrogen peroxide in a solution is determined by redox titration with potassium permanganate.

## Appliances

Titrator: TitroLine 6000/7000/7750  
 Burette: WA 20

## Electrodes

Electrode : Pt 62, Pt 6280 or Pt 62 RG  
 Electrolyte : KCl (3 mol/L)

## Reagents

Titration agent: potassium permanganate ( $KMnO_4$ ) 0,02mol/L (0.1 N)  
 Standardization: with ferrous sulphate (Fe(II))  
 other reagents: sulfuric acid ( $H_2SO_4$ ) ca. 25%


## Description

### Preparation of the sample

The sample is weighed in accurately with four decimal places into a 100 mL beaker and filled up to 60 mL with distilled water. After that 5 mL of a sulfuric acid (25%) are added.

### Titration

Hydrogen peroxide is determined by titration with potassium permanganate. The reaction runs in an acid solution and is based on the following chemical equation:



The chemical equation shows that the ratio of permanganate and hydrogen peroxide is 2:5, what is allowed for the calculation below.

# Application

calculation:

$$H_2O_2[\%] = \frac{5 \cdot c(KMnO_4) \cdot V(KMnO_4) \cdot t(KMnO_4) \cdot M(H_2O_2) \cdot 100}{2 \cdot m(H_2O_2) \cdot 1000}$$

c(KMnO<sub>4</sub>): concentration of the measure solution [mol/L] (here: 0,02 mol/L)  
V(KMnO<sub>4</sub>): value of the measure solution [mL]  
t(KMnO<sub>4</sub>): titre of the measure solution  
M(H<sub>2</sub>O<sub>2</sub>): molar mass of H<sub>2</sub>O<sub>2</sub>  
(34,0146 g/mol)  
m(H<sub>2</sub>O<sub>2</sub>): amount of the sample

That means that 1 ml 0.02 mol/l KMnO<sub>4</sub> = 1,701 mg H<sub>2</sub>O<sub>2</sub>

Please take the titration parameter out of the suitable method. To get proper results you have to make a titre for the potassium permanganate with ferrous sulfate.

# Application

Julabo



## Method

### Method data

Method name:	H2O2	Created at:	02/28/13 15:52:24
Method type:	Automatic titration	Last modification:	03/05/13 17:29:46
Measured value:	mV	Damping settings:	None
Titration mode:	Dynamic	Documentation:	GLP
Dynamic:		Average	
Measuring speed / drift:	User-defined:	minimum holding time:	05 s
		maximum holding time:	12 s
		Measuring time:	03 s
		Drift:	50 mV/min
Initial waiting time:	0 s		
Titration direction:	Increase		
Pretitration:	1.000 ml	Delay time:	20 s
End value:	Off		
EQ:	On (1)		
Slope value:	Steep	Value:	700

### Dosing parameter

Dosing speed:	15.00 %	Filling speed:	30 s
Maximum dosing volume:	50.00 ml		

### Calculation formula


H2O2:	$(EQ1-B)*T*M*F1/(W*F2)$	Mol (M):	1.70100
Unit:	%	Decimal places:	2
m-value:	$EQ2*T*M*F1/(W*F2)$	Mol (M):	1.00000
Unit:	mmol/l	Decimal places:	2
Blank value (B):	0.0000 ml	Titre (T):	1.00000000
Factor 1 (F1):	0.1000	Weight (W):	man
Factor 2 (F2):	1.0000	Statistics:	Off

# Application

## Example

sample titration:

Titrationsdiagramm



Methodendaten

Methodenname:	H <sub>2</sub> O <sub>2</sub> mit KMnO <sub>4</sub>	Titrationsdauer:	3 m 28 s
Enddatum:	17.07.12	Endzeit:	16:42:11

Titrationsdaten

Proben ID:	ohne	Einwaage:	1.0000 g
Start mV:	462.8 mV	End mV:	932.2 mV
EQ:	10.899 ml / 788.0 mV	result:	10.90 ml

## Hints

If you have any questions concerning the application, you are welcome to contact us.

优莱博技术（北京）有限公司

地址：北京市朝阳区酒仙桥东路1号院M8号楼C厅3层301室

■ Tel. : 4008092068 ■ Fax : 4008092068-112

■ E-Mail : info@julabo.cn

■ Web: www.julabo.cn