

Titration of Indigo and Sodium hydrosulfite in Indigo Vat 40

date: 14.5.2013 page 1 of 6



Use

Indigo is a so-called vat dye, which means that it needs to be reduced to its water soluble leuco-form before dying. The reduced form is absorbed into the fibres, and when oxidized back to its blue form it stays within the fibre. Earlier the reduction and dyeing was done with fermentation. Nowadays, the most of the reduction has been done chemically by sodium dithionite.

Concentrations of indigo and dithionite (hydrosulfite) are titrated simultaneously automatic titration using potassium ferricyanide as titrant

Appliances

Titrator: 7000 M1/20 consists of

- Basic device
- Magnetic stirrer TM 235
- 20 mL exchange unit WA 20, with brown glass bottle for titrant complete

Accessories:

- Titration head Z 306
- Glass beaker 100 ml without sprout, tall form
- Nitrogen gas for purging

Electrodes

Electrode: SA Pt 7780 or Pt 62

Cable: L 1 A

date: 14.5.2013 page 2 of 6



Reagents

- Titrant: Potassium ferricyanide 0.05 mol/l
- NaOH 0.1 mol/l
- Dispersing agent solution such as Sera Sperse Setamol WS 5 %

Description

Potassium ferricyanide 0.05 mol/l solution

4.11 g $K_3[Fe(CN)_6]$ are weighed in a 250 ml volumetric flask and filled up with distilled water to the mark. The solutions have to be stored in a dark bottle. The solution can be used for 1 week.

Sample titration

In the 100 mL glass beaker (tall form) are added 50 ml of the NaOH 0.1 m solution + 1 ml of a 5 % dispersing agent solution. The titration head is placed on the beaker and the solution is then purged with nitrogen for 5 minutes. To this solution are then pipetted 5 or 10 ml of the Indigo sample solution. Load the suitable method and start the titration.

Important:

The solution must be yellow. If the solution is brown or even blue then the purging with nitrogen was not enough. Maybe the addition of some Hydrosulfite turns the colour of the sample solution into yellow.




date: 14.5.2013 page 3 of 6



This example was made with 5 ml sample volume:

Standard documentation



Method data

Method name: Hydro-Indigo 5 ml a End date: 25.01.12 Titration duration: 6 m 12 s

End time:

15:52:52

Titration data

Sample ID: Test 3

Start mV -733.2 mV End mV 165.4 mV

EQ1: 2.280 ml/-677.8 mV Hydrosulfit: 2.234 g/l

EQ2: 5.092 ml/-194.5 mV Indigo: 3.824 g/l

Calculation formula

Hydrosulfit: (EQ1*F1)-F2
Mol (M): 1.00000
Indigo: (EQ2-EQ1)*F3
Mol (M): 1.00000

Factor 1 (F1): 0.9800 Factor 2 (F2): 0.0000

Factor 3 (F3): 1.3600

date: 14.5.2013 page 4 of 6



Method

Method data

Method name: Hydro-Indigo 5 ml Method type: Automatic titration

Measured value: mV
Titration mode: Linear
Linear steps: 0.100 ml

Created at: 05/14/13 14:39:54 Last modification: 05/14/13 14:39:54

Damping settings: None Documentation: GLP

Measuring speed / drift: Normal: minimum holding time: 02 s

maximum holding time: 15 s

Measuring time: 02 s

Drift; 20 mV/min

 Initial waiting time:
 0 s

 Titration direction:
 Increase

 Pretitration:
 Off

 End value:
 300.0 mV

 EQ:
 On (1)

 Slope value:
 User-defined

Value: 1350

Dosing parameter

Dosing speed: 100.00 % Filling speed: 30 s Maximum dosing volume: 20.00 ml

Calculation formula

Hydrosulfit: (EQ1*F1)-F2

 Unit:
 g/l
 Decimal places:
 3

 Indigo:
 (EQ2-EQ1)*F3

Unit: g/l Decimal places; 3

 Factor 1 (F1):
 0.9800
 Factor 2 (F2):
 0.1000

 Factor 3 (F3):
 1.3600
 Statistics:
 Off

date: 14.5.2013 page 5 of 6



Notes

If you have any questions on the application, you can feel free to contact us..

优莱博技术(北京)有限公司

地址:北京市朝阳区酒仙桥东路1号院M8号楼C厅3层301室

■E-Mail : info@julabo.cn
Web: www.julabo.cn



date: 14.5.2013 page 6 of 6