

KRC50 Refrigerators for Chemicals

for storing and cooling of chemicals and hazardous substances

JULABO refrigerators for chemicals are designed for storing and cooling hazardous substances. The spark free interior prevents damage caused by spilled or evaporating chemicals.

Your advantages

- · With spark free interior
- Storing and cooling of hazardous substances
- Self-protecting control circuit
- Digital temperature display (LED)
- · Overload protection for cooling compressor with test button
- · Cut-off in case of disturbance with optical alarm signal
- · Working and safety sensors are protected against short circuits and disconnection



Technical Data

Technical Data	
Order No.	8800705
Category	Refrigerators for chemicals
Working temperature range (°C)	-2 +1 2
Temperature stability (°C)	±1
Setting / display resolution	0.1 °C
Temperature Display	LED
Ambient temperature	035
Weight (kg)	32
Included with each unit	1 glass insert, 2 sliding drawers
Power requirement V / Hz / A	230/50
Compressor overheating protection °C	105
Alarm signal	optical
Safety sensor control	PTC
Working sensor	PTC
Temperature selection/display	Analog, LED
Volumetric capacity liters	68
Inner dimensions W x L x H cm	42x29x44
Outer dimensions W x L x H cm	55x64x63
Available voltage versions	220 - 240V / 50 Hz